i-MyRun II

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Do not pour liquid onto Power supply nor electrode connectors, as electrical shock, shortout, or a fire could result.

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If moving Gel electrophoresis tank while containing buffer in it, separate the Power supply first before moving, unless hold the Power supply and the tank firmly with both hands, as electrode connectors will be stressed, equipment damage or dropping of product could result.

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Unplug the power plug of the Power cable from the outlet when product is not in use.

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Crossover function activates when in cases as buffer concentration becomes too high, or overcurrent event occurs due to buffer temperature rose over proper level (over 40 degrees Celsius [°C]/104 degrees Fahrenheit [°F]).

* Crossover function: Automated function to suppress sample aggregation. Sudden current change will be detected and E1 error (overcurrent error) will show in the display.

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Do not leave buffer in Gel electrophoresis tank for prolonged period, as degradation of electrode could result.

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Gel electrophoresis tank and Gel casting tray both have UV (Ultraviolet) light transparency, though prolonged exposure to UV light may cause discoloration or decrease of UV light transparency level. Always wear a visor or approved goggles to help vision and to protect your eyes and skin from UV light.

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Always wear suitable protective gloves and avoid skin contact with nucleic acid stain, as it may be hazardous to human body.

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
Gel electrophoresis tank is washable, though dry well after the wash, as electrical shock, a fire or equipment damage could result when using it wet.

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
How to attach / detach Power supply and Gel electrophoresis tank

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
When attaching Power supply, slide in a little tilted as shown in the figure.

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Power supply will detach from the Gel electrophoresis tank by rotating towards it.

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* If both Safety switches are not functioning, "E2" error will show up in the display.

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When timer is set to 0 (zero), "--" will be displayed for continuous run mode.

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Switches voltage in order of 50→75→100→120→135 V as button is pressed.

e.g. When setting to 50 V, only 50 V will light up. When setting to 120 V, all 50, 75, 100 and 120 V will light up.

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
When pressed, electrophoresis will start as programmed. During the run, lamp will light up.

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Place Gel casting stand in a horizontal position. Install Gel casting tray and Comb(s) on Gel casting stand.

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Mix necessary amounts of buffer and agarose, heat to dissolve agarose completely.

* Use of 0.5 × TBE or 0.5 × TAE buffer is recommended.

* For preparation of 4 mm thickness agarose gel, amount of 70 - 80 mL gel will be required. For 4 mm thickness agarose gel, 9 μ L (26 well comb) or 18 μ L (13 well comb) of sample can be applied.

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*Pouring in over 70 °C (158 °F) agarose gel into Gel casting system may cause deforming.

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Leave under room temperature for 15 - 20 minutes until agarose gel solidifies.

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After agarose gel has solidified completely, draw out Comb(s) gently.

❗ If you wish to store agarose gel you have prepared, be sure to cover the gel with buffer after gently drawing out Comb(s). Wrap the gel in plastic wrap and store it under low temperature.

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

Pour in buffer into Gel electrophoresis tank that has same concentration of what was used to prepare agarose gel. For 4 mm gel, c.a. 350 mL of buffer is required to perform gel electrophoresis.

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Remove Gel casting tray with agarose gel from Gel casting stand and set it in the center of Gel electrophoresis tank. Place Gel casting tray in such way that the negatively charged nucleic acid migrates in the direction of the "+" marked on the tank. Be careful not to let air bubbles enter the space between Gel electrophoresis tray and the tank.

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Put on Gel electrophoresis tank cover and press "PUSH" mark, which is located at the top right corner of the cover.

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After setting timer and voltage, press "RUN" button to start gel electrophoresis. Please double-check that "RUN" button lights up and air bubbles forming from platinum wires at both ends of the tank.

❗ Performing gel electrophoresis for prolonged time under condition of 100 V or more will raise the buffer temperature, as sample distortion could result.

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
When gel electrophoresis completes, turn off Power supply and unplug the power plug of the Power cable from the outlet.

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Due to UV (Ultraviolet) light transparency of Gel electrophoresis tank and Gel casting tray, the result can be observed using transilluminator without taking out the gel from the tank, if staining dye has been added into the agarose gel when preparing it.

❗ Always wear a visor or approved goggles to help vision and to protect your eyes and skin from UV light.

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
If staining dye has not been added into the agarose gel when preparing it, take out the gel from Gel casting tray and stain it for necessary time length.

❗ When using ethidium bromide for DNA staining, always wear gloves to avoid skin contact.

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When staining is completed, observe the result using transilluminator.

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電源 Input voltage and power	AC 100 - 240 V, 90 VA, 50 - 60 Hz
使用温度条件 Temperature of operation	0 - 35 °C (ただし結露がないこと) from 0 - 35 °C / 32 - 95 °F (avoid condensation to form at lower temperatures)
直流出力電源 Output constant voltage	DC 50, 75, 100, 120, 135 V, 54 VA
タイマー Timer Settings	0 - 99 分 (0分に設定した場合、「--」と表示され連続運転モードに切り替わります。) Timer settings from 0 - 99 min (When timer is set to 0 [zero], "--" will be displayed for continuous run mode)
外寸 (幅、高さ、奥行き) Dimensions (width, height, length)	215mm (W) × 65mm (H) × 240 mm (L)
材質 Materials	電気泳動槽、ゲルキャスティングトレイ：紫外線透過性樹脂 蓋：紫外線非透過性樹脂 Gel electrophoresis tank : UV light transparent resin Gel casting tray : UV light transparent resin Gel electrophoresis tank cover : UV light un-transparent resin

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If an error occurs, the error message will be displayed during execution depending on its error level.

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COSMO BIO Co., LTD.
COSMO BIO CO., LTD.