

SOD Activity Assay

KB-03-011

100/200/400 test (96 well plate)

BOCKit

A brand of  **BioQuoChem**

Index

Introduction	Pag. 1
Materials	Pag. 2
Assay Principle	Pag. 4
Reagent Preparation	Pag. 5
Assay Protocol	Pag. 6
Data Analysis	Pag. 7
Warranties and Limitation of Liability	Pag. 8



All chemicals should be handled with care



➤ This kit is for R&D use only

Introduction

Superoxide dismutases (SODs) are metallo enzymes that catalyse the dismutation of the superoxide anion to molecular oxygen and hydrogen peroxide and thus form a crucial part of the cellular antioxidant defense mechanism.

Excessive reactive oxygen species, especially superoxide anion ($O_2^{\bullet-}$), **play important roles in the pathogenesis of** many cardiovascular diseases, including hypertension and atherosclerosis. Superoxide dismutases (SODs) are the major antioxidant defense systems against $O_2^{\bullet-}$, **which** consist of three isoforms of SOD in mammals: the cytoplasmic Cu/ ZnSOD (SOD1), the mitochondrial MnSOD (SOD2), and the extracellular Cu/ZnSOD (SOD3), all of which require catalytic metal (Cu or Mn) for their activation.

Superoxide Dismutase Activity Assay Kit (Colorimetric) is a sensitive kit using WST-1 that produces a water-soluble formazan dye upon reduction with superoxide anion. The rate of the reduction with a superoxide anion is linearly related to the xanthine oxidase (XO) activity, and is inhibited by SOD. Therefore, the inhibition activity of SOD can be determined by a colorimetric method.

Materials

BQCKit SOD Activity Assay kit *KB03011-100 tests* contains:

Product	Quantity	Storage
Reagent A	2 mL	4°C
Reagent B	20 µL	4°C
Reagent C	40 mL	4°C
Reagent D	10 mL	4°C
SOD Standard (3U/mL)	50 µL	4°C

BQCKit SOD Activity Assay kit *KB03011-200 tests* contains:

Product	Quantity	Storage
Reagent A	4 mL	4°C
Reagent B	40 µL	4°C
Reagent C	80 mL	4°C
Reagent D	20 mL	4°C
SOD Standard (3U/mL)	100 µL	4°C


Materials

BQCKit SOD Activity Assay kit *KB03011-400 tests* contains:

Product	Quantity	Storage
Reagent A	8 mL	4°C
Reagent B	80 µL	4°C
Reagent C	160 mL	4°C
Reagent D	40 mL	4°C
SOD Standard (3U/mL)	200 µL	4°C

Assay Principle

Bioquochem Superoxide Dismutase Activity Assay kit utilizes a tetrazolium salt for detection of superoxide radicals generated by xanthine oxidase and hypoxanthine (Scheme 1). One unit of SOD is defined as the amount of enzyme needed to exhibit 50% dismutation of the superoxide radical.



Scheme 1. SOD inhibition assay mechanism

XOD and SOD Antagonism in the Generation of Formazan Dye. The conversion of xanthine and O_2 to uric acid and H_2O_2 by XOD generates superoxide radicals. The superoxide anions reduce a tetrazolium salt (WST-1) to a colored formazan product (WST-1 formazan) that absorbs light. SOD scavenges superoxide anions, thereby reducing the rate of formazan dye formation.

Reagent Preparation

For 100 assays*:

*For a number of assays different from 100, recalculate Reagents volumes.

Working Reagent:


Mix 2 mL of Reagent A with 38 mL of Reagent C.

Enzyme Solution:

Centrifuge Reagent B in a microcentrifuge, and then mix by pipetting. Dilute it (12 μL with 2 mL of Reagent D).

Standard:

Prepare in 1.5 mL tubes, the following SOD standard Solutions with Phosphate Buffer pH 7.4 as diluent (not included): 1U/mL; 0.1U/mL; 0.05U/mL; 0.01U/mL and 0.001U/mL.



Scheme 2. Standard curve preparation

Assay Protocol

Considering a 96 well plate:

1. Pipette 20 μL of sample to each sample and Blank C well (you must prepare one Blank C well for each sample or standard).
2. Pipette 20 μL of ddH₂O to Blank A and Blank B wells.
3. Add 200 μL of the Working Solution (previously prepared) to each well.
4. Add 20 μL of Reagent D to Blank B and to each Blank C.
5. Add 20 μL of the Enzyme Solution (previously prepared) to Blank A well and to each sample well (Is important to add Enzyme Solution at the same time in all wells, using for example a multichannel pipette). Then mix the plate thoroughly.
6. Incubate the plate at 37°C for 20 min.
7. Read the absorbance at 450 nm using a microplate reader.

	Blank A	Blank B	Sample/Standard	Blank C
Sample/Standard			20 μL	20 μL
Working Solution	200 μL	200 μL	200 μL	200 μL
Enzyme Solution	20 μL			
Reagent D		20 μL		20 μL
ddH ₂ O	20 μL	20 μL		

Scheme 3. Instructions to fill the 96 well plate

Data Analysis

1. Plot the % of inhibition at 450 nm of standards as function of their final concentrations (Scheme 2).

$$\text{SOD Activity (\% inhibition)} = \left\{ \frac{[(\text{ABlank A} - \text{ABlank B}) - (\text{Asample/standard} - \text{ABlank C})]}{(\text{ABlank A} - \text{ABlank B})} \right\} \times 100$$

2. Create a standard curve by plotting the % inhibition at 450 nm for each standard against the concentration of SOD standards.
3. Calculate the SOD activity of the samples using the equation obtained from the linear regression of the standard curve replacing the % inhibition values for each sample.

Warranties and Limitation of Liability

Bioquochem shall not in any event be liable for incidental, consequential or special damages of any kind resulting from any use or failure of the products, even if Bioquochem has been advised of the possibility of such damage including, without limitation, liability for loss of use, loss of work in progress, down time, loss of revenue or profits, failure to realize savings, loss of products of buyer or other use or any liability of buyer to a third party on account of such loss, or for any labor or any other expense, damage or loss occasioned by such product including personal injury or property damage is caused by Bioquochem's gross negligence. Any and all liability of Bioquochem hereunder shall be limited to the amounts paid by buyer for product.

Buyer's exclusive remedy and Bioquochem's sole liability hereunder shall be limited to a refund of the purchase price, or the replacement of all material that does not meet our specifications.

Said refund or replacement is conditioned on buyer giving written notice to Bioquochem within 30 days after arrival of the material at its destination.

Expiration date: 1 year from the date of delivery

For further details, please refer to our website www.bqckit.com.

