

Cell Growth Determination Kit (MTT)

KF-03-001 A/B

500/2000/5000 test (96 well plate)

BOCKit

A brand of  **BioQuoChem**

Index

Introduction	Pag. 1
Materials	Pag. 2
Assay Principle	Pag. 3
Sample Preparation	Pag.4
Assay Protocol	Pag. 5
Warranties and Limitation of Liability	Pag. 7



All chemicals should be handled with care



This kit is for R&D use only

Introduction

Cell proliferation has been shown to have multiple functions in development and pattern formation, including roles in growth, morphogenesis, and gene expression.

Methods commonly used for this purpose are hemocytometer counting, determination of protein content, wet or dry weight measurement, and determination of the optical density (OD). While hemocytometer counting, and protein determination have the disadvantage of being time-consuming and tedious, the measurement of wet or even dry weight is not practical for very small culture volumes.

An alternative method is based on the transformation and colorimetric quantification of MTT [3-(4,5-dimethylthiazol-2-yl)-2,5-diphenyltetrazolium bromide]. The respiratory chain and other electron transport systems reduce MTT and other tetrazolium salts and thereby form non-water-soluble violet formazan crystals within the cell. The amount of these crystals can be determined spectrophotometrically and serves as an estimate for the number of mitochondria and hence the number of living cells in the sample. These features can be taken advantage of in cytotoxicity or cell proliferation assays, which are widely used in immunology, toxicology, and cellular biology.

Materials

BQCKit MTT Assay kit *KF03001 A* contains:


Product	Quantity	Storage
MTT Solution		
5000 test	1 bottle	4°C
2000 test	1 bottle	
500 test	5 vials	
MTT Solvent		
5000 test	1 bottle	RT
2000 test	2 bottles	
500 test	1 bottle	

BQCKit MTT Assay kit *KF03001 B* contains:

Product	Quantity	Storage
MTT Solution		
5000 test	1 bottle	4°C
2000 test	1 bottle	
500 test	2 vials	
MTT Solvent		
5000 test	1 bottle	RT
2000 test	2 bottles	
500 test	1 bottle	

Assay Principle

The principle of the MTT assay is based in mitochondrial activity. For viable cells, mitochondrial activity is constant and thereby an increase or decrease in the number of viable cells is linearly related to mitochondrial activity. The mitochondrial activity of the cells is reflected by the conversion of the tetrazolium salt MTT into formazan crystals, which can be solubilized for homogenous measurement. Thus, any increase or decrease in viable cell number can be detected by measuring formazan concentration spectrophotometrically using a plate reader at 570 nm vs. 690 nm.



Scheme 1. *Reduction of MTT to formazan*

Sample Preparation

This protocol is for a 96 well format. Volumes of culture cells, media and reagents may differ from the format described below.

1. Add 100 μL of culture cells to each well at an appropriate density. Include one set of wells with medium but no cells (control).
2. Incubate the cells overnight.
3. Treat cells on day two with agonist, inhibitor or drug ($V_f = 100 \mu\text{L}$) or change culture media if no treatments are required.
4. After the incubation time (drug and cell-dependent), follow the protocol described in Assay Protocol section.

Assay Protocol

Performing the assay

After sample preparation, follow next steps:

NOTE: These volumes are for a 96 well plate. For other sizes extrapolate the reagents volume.

1. Add 10 μ L of MTT solution to each well (10% of the culture media volume).
 2. Incubate for 4 hours at 37°C in a culture hood. The optimal incubation time may differ in each assay.
- From here on, there are two different protocols depending on the reference:

KF03001 A:

3. After the incubation period, remove culture medium from the culture hood and dissolve the resulting formazan crystals^{a,b}.
4. Cover and agitate 96 well plate on an orbital shaker for 15 minutes. (Within 1 hour of MTT solvent addition).
5. Read absorbance at 570 nm.
 - a. If cells are attached to culture vessels growth surface, remove and discard the culture media. Add MTT solvent in an amount equal to the original culture volume.
 - b. If cells are not attached, add MTT solution directly to the culture media in an amount equal to the original culture volume.

Assay Protocol

© 2015 ThermoFisher Scientific. All rights reserved. For research use only. For more information, visit www.thermo.com.

KF-03-001 B:

3. After the incubation period, dissolve the resulting formazan crystals with the MTT solvent in an amount equal to the culture volume.
4. Cover and agitate 96 well plate on an orbital shaker for 15 minutes. (Within 1 hour of MTT solvent addition).
5. Read absorbance at 570 nm.

Warranties and Limitation of Liability

Bioquochem shall not in any event be liable for incidental, consequential or special damages of any kind resulting from any use or failure of the products, even if Bioquochem has been advised of the possibility of such damage including, without limitation, liability for loss of use, loss of work in progress, down time, loss of revenue or profits, failure to realize savings, loss of products of buyer or other use or any liability of buyer to a third party on account of such loss, or for any labor or any other expense, damage or loss occasioned by such product including personal injury or property damage is caused by Bioquochem's gross negligence. Any and all liability of Bioquochem hereunder shall be limited to the amounts paid by buyer for product.

Buyer's exclusive remedy and Bioquochem's sole liability hereunder shall be limited to a refund of the purchase price, or the replacement of all material that does not meet our specifications.

Said refund or replacement is conditioned on buyer giving written notice to Bioquochem within 30 days after arrival of the material at its destination.

Expiration date: 1 year from the date of delivery

For further details, please refer to our website www.bqckit.com.

