

For research use only. Not for clinical diagnosis.

Catalog No. BAM-02-011-EX

Taq DNA Polymerase, Economy

BACKGROUND

Thermus aquaticus DNA polymerase (Taq DNA polymerase) gene was expressed in **E.Coli** in large quantities and highly purified. The enzyme has thermostable DNA polymerase activity and the MW is 94 kDa, same as that of the natural enzyme.

■ This enzyme is suitable for PCR reactions; capable of amplifying DNA with various primers.

Applications: 1) High-throughput PCR

2) Colony PCR

3) Incorporation of dUTP, dITP, and fluorescence-labeled nucleotides

4) Primer extension

5) Addition of a single nucleotide (adenosine) at the 3'-blunt ends

Size: 200 U (5U/μl)

Concentration: 5 units/ul, where one unit is defined as the amount of enzyme that can incorporate 10

nmols of total dNTPs into an acid-insoluble material in 30 minutes at 74°C when

activated salmon sperm DNA was used as template/primer.

Form: 20mM Tris-HCl (pH 8.0), 100mM KCl, 0.1mM EDTA, 1mM DTT, 50% glycerol, 0.5%

Tween20, 0.5% Igepal CA-630

Quality Assurance: Greater than 95% of protein determined by SDS-PAGE (CBB staining) (Fig.1) The

absence of endonucleases and exonucleases was confirmed.

PCR Test: Good amplification result was obtained in PCR reaction using λDNA as a template

(Fig.2).

Reagents Supplied

with Enzyme:

10 x Reaction Buffer (Taq): 100mM Tris-HCI (pH 8.3), 500mM KCI, 15mM MgCl₂

Storage: Store at -20°C

References:

Related Products

| BAM-02-001-EX | Taq DNA Polymerase(+dNTPs) | |
|---------------|----------------------------|--|
| BAM-02-021-EX | Pfu DNA Polymerase(+dNTPs) | |





| General composition of PCR reaction mixture (total 50ul) | | | |
|--|---------------------------|--|--|
| Taq DNA polymerase (5 units | /ul) *0.25 ul | | |
| 10 x Reaction Buffer (<i>Taq</i>) | 5 ul | | |
| 2.5mM (each) dNTPs | 4ul | | |
| Template | <500ng | | |
| Primer 1 | 0.2 ~ 1.0uM (final conc.) | | |
| Primer 2 | 0.2 ~ 1.0uM (final conc.) | | |
| Sterile distilled water | up to 50ul | | |
| *Use of excess amount is not recommended. | | | |

| PCR c | condition | | | |
|---------------------------------|-----------|-----------|--|--|
| 98°C | 10 sec |) | | |
| 57°C | 30 sec | 25 cycles | | |
| 72°C | 8 min | J | | |
| (2 min in the case of 2 kb DNA) | | | | |

1: 2 kb 2: 4 kb 3: 6 kb 4: 8 kb




Fig.1SDS-PAGE of Taq DNA polymerase

Typical other BioAcademia supplier ${
m M}$




Fig.2 Amplification of λ DNA

For research use only. Not for clinical diagnosis.

Manufactured by BioAcademia,Inc.



COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

Phone: +81-3-5632-9617 FAX: +81-3-5632-9618