



cDNA Library, *S. cerevisiae*, Log Phase

BACKGROUND

This cDNA library (plasmid DNA) is constructed from *Saccharomyces cerevisiae*, strain S288C-derived poly(A)⁺ RNA at the log phase by the Linker-Primer method (Ref. 1) by Prof. H. Nojima of Osaka University. This library is unidirectionally cloned by using the oligo (dT)₁₈ linker primer which contains the restriction enzyme site of *NotI*, and *BamHI* (*BglI*)-*SmaI* adaptor.

The pLZ3 vector (shown below) used in this library can not replicate in *S. cerevisiae* but contains pUCori for replication in *E. coli*.

Applications:

PCR screening of known or unknown gene: Prepare the primers for the known or unknown gene (cDNA) and amplify the gene by PCR from this library followed by cloning to an appropriate vector.

Standard amplifying conditions: 35 cycles of PCR reactions using 10-100 ng of cDNA as a template. (Change the quantity of template and the number of cycles depending on the expression rate of mRNA of the objective gene.)

Size:

500 ng (40 ng/ul, 13ul) in 10 mM Tris-HCl-1mM EDTA (pH 7.5)

Quality:

- 1) Number of independent clones: 3.6×10^6
- 2) Average insert size : longer than 1 kb

Storage:

Store at -20°C

References:


Construction of this library is described in Supplementary data of Ref.3

- 1) Kobori M *et al* "Large scale isolation of osteoclast-specific genes by an improved method involving the preparation of a subtracted cDNA library." *Genes Cells* 3: 459-475 (1998) PMID: [9753427](#)
- 2) Tanaka S and Nojima H "Nik1: a Nim1-like protein kinase of *S. cerevisiae* interacts with the Cdc28 complex and regulates cell cycle progression." *Genes Cells* 1, 905-921 (1996) PMID: [9077450](#)
- 3) Tougan T, Okuzaki D, [Nojima H](#). Chum -RNA allows preparation of a high-quality cDNA library from a single-cell quantity of mRNA without PCR amplification. *Nucleic Acids Res.*, 36(15):e92, (2008) PMID:[18603591](#)

Note

* This library is to be used only by the purchaser. It is not allowed to amplify and transfer the library to a third person.

* Related products: human tissue specific cDNA libraries and cDNA libraries of model organisms.



; MCS(pLZ3)

```

          CpoI (3)   SauI (b)   MluI (5)
PstI (3) SacI (3)   ApaI (3)
SseI (3) -----
NNNCTGCA CCTGCAGGAGCTCGGACCGGGCCCTTAGGACGGTAATACGACTCCTATAGGGAATTCGACGCTAGATCTTAAGGGCGCCAAAGGGTTGGCCA
NNNG  ACGTGGACGCTCCGAGCCTGGCCCGGGAATCCTGCGCATTATGCTGAGTGATATCCCTTAAGCTGCAGATCTAGAATTCGCGCGGTCCCAACCGGT

          AatII (3) BglII (5)   AscI (5)   BalI (b)
          -----
          EcoRI (5) XbaI (5) AflII (5)   BstXI (5)
          -----
          T7 Promoter
          -----
          BstEII (5)
          -----
          NheI (5)
          -----
          SceI (3)
          -----
          NotI (5) T3 promoter
          -----
          SmaI (3)   NruI (b)   SacII (3)
          -----
          SplI (5)
          -----
          PacI (3)
          -----
          SacI (3)
          -----
CGTGGTAACCAACGGGGTGGCTAGCTAGGGATAACAGGGTAATATAGCGCCGCCCTTATAGTGGGGTTAATTTAAATCGTACGTCGCGGATTAATTAACCGCGGTGGAGCT CAAT
GCACCATTTGGTCCCAACCGATCGATCCCTATTGTCCCATTTATATCGCCGGCGGGAATCACTCCCAATTAATTTAGCATGCAGCGCTAATTAATTGGCGCCACC TCGACTTA

TCGCCCTATAGTGAGTCGTATTA -3'
AGGGGATATCACTCAGCATAAT -5'

```

Fig. Structure of pLZ3 and the restriction sites.

For research use only. Not for clinical diagnosis.

Manufactured by BioAcademia, Inc.



COSMO BIO Co., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

http://www.cosmobio.co.jp/index_e.asp E-mail: export@cosmobio.co.jp

Phone : +81-3-5632-9617

F AX : +81-3-5632-9618