Easy-WESTERN Super

Primary Antibody Detecting Reagent for Western Blotting

(For high sensitivity and strong signal detection)

User 's Manual

Beacle, Inc.
OKAYAMA JAPAN

Storage condition and important information before use

1. The Multi-Antibody Detection (MAD) reagent should be stored at -20 °C immediately after receiving, and should be returned to -20 °C freezer as soon as possible after the use in order to prevent inactivation.
2. Dilution Buffer (20 x concentrate) should be stored at 4 °C.
3. Marker Detection Reagent should be stored at 4 °C.

---- Contents ----

(1) Introduction	2
(2) Contents	2
(3) Related Products	2
(4) Preparation and the method for attached reagents	2
(5) Protocol for the use as a substitute for 2nd antibody	3
(6) Protocol for the use as reprobing enhancer	3
(7) Protocol for the use as simple signal enhancer	3
(8) Trouble shootings	4
(9) Contact infromations	4

Important Notices

1. The reagent should be used only for research purpose, and not for diagnostic or medical purpose.
2. It is important to keep these reagents in appropriate condition (as indicated above) to prevent their inactivation

(1)INTRODUCTION

Easy-WESTERN (EZW) is a primary antibody detecting reagent kit for western blotting. The kit is based on our original Multi-Antibody Detection (MAD) technology. EZW enables you to detect primary antibody without using secondary antibody, and gives you many advantages, such as high sensitivity detection, signal enhancement, quick one-step detection and simultaneous multi antigen detection and more that are not experienced by ordinary western blotting technique.

Two types of Easy-WESTERN kit are available, Easy-WESTERN Super and Easy-WESTERN Quick. Both kits are the same in their composition except for the manual attached. The manual in Easy-WESTERN Super is written such a way that users can easily perform detections with higher sensitivity or signal enhancement, while that in Easy-WESTERN Quick with quick operation and multi-antigens. Both manuals can be downloaded from our webpage.

(Principle): The main component of the kit is MAD reagent that is nano-size protein particles with high affinity to antibodies. Each particle is composed of about 100 antibody-binding proteins and is labeled with HRP of 50 molecules. Because of these properties, MAD reagent enables high sensitivity and quick detection.

Advantages

- 1.No need of secondary antibody: MAD reagent can detect most types of primary antibodies¹⁾, and you do not have to prepare many types of secondary antibodies.
- 2.Can get higher signal as compared ordinary method using 2nd antibody²⁾: You can detect weakly expressed antigen or reduce amount of precious primary antibody needed.
- 3.Can enhance weak signals detected by ordinary 2nd antibody method: You can enhance a weak signal once observed by 2nd antibody, just by processing the same membrane with EZW.
- 4.Easy to use: MAD reagent, the main component of EZW, can be used easily. Just exchange the secondary antibody to MAD reagent and apply to membranes.
 - 1) MAD reagent may not work well with mouse IgG1 and goat IgG. For Mouse IgG1 please use Mouse IgG enhancer that can cover the less sensitivity.
 - 2) EZW's performance depends on the type of antibody, and we do not warrant the higher sensitivity in all cases.

(2)PRODUCT LINEUP AND COMPOSITION (the manual applies to below products)

Product #	Product name	Content
BCL-EZS01	Easy-WESTERN Super basic	MAD reagent, Dilution buffer, for 50 tests
BCL-EZS02	Easy-WESTERN Super Marker detector set	Basic plus Marker detector, for 50 tests
BCL-EZS04	Easy-WESTERN Super Mouse enhancer set	Basic plus Mouse IgG enhancer, for 50 tests
BCL-EZS03	Easy-WESTERN Super full set	Basic plus Marker detector and Mouse IgG enhancer, for 50 tests

(3)Related products

Product #	Product name	description
BCL-EZM01	Marker detector	for Easy-WESTERN (for 50 test)
BCL-EZE01	Mouse IgG enhancer	for Easy-WESTERN (for 50 test)
BCL-EZB01	Dilution buffer	for Easy-WESTERN (x20, conc., 50mL)
BCL-125A	Signal Booster Solution A (250mL)	Enhancer for antibody-antigen reaction

(4) PREPARATION AND METHODS FOR ATTACHED REAGENTS

1. Prepare dilution buffer by diluting 1 part of 20 x concentrate with 19 parts of DW. Before use, add skim milk at 5% to the diluted buffer (5% skim milk/dilution buffer)
2. Prepare TBS-T(150mM NaCl, 10mM Tris-HCl, 0.1% Tween-20, pH 7.6).
3. For membrane blocking we recommend skim milk (not included in the kit), and this manual describes the method using skim milk, but any other blocking reagent may be possible to use.

4. Marker detector is used to detect molecular markers, that are designed to be detected by secondary antibodies (such as MagicMark XP, Invitrogen), by MAD reagent.
 - To use, add 1/10000 volume of Marker detector into primary antibody reaction mixture.
 - The remaining procedure is the same as indicated in the below instructions.
5. Mouse IgG enhancer is used to enhance weak signals detected by using mouse IgGs such as IgG1
 - To use add Mouse IgG enhancer (dilution factor 2000) into 5% skim milk-containing dilution buffer with which MAD reagent is diluted.
 - The remaining procedure is the same as indicated in the below instructions.
 - We recommend pre-examination when you use this reagent in methods other than general method where MAD reagent is merely used as a substitute of secondary antibody.

(5) METHOD TO USE AS A SUBSTITUTE FOR SECONDARY ANTIBODY

This is the most general way of usage and best suited to get higher sensitivity and strong signal.

1. Separate protein sample(s) by general SDS-PAGE procedure.
2. Transfer protein to PVDF membrane from gel by general protocol.
3. Block with 5 % skim milk in TBS-T for 1 hour at room temperature.
4. Wash membrane with TBS-T (5 min x 3 times)
5. Incubate the membrane with primary antibody diluted in 1% skim milk/TBS-T for 1 hr at RT. (The primary antibody's dilution should follow the supplier's recommendation at the first trial. Appropriate dilution should be determined then. Signal Booster, solution A, can be used to further strength the signal)
6. Wash membrane with TBS-T (5 min x 3 times)
7. Dilute MAD reagent to 1/2000 in 5% skim milk/dilution buffer and incubate membrane in it for 1 hr at RT. MAD reagent can be used at 1/4000 to 1/1000 dilution.
8. Wash membrane with TBS-T (5 min x 3 times)
9. Detect signal by commercial HRP detection reagent.

(6) METHOD FOR ENHANCING SIGNALS BY REPROBING

This is a method to enhance weak signals once detected by ordinary detection using the same membrane.

1. Membrane with weak signal detected by ordinary 2nd antibody method. (Membrane have to be still wet with buffer solutions, dried membrane cannot be used)
2. Wash membrane with TBS-T (5 min x 3 times)
3. Incubate membrane with MAD reagent for 1 hr at RT as described in (5)-7.
4. Wash membrane with TBS-T (5 min x 5 times)
5. Detect signal by commercial HRP detection reagent.

(7) METHOD FOR ENHANCING SIGNALS DETECTED BY 2ND ANTIBODY

This is a method to enhance 2nd antibody-detected signals as the 3rd reaction

2. Perform the same procedures as described in (5) 1 to 6.
3. Incubate membrane with HRP-conjugated secondary antibody.
4. Wash membrane with TBS-T (5 min x 3 times)
5. Incubate membrane with MAD reagent for 1 hr at RT as described in (5)-7.
6. Wash membrane with TBS-T (5 min x 5 times)
7. Detect signal by commercial HRP detection reagent.

(8) TROUBLE SHOOTINGS

Trouble	Causes and Measures
Weak signal	1. Low amount of antigen: Increase the concentration of antigen
	2. Low primary antibody concentration: increase the concentration of primary antibody.
	3. Low protein transfer: Increase the electric current or transfer time.
	4. Strong blocking: over blocking reduced the signal intensity. Reduce the time or lower the concentration of blocking agents
	5. Primary antibody is either mouse IgG ₁ or Goat IgG: Consider the usage of EZW for enhancing signal method.
	6. Adherence of MAD reagent: When dilute MAD reagents in buffer without blocking agents, use low protein binding tubes.
whited out of signal band	7. Too much antigen or antibody: Too much signal inhibits luminescence, reduced the antigen or antibody to appropriate concentrations.
Many extra-bands	8. Too much antibodies: excess antibody resulted in increases non-specific signals. Reduced the antibody to appropriate concentration.
	9. Too much protein: Reduces the amount of protein in electrophoresis.
	10. not enough blocking: block the membrane with 5% skim milk /TBS-T for over 1 hrs.
	11. MAD reagent inactivated: Inappropriate storage of MAD reagents causes inactivation, and sometimes produces non-specific signal. Reduce the concentration of MAD reagents or replace with new one.
High background	12. not enough washing: increase the number and the duration of washing.
	13. High antibody concentration or long incubation: When high background with enough signal band visible, reduce the antibody concentration or reaction time with antibody.
	14. Too much MAD reagent: Reduce the concentration of MAD reagents.
	15. When using antigen-antibody reaction enhancer, insufficient washing causes high background: Increase the number and the duration of washing.
By simultaneous detection, one band shown weak signal.	16. One primary antibody weakly binds to antigen or MAD reagent: increase the concentration of the antibody with weak signal.
	17. Insufficient washing of primary antibody: Increase the number and the duration of washing.
	18. Washing of multiple primary antibodies conducted in the same chamber: A trace amount of antibody by reacting with MAD reagent resulted in weaker signal in non-balanced signals one, Please wash membrane separately.
No signal with one-step method	19. With not purified antibody, such as serum, one-step method may not work due to contaminated IgGs. Use purified antibodies or detect by two-step method.

The Q&A contains items related to Easy-WESTERN Quick

(9) CONTACT

Beacle Inc. [Manufacturer]

5303 Haga, Kita-Ku, Okayama 701-1221, Japan

Tel/FAX: +81-86-268-8091

E-mail: technical-support@beacle.com

Website: www.beacle.com

【Distributor】

COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN

URL: <http://www.cosmobio.com>

e-mail: export@cosmobio.co.jp

[International/Outside Japan] Phone: +81-3-5632-9617 [国内連絡先] Phone: 03-5632-9610

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