## 1. Introduction

CellEase®Tissue II is the reagent kit for extraction of genomic DNA from tissue cells rapidly and efficiently. You can prepare PCR (Polymerase Chain Reaction) grade template DNA with simple protocol (mix with reagent and test samples and then incubate without any purification steps).

## 2. Contents

Reagent A	: 1.0ml
Reagent B	: 1.0ml
Instruction manual	: 1
Volume	: 50 reactions
Store at	: 4°C

## 3. Principle

The reagent A disrupts the cells and stabilizes the genomic DNA from the samples. The reagent B degrades the cell extracts quickly. You can prepare the DNA from the cells efficiently in short time. These reagents don't have any inhibitor for PCR and also inactivate the PCR inhibitors from the cell debris. Thus, you can use the DNA samples with CellEase® Tissue II directly to PCR. These CellEase® Tissue II reagents don't include any toxic or harmful chemical compounds.

## 4. Use

Mouse tail, animal tissue, cultured cells, etc.

## 5. Preparation

Prepare the CellEase® mixture before use. Mix well the reagent A, B and distilled water (1:1:3). The mixture is not able to store for long time.

**Table 1. Preparation of CellEase® mixture**

Number of samples	Reagent A (μl)	Reagent B (μl)	Distilled water (μl)
1	20	20	60
5	100	100	300
10	200	200	600
25	500	500	1500
50	1000	1000	3000

## 6. Standard protocol (DNA template for PCR)

1) Cut samples into appropriate size and put in the micro-tube (The optimal sample amount will vary with the type of samples).

- 2) Add CellEase® mixture into the test samples.
- 3) Incubate the samples at 72°C for 6 minutes.
- 4) Then, continually incubate the samples at 94°C for 3 minutes.
- 5) Take appropriate amount of samples and use as a template DNA of PCR.



Fig. 1. The Schematic diagram for the extraction of DNA

## 5. Caution

- 1) Please use samples without including too much skin, hair, bone etc.
- 2) By using finely cut samples, genomic DNA can be more efficiently extracted.
- 3) You can adjust the volume of CellEase® mixture as amount of sample volume.
- 4) If it isn't enough DNA extracts with standard protocol. You can use a little longer incubation time (up to 60min) on the step of 72°C.
- 5) The rate of the mixture of reagent A, B and distilled water should be keep correctly.
- 6) The DNA extracts should be added less than 20% of the total volume of PCR.
- 7) The DNA extracts can be kept at 4°C for a few days. Alternatively it would be stored at -20°C until use.
- 8) If necessary, remove debris in DNA extraction fraction by using centrifugation.
- 9) CellEase®Tissue II is for research use only. Please don't use CellEase®Tissue II to the other purpose.

M	Marker (100bp ladder)	
1	DNA extract by using CellEase Tissue II	conc.
2		× 10
3		× 100
4	DNA extract by using conventional CellEase	conc.
5		× 10
6		× 100

※ The protocol of conventional CellEase and CellEase II kit were followed by the instruction manual respectively.

※ The thickness of the DNA bands were depending upon the amount of test samples and a parts of tissue including bone, skin or fat.

The clear DNA bands were detected from more than × 10 dilution of DNA extracts by using CellEase Tissue II


**Protocol of CellEase Tissue II**

- 1) Cut 5mm of the mouse tail (c.a.10mg) and put it in the micro test tube (ordinary use 0.2ml or 0.5ml tubes for PCR).
- 2) Mix the Reagent A, B and distilled water.  
(20 ul CellEase A, 20 ul CellEase B, 60 ul distilled water)
- 3) Add 100 ul of the mixture to the samples.
- 4) Incubate at 72°C for 6 minutes. Then incubate at 94°C for 3 minutes.
- 5) Transfer 6ul of extracts to PCR reaction mixture and amplify the target DNA fragment.

**PCR reagent**

6.0 ul	Test sample
5.0 ul	× 10 buffer (+Mg <sup>2+</sup> )
5.0 ul	dNTPs
1.0 ul	Forward Primer (10pmol/ul)
1.0 ul	Reverse Primer (10pmol/ul)
0.5 ul	Ex Taq DNA polymerase (5 U/ul, Takara, Japan)

Fill up to 50ul by distilled water

※ We can provide application data for the other kind of samples. Please don't hesitate to contact us.

※ CellEase® is the registered trademark of Biocosm Inc.

*For research use only. Not for clinical diagnosis.*

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