





New Tool For Protein Research

CloverDirect™

tRNA Reagents for Site-Directed Protein Functionalization

Expression of proteins with unnatural amino acids using four-base codon (CGGG) or amber stop codon (UAG)

CloverDirect™ tRNA Reagents for Site-Directed Protein Functionalization allow the incorporation of unnatural amino acids at defined positions of proteins using *in vitro* translation system. Various unnatural amino acids are available for fluorescence labeling, biotin labeling, post-translational modification, and unnatural mutagenesis of proteins.



Principle

If CGGG or UAG codon is recognized by unnatural aminoacyl-tRNA, full-length protein with unnatural amino acid is successfully synthesized. On the contrary, if CGGG or UAG codon is recognized by natural translation factor (arginyl-tRNA_{CCG} or release factor 1), the protein synthesis is terminated. Therefore, the translation product obtained as a full-length protein contains the unnatural amino acid at 100% efficiency.

Features

Accurate and Quantitative:

Incorporation position of unnatural amino acids is defined by CGGG four-base codon or UAG amber codon. An unnatural aminoacyl-tRNA recognizes the CGGG or UAG codon during translation to introduce unnatural amino acids into proteins in a site-directed and quantitative fashion.

Fast and Easy:

Proteins with unnatural amino acids can be obtained within a few hours just by adding CloverDirect™ reagents and DNA or mRNA having CGGG or UAG codon to an *in vitro* translation system.

Flexible:

Various unnatural amino acids containing fluorescent groups, biotin, PEG, photo-cross-linker, etc. are available.

Product Contents

- Unnatural aminoacyl-tRNA X 1
- tRNA dissolving buffer X 1

One tube contains unnatural aminoacyl-tRNA sufficient for 300 μL of *in vitro* translation reaction.



COSMO BIO CO., LTD.

www.cosmobio.com


Applications of site-directly labeled proteins

Site-directly fluorescent-labeled proteins are available to the following analyses.

- Protein interaction analysis using Single Molecule fluorescence detection system (MF20 / FluoroPoint-Light; OLYMPUS).
- Conformation analysis of protein by inter- or intra-molecular fluorescence resonance energy transfer (FRET).
- Functional analysis by cell imaging.
- Interaction analysis by protein array.
- Expression analysis by fluorescent detection of SDS-PAGE.

Expression of site-directly labeled proteins

Fluorescent- and biotin-labeled unnatural amino acids are incorporated into eight prokaryote and eukaryote proteins. The site-directly fluorescent-labeled proteins can be visualized on SDS-PAGE using a laser-based fluorescence scanner. The proteins are also detectable by an antibody against tag peptide or biotin.



2UAG : UAG codon is inserted after initiator AUG codon.
 ProX tag : ProX tag is fused to the N-terminus.
 Applied volume: 0.25 μ L of translational reaction mix
 Fluorescence images (Top) are visualized with Ex and Em wavelengths listed below:
 HiLyte Fluor488 Ex : 488 nm / Em : 520 nm
 TAMRA Ex : 532 nm / Em : 580 nm
 ATTO633 Ex : 635 nm / Em : 670 nm
 Western blotting (Bottom) are analyzed by anti-His tag antibody (for fluorescent amino acids) and anti-biotin antibody (for biotin)

References

- 1) Daisuke Kajihara, Ryoji Abe, Issei Iijima, Chie Komiyama, Masahiko Sisido, Takahiro Hohsaka *Nature Methods*, **3**, 923-929 (2006).
- 2) Takayoshi Watanabe, Norihito Muranaka, Issei Iijima, Takahiro Hohsaka *Biochem. Biophys. Res. Commun.*, **361**, 794-799 (2007).
- 3) Hikaru Taira, Yosuke Matsushita, Kenji Kojima, Kaori Shiraga, Takahiro Hohsaka *Biochem. Biophys. Res. Commun.*, **374**, 304-308 (2008).
- 4) Takahiro Hohsaka, Daisuke Kajihara, Yuki Ashizuka, Hiroshi Murakami, Masahiko Sisido *J. Am. Chem. Soc.*, **121**, 34-40 (1999).

Product List

Site-Directed Fluorescence Labeling

Catalog No.	Description	Size
PRX-CLD1001	CloverDirect™ CR110-X-AF (amber)	300 μ L
PRX-CLD2001	CloverDirect™ CR110-X-AF (amber)	5×300 μ L
PRX-CLD2002	CloverDirect™ CR110-X-AF (CGGG)	5×300 μ L
PRX-CLD01	CloverDirect™ HiLyte Fluor 488-AF (amber)	300 μ L
PRX-CLD05	CloverDirect™ HiLyte Fluor 488-AF (amber)	5×300 μ L
PRX-CLD2004	CloverDirect™ HiLyte Fluor 488-AF (CGGG)	5×300 μ L
PRX-CLD02	CloverDirect™ TAMRA-X-AF (amber)	300 μ L
PRX-CLD06	CloverDirect™ TAMRA-X-AF (amber)	5×300 μ L
PRX-CLD03	CloverDirect™ ATTO633-AF (amber)	300 μ L
PRX-CLD07	CloverDirect™ ATTO633-AF (amber)	5×300 μ L
PRX-CLD2008	CloverDirect™ ATTO633-AF (CGGG)	5×300 μ L
PRX-CLD1009	CloverDirect™ ATTO655-X-AF (amber)	300 μ L
PRX-CLD2009	CloverDirect™ ATTO655-X-AF (amber)	5×300 μ L
PRX-CLD1010	CloverDirect™ ATTO655-X-AF (CGGG)	300 μ L
PRX-CLD2010	CloverDirect™ ATTO655-X-AF (CGGG)	5×300 μ L

Site-Directed Biotin Labeling

Catalog No.	Description	Size
PRX-CLD2101	CloverDirect™ Biotin-AF (amber)	5×300 μ L
PRX-CLD2102	CloverDirect™ Biotin-AF (amber)	5×300 μ L
PRX-CLD2103	CloverDirect™ Biotin-X-AF (amber)	5×300 μ L
PRX-CLD2104	CloverDirect™ Biotin-X-AF (CGGG)	5×300 μ L
PRX-CLD04	CloverDirect™ Biotin-XX-AF (amber)	300 μ L
PRX-CLD08	CloverDirect™ Biotin-XX-AF (amber)	5×300 μ L
PRX-CLD2106	CloverDirect™ Biotin-XX-AF (CGGG)	5×300 μ L

Site-Directed Post-Translational Modification

Catalog No.	Description	Size
PRX-CLD2201	CloverDirect™ Tyr(PO3H2) (amber)	5×300 μ L
PRX-CLD2202	CloverDirect™ Tyr(PO3H2) (CGGG)	5×300 μ L
PRX-CLD2203	CloverDirect™ Lys(Me) (amber)	5×300 μ L
PRX-CLD2204	CloverDirect™ Lys(Me) (CGGG)	5×300 μ L
PRX-CLD2205	CloverDirect™ Lys(Me2) (amber)	5×300 μ L
PRX-CLD2206	CloverDirect™ Lys(Me2) (CGGG)	5×300 μ L
PRX-CLD2207	CloverDirect™ Lys(Ac) (amber)	5×300 μ L
PRX-CLD2208	CloverDirect™ Lys(Ac) (CGGG)	5×300 μ L

Site-Directed Unnatural Mutagenesis

Catalog No.	Description	Size
PRX-CLD2301	CloverDirect™ PEG4-AF (amber)	5×300 μ L
PRX-CLD2302	CloverDirect™ PEG4-AF (CGGG)	5×300 μ L
PRX-CLD2303	CloverDirect™ PEG8-AF (amber)	5×300 μ L
PRX-CLD2304	CloverDirect™ PEG8-AF (CGGG)	5×300 μ L
PRX-CLD2305	CloverDirect™ PEG12-AF (amber)	5×300 μ L
PRX-CLD2306	CloverDirect™ PEG12-AF (CGGG)	5×300 μ L
PRX-CLD2321	CloverDirect™ BPA (amber)	5×300 μ L
PRX-CLD2322	CloverDirect™ BPA (CGGG)	5×300 μ L
PRX-CLD2323	CloverDirect™ AcPhe (amber)	5×300 μ L
PRX-CLD2324	CloverDirect™ AcPhe (CGGG)	5×300 μ L
PRX-CLD2331	CloverDirect™ AzoAla (amber)	5×300 μ L
PRX-CLD2332	CloverDirect™ AzoAla (CGGG)	5×300 μ L

We provide custom services for the synthesis of unnatural aminoacyl-tRNAs and the expression of proteins with unnatural amino acids.

manufacturer:



COSMO BIO Co., LTD.

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
 KOTO-KU, TOKYO 135-0016, JAPAN
 TEL : (81)3-5632-9617
 FAX : (81)3-5632-9618
 e-mail : export@cosmobio.co.jp
 URL : www.cosmobio.com

www.cosmobio.com