

Anti Human Collagen IV α Chain Monoclonal antibodies

For staining human renal and skin biopsy sections easily and rapidly!


Epitope-defined monoclonal antibodies against the NC1 domains of $\alpha 1 - \alpha 6$ (IV) chains of type-IV collagen are against the triple-helical domain. These epitope-defined monoclonal antibodies against different site of the collagen molecule are useful for research on type-IV collagen and diagnosis of hereditary disease related to this collagen because of their reliability.

Line UP : Anti Human collagen IV, Rat monoclonal

Clone	Subclass	Specificity	$\alpha 1$	$\alpha 2$	$\alpha 3$	$\alpha 4$	$\alpha 5$	$\alpha 6$	Epitope
H53	IgG2a κ	imperfection III	—	—	—	—	+	—	DVEF with an amino acid sequence ¹⁾
B51	IgG2a	NC1 domain ³⁾	—	—	—	—	+	—	unknown
H25	IgG1 κ	imperfection XIII	—	+	—	—	—	—	EAIQP with an amino acid sequence ¹⁾
H52	IgG2b κ	NC1 domain ⁴⁾	—	—	—	—	+	—	SKPQSETL with an amino acid sequence


The chromatic figure of normal human kidney stained with the product.

H53




Applicable to staining for basement membrane of glomerular, Bowman's capsular and partial renal tubular.

B51




H25



Applicable to staining for basement membrane of glomerular, Bowman's capsular, renal tubular and capillary vessel.


H52



Human cryostat section stained with monoclonal antibody H52. Glomerular basement membrane (BM), part of tubular BM, and Bowman's capsular BM are stained.

H53, B51, H25 Cocktail


This product can be used for staining human renal and skin biopsy sections easily and rapidly as two monoclonal antibodies are conjugated with two different fluorochromes. Cocktail (H53+B51+H25) is useful especially for research of Alport syndrome*



Double exposure of normal human kidney section stained with the product.

FITC fluorescence ($\alpha 5$ chain) is observed in the GBM, part of the TBM and Bowman's capsular BM. BMs were two fluorochromes are present look orange to yellow.

*Alport syndrome, an inherited disease, shows the absence or reduction of the $\alpha 5$ (IV) chain in the glomerular, tubular and Bowman's capsular basement membranes.



Double exposure of Alport human kidney section stained with the product.

Because no FITC fluorescence ($\alpha 5$ chain) is observed in the kidney of the patient with X-linked Alport syndrome, only Texas Red fluorescence ($\alpha 2$ chain) is confirmed.



Anti Human Collagen IV α Chain Monoclonal antibodies

Anti Human Collagen IV α Chain Monoclonal antibody series

Specification	Clone	Labeling	Cat. No.	Specificity	Size	Acid treatment	Storage	Application
Cocktail*	H53	FITC	SGE-CFT45325	α 5	1 ml	unnecessary	4°C in dark	Immunofluorescence (Frozen section)
	B51	FITC		α 5				
	H25	Texas Red		α 2				
Single	H53	FITC	SGE-CFT453	α 5	500 μ l	unnecessary	4°C in dark	
		No labeled	SGE-C453				-20°C	
	B51	FITC	SGE-CFT451	α 5		unnecessary	4°C in dark	
		No labeled	SGE-C451				-20°C	
	H25	Texas Red	SGE-CFT425	α 2		unnecessary	4°C in dark	
		No labeled	SGE-C425				-20°C	
	H52	No labeled	SGE-C452	α 5		necessary	-20°C	WB Immunofluorescence (Frozen section)

* Useful for a study of the Alport syndrome

【Appearance】 Labeled antibody : affinity purified from supernatant.
No labeled antibody : Collected from supernatant by high-density culture and diluted with serum GIF medium. (Containing 10% FBS).

Citation

- Samar M. Said Brooke McCann *et al.* Negative Staining for COL4A5 Correlates With Worse Prognosis and More Severe Ultrastructural Alterations in Males With Alport Syndrome. *Kidney Int Rep.* 2 (1): 44-52. (2017)
- Malone AF, Miner JH *et al.* Functional assessment of a novel COL4A5 splice region variant and immunostaining of plucked hair follicles as an alternative method of diagnosis in X-linked Alport syndrome. *Pediatr Nephrol.* 32 (6): 997-1003. (2017)
- Nozu K, Makita J *et al.* X-linked Alport syndrome caused by splicing mutations in COL4A5. *Am Soc Nephrol.* 9 (11): 1958-64. (2014)
- Matsubara S, Suzuki M.J *et al.* Pregnancy complicated with Alport syndrome: a good obstetric outcome and failure to diagnose an infant born to a mother with Alport syndrome by umbilical cord immunofluorescence staining. *Obstet Gynaecol Res.* 35 (6): 1109-14. (2009)
- Patey-Mariaud de Serre N, Knebelmann B *et al.* Collagen alpha5 and alpha2(IV) chain coexpression: analysis of skin biopsies of Alport patients. *Kidney Int.* 72 (4): 512-6. (2007)
- Kharrat M, Saudi J *et al.* Autosomal dominant Alport's syndrome: study of a large Tunisian family. *Kidney Dis Transpl.* 17 (3): 320-5. (2006)

References

- Kagawa M *et al.* Epitope-defined monoclonal antibodies against type-IV collagen for diagnosis of Alport's syndrome. *Nephrol. Dial. Transplant.* 12: 1238-1241. (1997)
- Naito I *et al.* Relationship between COL4A5 gene mutation and distribution of type IV collagen in male X-linked Alport syndrome. *Kidney Int.* 50: 304-311. (1996)
- Borza DB *et al.* The NCI domain of collagen IV encodes a novel network composed of the alpha 1, alpha 2, alpha 5, and alpha 6 chains in smooth muscle basement membranes. *J. Biol. Chem.* 276: 28532-28540. (2001)
- Sado *et al.* Establishment by the rat lymph node method of epitope-defined monoclonal antibodies recognizing the six different α chains of human type IV collagen. *Histochem. Cell Biol.* 104: 267-275. (1995)

Producer : Shigei Medical Research Institute

All antibodies were prepared by the rat lymph node method developed by Shigei Medical Research Institute.



COSMO BIO Co., LTD.

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
TEL : (81)3-5632-9617
FAX : (81)3-5632-9618
e-mail : export@cosmobio.co.jp
URL : www.cosmobio.com



10193