



Anti-8-Nitro-cGMP

Background

8-Nitroguanosine is a nitrated nucleic acid which is formed from the reactions of guanosine with peroxyxynitrite, myeloperoxidase, nitrite, and peroxide. It is known that the nitration of guanine is enhanced in virus infection^{1,2}, bacterial infection^{3,4}, inflammatory disease⁵, cancers⁵, and diseases associated with smoking⁶. 8-Nitroguanosine is thought to be one of the makers of DNA damage caused by oxidative stress. Cyclic GMP (cGMP) is a second messenger that activates protein kinase G. 8-Nitro-cGMP (nitrated cGMP) has been identified as a nitrated cGMP formed in vivo³. 8-Nitro-cGMP can act as a unique second messenger distinct from cGMP to induce antioxidative adaptive responses^{3,7-9}. Mode of actions of 8-nitro-cGMP mainly relies on its adduction to protein cysteine residues called "protein S-guanylation", as a posttranslational modification.

Product type	Primary Antibody
Immunogen	8-Nitro-cGMP-conjugated bovine serum albumin
Raised in	Mouse
Clone number	1G6
Isotype	IgG1
Source	Ascites
Purification	Ion-Exchange Chromatography
Buffer	10 mM Tris-HCl (pH 7.4) + 50 mM NaCl (no preservatives added)
Concentration	1 mg/mL
Volume	50 ug
Label	Unlabeled
Specificity	8-Nitro-cGMP
Storage	Store at -70°C Aliquot to avoid cycles of freeze/thaw.

Recommended Dilutions	Immunocytochemistry, 1 ug/mL in Can get signal solution I Other applications have not been tested or not reactive. Optimal dilutions/concentrations should be determined by the end user.
-----------------------	---




Fig 1. Immunocytochemistry of 8-nitro-cGMP in C6 cells treated with authentic 8-nitro-cGMP. C6 cells (0.25×10^5 cells in 0.5 mL/well containing circle glass) were treated with 8-nitro-cGMP ($100 \mu\text{M}$) or untreated for 3 h in FBS-free DMEM. Cells were fixed with Zamboni solution (4% paraformaldehyde and 10 mM picric acid in 0.1 M phosphate buffer, pH 7.4) at room temperature (RT) for 1 h. After washing with PBS, the preparations were permeabilized with PBS containing 0.5% TritonX-100 at RT for 15 min, washed, and incubated in BlockAce at RT for 1 h. After washing with TBST three times, the preparations were incubated with 1G6 ($1 \mu\text{g}/\text{mL}$ in Can get signal solution I) at 4°C over night. After washing with TBST three times, the preparations were incubated with goat anti-mouse IgG (H+L), highly cross-adsorbed, HiLyte™ Fluor 555-labeled (AnaSpec, Inc, San Jose, CA) ($0.2 \mu\text{g}/\text{mL}$, 1:5000 in Can get signal solution II) at RT for 1 h. After washing with TBST three times and with distilled water one time, the preparations were mounted on glass slides using Prolong® Diamond Antifade Mountant with DAPI (Invitrogen; Carlsbad, CA).



Fig 2. Immunocytochemistry of 8-nitro-cGMP formed in immunologically stimulated C6. C6 cells (0.25×10^5 cells in 0.5 mL/well containing circle glass) were treated with LPS $10 \mu\text{g}/\text{mL}$ + Cytokines ($10 \text{ ng}/\text{mL}$ IL- 1β , $100 \text{ U}/\text{mL}$ TNF- α and $100 \text{ U}/\text{mL}$ IFN- γ) (all cytokines from R & D Systems) for 30 h. Cell fixation and immunocytochemistry protocols were same as that for Fig. 1.

References:

- 1) T. Akaike, S. Okamoto, T. Sawa, J. Yoshitake, F. Tamura, K. Ichimori, K. Miyazaki, K. Sasamoto and H. Maeda, 8-nitroguanosine formation in viral pneumonia and its implication for pathogenesis, *Proc. Natl. Acad. Sci. USA*, 100, 685-690 (2003).
- 2) J. Yoshitake, T. Akaike, T. Akuta, F. Tamura, T. Ogura, H. Esumi, and H. Maeda, Nitric oxide as an endogenous mutagen for Sendai virus without antiviral activity, *J. Virol.*, 78, 8709-8719 (2004).
- 3) T. Sawa, M. H. Zaki, T. Okamoto, T. Akuta, Y. Tokutomi, S. Kim-Mitsuyama, H. Ihara, A. Kobayashi, M. Yamamoto, S. Fujii, H. Arimoto, and T. Akaike, Protein S-guanylation by the biological signal 8-nitroguanosine 3',5'-cyclic monophosphate, *Nat. Chem. Biol.*, 3, 727-735 (2007).
- 4) M. H. Zaki, S. Fujii, T. Okamoto, S. Islam, S. Khan, K. A. Ahmed, T. Sawa, and T. Akaike, Cytoprotective function of heme oxygenase 1 induced by a nitrated cyclic nucleotide formed during murine salmonellosis, *J. Immunol.*, 182, 3746-3756 (2009).
- 5) Y. Terasaki, T. Akuta, M. Terasaki, T. Sawa, T. Mori, T. Okamoto, M. Ozaki, M. Takeya and T. Akaike, Guanine nitration in idiopathic pulmonary fibrosis and its implication for carcinogenesis, *Am. J. Respir. Crit. Care. Med.*, 174, 665-673 (2006).
- 6) M. Feelisch, Nitrated cyclic GMP as a new cellular signal, *Nat. Chem. Biol.*, 3, 687-688 (2007).
- 7) K. A. Ahmed, T. Sawa, T. Akaike, Protein cysteine S-guanylation and electrophilic signal transduction by endogeneous nitro-nucleotides, *Amino Acids*, 41, 123-130 (2011).
- 8) T. Sawa, H. Arimoto and T. Akaike, Regulation of redox signaling involving chemical conjugation of protein thiols by nitric oxide and electrophiles, *Bioconjugate Chem.*, 21, 1121-1129 (2010).

For research use only, Not for diagnostic use.



COSMO BIO CO., LTD.

【JAPAN】
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU, TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】
2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobiousa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobiousa.com