




## Anti TGF- $\beta$ 1 LAP-D (R58)

### (LAP Degrades C-Terminus side cut end R58)

#### Background :

TGF- $\beta$  is produced as a latent form in which 25 kD active TGF- $\beta$  is trapped by its pro-peptide called Latency Associated Protein (LAP). Upon receiving certain stimuli, a conformational change is induced in a latent complex to release the active TGF- $\beta$  from the complex. The resultant TGF- $\beta$  binds to cognate signaling receptors and exerts various physiological and pathological activities. This reaction is called TGF- $\beta$  activation reaction, which is known to be induced by binding of the latent complex to cell adhesion proteins such as thrombospondin and integrins, and/or by being cleaved by the action of proteases such as serine proteases, cysteine proteases, and MMPs in an organ and context-depending manner.

Kojima and his colleagues in Cellular Molecular Pathology Research Unit (currently, Center for Integrative Medical Sciences, Liver Cancer Prevention Research Unit), RIKEN, Japan identified that a serine protease, plasma kallikrein induces release and activation of TGF- $\beta$  by cleaving between 58Arg-59Leu within LAP and thereby participates in the pathogenesis of the liver diseases. The anti-TGF- $\beta$ 1 LAP-degradates (LAP-D) antibodies are useful to investigate the molecular mechanism of TGF- $\beta$  activation and its related diseases including liver fibrosis/cirrhosis and liver degeneration as tools to detect LAP-D.



Host Species:	Mouse
Form:	Liquid, PBS (pH 7.4), 0.05% NaN <sub>3</sub>
Volume:	100 $\mu$ g (1 mg/mL)
Specificity:	Recognizes C-terminus cut end of LAP degradates (LAP-D) R58 when latent TGF- $\beta$ is digested with Plasma Kallikrein (PLK).
Antigen:	R58 peptide [CGQILSKLR]
Clonality:	Monoclonal (clone # 18F9-16)
Isotype:	IgG3
Applications:	Immunohistochemistry (IHC): 10 $\mu$ g/mL * Optimal dilutions/concentrations should be determined by each researcher.
Purification method:	Purified from cell culture of serum-free medium by affinity column (Protein A)
Conjugation:	none
Storage condition:	Store below -20°C (below -70°C for prolonged storage) *Aliquot to avoid cycles of freeze/thaw

#### References:

1. LAP degradation product reflects plasma kallikrein-dependent TGF- $\beta$  activation in patients with hepatic fibrosis, Hara M., Kirita A., Kondo W. et al. Springerplus. May 1; 3: 221. PMCID: PMC4033717 (2014)

\* Anti TGF- $\beta$ 1 LAP-D (R58) was generated & licensed under RIKEN, Japan.

## Example Assay Data:

### 1. Immunohistochemistry (IHC) Staining




Figure 1. Immunohistochemistry (IHC) staining of liver section from a hepatitis model mice

【sample】 liver, paraffin sections (4 μm) from hepatitis model mice (positive control, treated) and sham operated mice (negative control, sham)

Primary Antibody: anti TGF-β1 LAP-D (R58) antibody (10 μg/mL)

#### 1. Sectioning

Paraffin section (4 μm) / adhesive coated glass slide

#### 2. Deparaffinizing

Xylene 5 min x 3

100% Ethanol 3 min x 2

90% Ethanol 3 min

80% Ethanol 3 min

70% Ethanol 3 min

50% Ethanol 3 min

MilliQ water

#### 3. Staining

epitope retrieval: Target Retrieval Solution pH 9.0 [Agilent Technologies, S236784-2], microwave 98 °C, 10 min

wash: PBST (0.1% tween 20), 5 min x 3

blocking: 10% donkey serum [ImmunoBioScience, IHR-8135] / PBST (0.1% tween 20), room temperature, 1 hr

primary antibody: anti TGF-β1 LAP-D (R58) antibody (10 μg/mL)

10% donkey serum / PBST (0.1% tween 20), 4 °C , overnight

wash: PBST (0.1% tween 20), 5 min x 3

secondary antibody: donkey anti-mouse IgG alexa555 [Thermo Fisher Scientific, A-31570], 1/1000, room temperature, 2 hrs

wash: PBST (0.1% tween 20), 5 min x 3

counter staining: DAPI

wash: PBST (0.1% tween 20), 5 min x 3

mounting

For research use only, Not for diagnostic use.



COSMO BIO Co., LTD.

2-20, Toyo 2-Chome, Koto-ku, Tokyo 135-0016, JAPAN

T: +81-3-5632-9617

F: +81-3-5632-9618

E: [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

W: [www.cosmobio.com](http://www.cosmobio.com)

日本のお客様

T: 03-5632-9610

F: 03-5632-9619

W: [www.cosmobio.co.jp](http://www.cosmobio.co.jp)