

Alginate 3D Cell Culture Kit

A three-dimensional (3D) cell culture kit using alginate gel.

-  **Easy to produce alginate gel beads and harvest the cultured cells from gel.**
-  **Detect the ability of cells to undergo anchorage-independent growth like as soft-agar culture.**
-  **Suitable for 3D culture including tumor cells and chondrocytes.**



Kit Components:


Sodium alginate solution (25 mL)
Calcium chloride solution (100 mL x2 bottles)
Sodium citrate solution (100 mL)
Plastic flexible needle (4 pieces)
24-well plate (4 pieces, for suspension culture)

Materials and Instruments Required, but Not Included:

Physiological saline (sterile)
22G Hypodermic needle
5 mL Syringe

Preparation of Alginate Beads

- (1) Prepare cell suspended alginate solution in 5mL syringe.
- (2) Drop the cell suspension into each well containing calcium chloride solution and produce alginate beads.
- (3) Wash alginate beads.
- (4) Dispense culture medium and let start the culture.




Alginate 3D Cell Culture Kit

A three-dimensional (3D) cell culture kit using alginate gel.


Cell Recovery from Alginate Beads

- (1) Remove the medium from the well.
- (2) Add sodium citrate solution and dissolve the alginate beads.
- (3) Centrifuge and harvest the cells as a precipitated pellet.


Reference Data




Porcine Chondrocytes cultured in alginate beads



HepG2 cell cultured in alginate beads



Growth curve of HepG2 cells cultured in alginate beads (Mean $\pm$ SD, n=3)



Effect of 5-FU on the growth of HepG2 cells in alginate beads culture

Product name	Product code	Quantity	Remarks
Alginate 3D Cell Culture Kit	CSR-ABC-KIT	1 kit	sterile, store below 4° C, Do not Freeze
Sodium alginate solution	CSR-ABC-AL	25 mL	sterile, store below 4° C, Do not Freeze
Calcium chloride solution	CSR-ABC-CA	100 mL	sterile, store below 4° C, Do not Freeze
Sodium citrate solution	CSR-ABC-CI	100 mL	sterile, store below 4° C, Do not Freeze

For research use only. Not for diagnostic use.



COSMO BIO CO., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU, TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobiousa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobiousa.com



Application example of Alginate 3D Cell Culture Kit

1. Chondrocyte

Normal porcine chondrocytes prepared from articular cartilage of knee joint were cultured in alginate beads at 2×10^6 cells/mL (in DMEM/F-12 containing 10%FBS, 100ng/mL IGF-I and 25 μ g/mL L-ascorbic acid) in the presence or absence of oversulfated chondroitin sulfate. At end of culture, the alginate beads were dissolved with sodium citrate solution and the cell pellet was digested with Pronase solution (1mg/mL). The DNA and collagen contents of the digested sample were determined using the Hoechst 33258 fluorescent dye method and measurement of hydroxyproline content, respectively (Tissue Eng Part A. 2010 May; 16:1575-84).




Image 1. Normal porcine chondrocytes cultured in alginate beads.




Figure 1. The effect of oversulfated chondroitin sulfate on the collagen production of chondrocytes in alginate beads. (at Day 6, Mean $\pm$ SD, n=3)

2. HepG2 cells

A human liver carcinoma cell line, HepG2 cells, was cultured in alginate beads (5×10^5 cells/mL, 10 beads/well, DMEM containing 10% FBS) for 9 days.




Image 2. HepG2 cells cultured in alginate beads.

HepG2 cells were cultured in alginate beads (1×10^5 cells/mL, 10 beads/well, DMEM containing 10% FBS). At the end of the culture, the alginate beads were dissolved with sodium citrate solution and the cell pellet was digested with Pronase solution (1 mg/mL). The DNA content of the digested sample was determined using the Hoechst 33258 fluorescent dye method.




Figure 2. Growth curve of HepG2 cells cultured in alginate beads (Mean $\pm$ SD, n=3)

HepG2 cells were suspended in alginate beads and cultured (2×10^5 cells/mL, 10 beads/well, DMEM containing 10% FBS). After 24 hours, antitumor drug, 5-fluorouracil (5-FU), was added and further cultured for 6 days. At end of culture, the alginate beads were dissolved with sodium citrate solution and the cell pellet was digested with Pronase solution (1 mg/mL). The DNA content of the digested sample was determined using the Hoechst 33258 fluorescent dye method.




Figure 3. Effect of 5-FU on the growth of HepG2 cells in alginate beads.

3. Saos-2 cells

A human osteosarcoma cell line, Saos-2 cells, was cultured in alginate beads (2×10^5 cells/mL, 10 beads/well, DMEM/F-12 containing 10% FBS). At end of the culture, the alginate beads were dissolved with sodium citrate solution and the cell pellet was digested with Pronase solution (1 mg/mL). The DNA content of the digested sample was determined using the Hoechst 33258 fluorescent dye method.




Image 3. Saos-2 cells cultured in alginate beads for 14 days.




Figure 4. Growth of Saos-2 cells cultured in alginate beads (Mean±SD, n=3)

4. 3T3-L1 cells

A murine fibroblast cell line, 3T3-L1 cells, was cultured in alginate beads (2×10^5 cells/mL, 10 beads/well, DMEM containing 10% FBS). At end of the culture, the alginate beads were dissolved with sodium citrate solution and the cell pellet was digested with Pronase solution (1 mg/mL). The DNA content of the digested sample was determined using the Hoechst 33258 fluorescent dye method.




Figure 5. Growth of 3T3-L1 cells cultured in alginate beads (Mean±SD, n=3)
(3T3-L1 cells cannot growth in alginate beads)



COSMO BIO CO., LTD.
Inspiration for Life Science

Certification of Analysis

Product Name: Sodium alginate solution

Code Number: CSR-ABC-AL

Lot Number: 15A279

Size: 25mL


Shipping Date : 2015.04.22

Exp. Date : 2017.3.31

Remark: Specifications are to be modified without notice.

	Result	Specification	Method
Colony Formation	pass	Recognition of colony formation by HepG2 cells in alginate beads	microscopy
Solubility of Alginate Beads	pass	Solubilized with sodium citrate solution within 15min	visual observation

Signature :


Tatsuya Miyazaki



COSMO BIO CO., LTD.

【JAPAN】

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU, TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>



COSMO BIO USA

【Outside Japan】

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
email: info@cosmobiousa.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobiousa.com

Manufacturer

PG
Research

PG Research, 4-13-11, Hanakoganei,
Kodaira, Tokyo 187-0002, Japan
URL: <http://www.pg-r.com>