# Spin column based Antibody Purification Kit

Cat. No. CSR-APK-10A

CSR-APK-10G

Last Updated: 2020/02/28

[www.cosmobiousa.com](http://www.cosmobiousa.com)

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This spin-column type IgG antibody-purify kit consists of immobilized Protein A (#APK-10A) or Protein G (#APK-10G) on the solid phase of filter type monolithic silica with uniform continuous pores. Compared with conventional ones, this product can easily purify IgG antibodies for a shorter time (approximately 3 minutes). The column can be recycled up to 5 times<sup>\*1</sup> without diminishing of antibody binding capacity. Each can purify IgG antibodies up to 0.4 mg with Protein A and 0.3 mg with Protein G per column under more than 90% of recovery rate.

\*1: In case purified antibodies or culture supernatants are used, up to five replicates of the columns can be reused. However, if the sample contains a solid component such as protein in the serum, it may cause a clogging of the column, which may result in a decrease in the regeneration frequency and fluid permeability. Therefore, this product should be used after removal of the solid component by centrifugal removal or filter filtration. Also, the number of replays will vary depending on the sample type of the customer, so it is not responsible for our web.


## 【II】 Composition

	Quantity
Monolithic silica-based Spin Column	10 columns
Sampling tube	10 tubes
Binding Buffer	30 mL
Wash Buffer	30 mL
Elution Buffer	30 mL
Neutralization Buffer	15 mL
Regeneration Buffer	30 mL

### 【III】 Cautions

1. Use the spin column with the cover opened.  
The sample will not pass through the monolithic silica smoothly when the cover of the spin column is closed as there will be a negative pressure generated in the spin column. Make sure the cover of the spin column is opened.
2. The spin column is filled with 20% ethanol solution. When reusing, be careful not to dry the column by filling it with 20% ethanol and store.  
Drying may degrade column performance.
3. Do not drop or bump the spin columns. Subjecting the spin columns to shocks may cause the monolithic silica seriously damaged.

### 【IV】 Typical Procedure

1. Adjust the pH of the antibody sample to around 7.0 in advance.  
If the sample amount is small, adjust to 500 µL with PBS etc.  
\*Before use, centrifuge the ethanol solution at 6,000 rpm for 30 seconds and discard.
2. Add 500 µL of Binding Buffer to the column, centrifuge at 6,000 rpm for 30 seconds and remove it.
3. Add 500 µL of antibody sample to the column, centrifuge at 6,000 rpm for 30 seconds, discard flow-through (FT).<sup>\*2</sup>

\*2: When concentrating the antibody, centrifuge and remove FT each time, then add 500 µL of antibody sample to the column. Sequential runs may cause clogging, so wash the column with Wash Buffer every 2 to 3 addition of the sample.

4. Add 500 µL of Wash Buffer to the column, centrifuge at 6,000 rpm for 30 seconds, discard FT.
5. Add 50 µL Neutralization Buffer to the sampling tube and set the column. Add 500 µL<sup>\*3</sup> of Elution Buffer to the column and centrifuge at 6,000 rpm for 30 seconds to recover the antibody sample.

\*3: Elution Buffer can be reduced to 100 µL and eluted at high concentration.

In this case, add Neutralization Buffer to the sampling tube at the ratio; Neutralization Buffer : Elution Buffer = 1: 10.

However, it may not be able to elute all the adsorbed IgG antibodies.

6. To regenerate the used spin column, add 500 µL of Regeneration buffer to the used spin column and centrifuge for 30 seconds, and remove FT. Store the spin column with 20% ethanol and store it at 4°. Do not allow the resin (silica monolith) to dry. Otherwise, the performance and quality of spin column may be deteriorated.

Column type	Recommended sample	
Protein G	Monoclonal antibody	Human IgG1, IgG2, IgG3, IgG4
		Mouse IgG1, IgG2a, IgG2b, IgG3
	Polyclonal antibody	Rabbit, Goat
Protein A	Monoclonal antibody	Human IgG1, IgG2, IgG4
		Mouse IgG2a, IgG2b, IgG3
	Polyclonal antibody	Rabbit

COSMO BIO Co., LTD.

【JAPAN】  
TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>

COSMO BIO USA<sup>12077</sup>

【Outside Japan】  
2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobioua.com](mailto:info@cosmobioua.com)  
URL: [www.cosmobioua.com](http://www.cosmobioua.com)  
Phone/FAX: (+1) 760-431-4600

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# Spin column based Antibody Purification Kit

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内容	容量
Monolithic silica-based Spin Column	10 columns
Sampling tube	10 tubes
Binding Buffer	30mL
Wash Buffer	30mL
Elution Buffer	30mL
Neutralization Buffer	15mL
Regeneration Buffer	30mL

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