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- 1) 曾根勝, 知久幹夫, 吉田光敏, 番場公雄, 小笠晃. 各種希釈保存液を用いた豚液状精液の長期保存試験.
日豚会誌 1992; 29: 41-50.
- 2) Mito T, Yoshioka K, Nagano M, Suzuki C, Yamashita S, Hoshi H. Transforming growth factor- α in a defined medium during in vitro maturation of porcine oocytes improves their developmental competence and intracellular ultrastructure. Theriogenology 2009; 72: 841-850.
- 3) Suzuki C, Yoshioka K. Effects of amino acid supplements and replacement of polyvinyl alcohol with bovine serum albumin in porcine zygote medium. Reprod Fertil Dev 2006; 18: 789-795.
- 4) Suzuki C, Yoshioka K, Sakatani M, Takahashi M. Glutamine and hypotaurine improves intracellular oxidative status and in vitro development of porcine preimplantation embryos. Zygote 2007; 15: 317-324.
- 5) Yoshioka K, Suzuki C, Tanaka A, Anas IMK, Iwamura S. Birth of piglets derived from porcine zygotes cultured in a chemically defined medium. Biol Reprod 2002; 60: 112-119.
- 6) Yoshioka K, Suzuki C, Itoh S, Kikuchi K, Iwamura S, Rodriguez-Martinez H. Production of piglets derived from in vitro-produced blastocysts fertilized and cultured in chemically defined media: effects of theophylline, adenosine, and cysteine during in vitro fertilization. Biol Reprod 2003; 69: 2092-2099.
- 7) Yoshioka K, Suzuki C, Onishi A. Defined system for in vitro production of porcine embryos using a single basic medium. J Reprod Dev 2008; 54: 208-213.
- 8) Mito T, Yoshioka K, Noguchi M, Yamashita S, Hoshi H. Recombinant human follicle-stimulating hormones and transforming growth factor- α enhance in vitro maturation of porcine oocytes. Mol Reprod Dev 2013; 80: 549-560.
- 9) Mito T, Yoshioka K, Yamashita S, Suzuki C, Noguchi M, Hoshi H. Glucose and glycine synergistically enhance the in vitro development of porcine blastocysts in a chemically defined medium. Reprod Fertil Dev 2012; 24: 443-450.

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TEL: (03) 5632-9610 FAX: (03) 5632-9619

TEL: (03) 5632-9620

Complete Culture Kit for Porcine Embryos in vitro

Catalog No. : CK030

For research use only, Not for diagnostic use.

- Please read this manual thoroughly before use -

www.cosmobio.co.jp

This kit is a complete culture kit for porcine embryos in vitro containing porcine oocyte/embryo collection medium (POE-CM), basic medium for porcine oocyte maturation (POM), porcine fertilization medium (PFM), defined medium for porcine embryos (PZM-5), defined medium for late stage porcine embryos (PBM), dbcAMP concentrated solution-100X (dbcAMP-100X) and Replate.

This kit can efficiently produce transferable embryos (blastocysts) under a low oxygen condition (5%O₂/5%CO₂/90%N₂, 39°C) without the coculture of cumulus/granulosa cells. This kit can support to increase the blastomere numbers of blastocysts and the formation of hatched blastocysts.

《I. Kit components》

- 1) Porcine oocyte/embryo collection medium : POE-CM(Cat.# CK020), 100 ml x 3
- 2) Basic medium for porcine oocyte maturation : POM(Cat.# CK021), 25 ml x 2
- 3) Porcine fertilization medium : PFM(Cat.#CK023), 100 ml
- 4) Defined medium for porcine embryos : PZM-5(Cat.# CK024), 25 ml
- 5) Defined medium for late stage porcine embryos : PBM(Cat.# CK025), 25 ml
- 6) dbcAMP concentrated solution-100X : dbcAMP-100X(Cat.# CK027), 0.3 ml
- 7) Replate : (Cat.# CK028), 10 plates x 3

〈Necessary Equipments and Reagents〉

- 1) Sterile PBS (-) or saline solution (0.9%)
- 2) Equine chorionic gonadotropin (eCG; 1,000IU): Yell pharmaceutical co., Ltd.
- 3) Human chorionic gonadotropin (hCG; 5,000IU): Yell pharmaceutical co., Ltd.
- 4) Paraffin liquid reagent(Mineral oil): Merk chemicals (Cat. No.: 107162)
- 5) Percoll: GE Healthcare UK Ltd. (17-0891-02)
- 6) Sperm diluent for in vitro fertilization (SEM-5X): COSMO BIO Co., Ltd (Cat. No.: CK026)
- 7) Culture dish
 - 35 mm culture dish: Greiner bio-one (627160)
 - Petri Dish: FALCON (1112)
- 8) Centrifuge tube (15 ml, 50 ml)
- 9) Disposable syringe (10 ml)
- 10) Syringe needle (18 G)
- 11) Pipette
- 12) Glass capillary pipette
- 13) Aspiration tube
- 14) Straw cutter
- 15) Hemacytometers, Thoma
- 16) Stereomicroscope
- 17) Multi-gas incubator
- 18) Clean bench
- 19) Centrifuge
- 20) Thermostat bath
- 21) Micro thermo plate
- 22) Vortex mixer
- 23) Others

〈Storage〉

Stored at 4°C in a dark room. Do not freeze.

〈Expiration date〉

The expiration date is shown on the label of the box (4 months after the preparation).

www.cosmobio.co.jp

《 II . Reagent preparation 》

☆ Preparation of Gonadotropin stock solutions (100X) (eCG, hCG each 1,000 IU/ml)

1) 2,000 IU/ml eCG solution

Add and dissolve eCG powder (1,000 IU/tube) in 0.5 ml POM medium.

2) 2,000 IU/ml hCG solution

Add and dissolve hCG powder (5,000 IU/tube) in 2.5 ml POM medium.

3) eCG and hCG (each 1,000 IU/ml) mixed solution

Mix equivalent of 2,000 IU/ml eCG and 2,000 IU/ml hCG and then, dispense 70 ul each and store at - 30°C.

☆ Sperm diluent for in vitro fertilization (SEM-5X*)

Composition	SEM-5X
Glucose	137.50 g/L
Na-citrate · 2H ₂ O	34.50 g/L
Citric acid · H ₂ O	14.50 g/L
NaHCO ₃	5.00 g/L
EDTA · 2Na · 2H ₂ O	11.75 g/L
Tris (TRIZMA BASE)	28.25 g/L
Gentamicin	0.05 g/L

Dispense and store at -30°C after filter sterilization with 0.22 um filter.

*“Sperm diluent for in vitro fertilization (SEM-5X) (Cat.No. : CK026)” is ready to use.

☆ Preparation of isotonic Percoll solution

Composition	80% Percoll solution	50% Percoll solution
Percoll	1.6 ml	1 ml
MilliQ water	-	0.6 ml
SEM-5×	0.4 ml	0.4 ml

Prepare just before use.

《 III. Medium preparation 》

☆ Preparation of basic medium for porcine oocyte maturation (POM)

- 1st POM medium : POM + 1 mM dbcAMP + 10 IU/ml eCG + 10 IU/ml hCG

1st POM medium (Prepare just before use)

- 1) Add 50 ul dbcAMP-100X to 4.95 ml POM medium (1 mM dbcAMP solution, Total 5 ml).
- 2) Remove 50 ul of 1 mM dbcAMP solution (Total 4.95 ml)
- 3) Add 50 ul Gonadotropin stock solution to 4.95 ml, 1 mM dbcAMP solution (1st POM medium, Total 5 ml).

- ① Oocyte washing : 1st POM (200 ul) + Mineral oil (about 100 ul)
- ② Oocyte maturation : 1st POM (150 ul) + Mineral oil (about 150 ul)

• 2nd POM medium : POM medium (Basic medium)

- ③ Oocyte washing : POM (200 ul) + Mineral oil (about 100 ul)
- ④ Oocyte maturation : POM (150 ul) + Mineral oil (about 150 ul)

☆ Preparation of porcine fertilization medium (PFM)

- ⑤ Oocyte washing : PFM (200 ul) + Mineral oil (about 100 ul)
- ⑥ In vitro fertilization : PFM (50 ul) + Mineral oil (about 200 ul)

☆ Preparation of defined medium for porcine embryos (PZM-5)

- ⑦ Embryo washing : PZM-5 (200 ul) + Mineral oil (about 100 ul)
- ⑧ Embryo culture : PZM-5 (50 ul) + Mineral oil (about 200 ul)

☆ Preparation of defined medium for late stage porcine embryos (PBM)


- ⑨ Embryo washing : PBM (200 ul) + Mineral oil (about 100 ul)
- ⑩ Late stage embryo culture : PBM (50 ul) + Mineral oil (about 200 ul)

* Prepare the given number of medium drops in the Replate and set the gas-phase equilibrium in a incubator (39°C, 5%CO₂,5%O₂,90%N₂) before the culture.

《IV. In vitro production of porcine embryos》

1. Oocyte collection and oocyte maturation culture

- 1) Wash porcine ovaries several times with PBS (-) or saline at room temperature (about 25°C) and collect cumulus-oocyte complexes (COCs) by aspirating from ovarian follicles (diameter=3-6 mm).
- 2) Transfer the follicular fluid containing COCs to centrifuge tube and stand it for 20 min.
- 3) Remove the supernatant and dilute the precipitation with POE-CM medium at room temperature. Dispense the precipitate to Petri dish. Harvest and select the oocytes that the cumulus cells adhere to more than two or three layers (cumulus-oocyte complexes;COCs) under the stereomicroscope (Picture 1).
- 4) Wash the selected COCs twice with 1st POM medium (①) prepared in Reproplate. Put 15-20 COCs in one drop of 1st POM medium (②) and then, incubate for 20-22 hours.
- 5) Wash the COCs three times with 2nd POM medium (③) prepared in Reproplate after 20-22 hours of the maturation culture. Transfer the COCs into the drop of 2nd POM medium (④) and incubate for additional 20-24 hours.

Picture 1 : Cumulus-oocyte complexes (COCs) just after the collection

2. Sperm preparation and in vitro fertilization


- 1) Keep warm PFM medium (about 20 ml) and 80% and 50% isotonic Percoll solution in a thermostat bath at 38°C.
- 2) Thaw frozen sperm*1 in the straw with warm water at 38°C. Suspend them into 50% Percoll solution.
- 3) Overlay the sperm suspension gently on the top of 80% Percoll solution. Then, centrifuge at 700 x g for 20 min (35°C or room temperature).
- 4) Aspirate the supernatant and add about 7 ml PFM medium to it and gently mix. Centrifuge it at 500 x g for 5 min.
- 5) Repeat 4).
- 6) Aspirate the supernatant and add 0.5-1 ml PFM medium to it and mix. Count the numbers of sperm with a hemacytometer.
- 7) Dilute sperm suspension with PFM medium at the concentrations of $2-20 \times 10^6/\text{ml}^{*2}$.
- 8) Prepare the drops of gas-phase equilibrated PFM medium (⑥) in the Reproplate and add 50 ul of sperm suspension (final sperm concentrations: $1-10 \times 10^6/\text{ml}^{*2}$) in each drop.
- 9) Wash the COCs twice with PFM medium (⑤) in the Reproplate after the in vitro maturation. And then, transfer 15-20 COCs into the drop containing sperm suspension and incubate for 10~20 hours*2.

*1 Please contact us if you would like to use fresh spermatozoa.


*2 The optimal sperm concentration and the incubation time of insemination are variable for different sperm preparation. Please make sure to determine these appropriate conditions before the experiments.

3. Embryo culture in vitro

- 1) Wash fertilized oocytes once with POE-CM medium at 38°C after in vitro fertilization. Transfer them into 1ml POE-CM medium of 15 ml centrifuge tube and remove cumulus cells of them by mixing with vortex mixer for 4 min.
- 2) Collect the denuded embryos (Picture 2) from centrifuge tube. Remove the fertilized oocytes with strongly attached cumulus cells and/or with degenerated cytoplasm.
- 3) Wash the selected fertilized oocytes twice with POE-CM medium. Wash them twice with PZM-5 medium (⑦) in the Reproplate.
- 4) Transfer 20-30 embryos into each drop of PZM-5 medium(⑧) in the Reproplate and incubate them. They develop to morula and blastocyst after 4-5 days after in vitro fertilization (Picture 3).
- 5) Wash morula and blastocyst twice with PBM medium (⑨) in the Reproplate. Transfer 10-20 of those into each drop of PBM medium (⑩) and incubate them. They develop to expanded blastocysts and hatched blastocysts after an additional 1-2 days of culture.


Picture 2 : Denuded embryos

Picture 3 : Blastocysts

Process of porcine embryo culture in vitro

< Reference >

- 1) Mito T, Yoshioka K, Nagano M, Suzuki C, Yamashita S, Hoshi H. Transforming growth factor- α in a defined medium during in vitro maturation of porcine oocytes improves their developmental competence and intracellular ultrastructure. *Theriogenology* 2009; 72: 841-850.
- 2) Suzuki C, Yoshioka K. Effects of amino acid supplements and replacement of polyvinyl alcohol with bovine serum albumin in porcine zygote medium. *Reprod Fertil Dev* 2006; 18: 789-795.
- 3) Suzuki C, Yoshioka K, Sakatani M, Takahashi M. Glutamine and hypotaurine improves intracellular oxidative status and in vitro development of porcine preimplantation embryos. *Zygote* 2007; 15: 317-324.
- 4) Yoshioka K, Suzuki C, Tanaka A, Anas IMK, Iwamura S. Birth of piglets derived from porcine zygotes cultured in a chemically defined medium. *Biol Reprod* 2002; 60: 112-119.
- 5) Yoshioka K, Suzuki C, Itoh S, Kikuchi K, Iwamura S, Rodriguez-Martinez H. Production of piglets derived from in vitro-produced blastocysts fertilized and cultured in chemically defined media: effects of theophylline, adenosine, and cysteine during in vitro fertilization. *Biol Reprod* 2003; 69: 2092-2099.
- 6) Yoshioka K, Suzuki C, Onishi A. Defined system for in vitro production of porcine embryos using a single basic medium. *J Reprod Dev* 2008; 54: 208-213.
- 7) Mito T, Yoshioka K, Noguchi M, Yamashita S, Hoshi H. Recombinant human follicle-stimulating hormones and transforming growth factor- α enhance in vitro maturation of porcine oocytes. *Mol Reprod Dev* 2013; 80: 549-560.
- 8) Mito T, Yoshioka K, Yamashita S, Suzuki C, Noguchi M, Hoshi H. Glucose and glycine synergistically enhance the in vitro development of porcine blastocysts in a chemically defined medium. *Reprod Fertil Dev* 2012; 24: 443-450.

For research use only, Not for diagnostic use.

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COSMO BIO Co., LTD.