Fluoresceinamine-Labeled Sodium Chondroitin Sulfate E (Squid Cartilage)

Product Code (Outside Japan): CSR-FACS-E2(SqC)1, CSR-FACS-E2(SqC)3

Product Code (Japan): FACS-E2(SqC)1, FACS-E2(SqC)3

Labeled Content: 1 mg/vial (FACS-E2(SqC)1), 3 mg/vial (FACS-E2(SqC)3)

Appearance: yellow—orange amorphous block or powder (lyophilizate)

Source of Sodium Chondroitin Sulfate E: Squid Cartilage

Fluorescent Probe: Fluoresceinamine

CAS number of Fluorescent Probe: 3326-34-9

Overview: Chondroitin sulfate (CS) is a sulfated glycosaminoglycan composed of repeating disaccharide units of N-acetyl-D-galactosamine (GalNAc) and D-glucuronic acid (GlcUA). CS is abundant in cartilage, where it exists as unbranched polysaccharide chains covalently linked to the core proteins of proteoglycans.

The Fluoresceinamine (FA)-Labeled Sodium Chondroitin Sulfate E (FACS-E) in this product is prepared by forming an amide bond between the fluorescent moiety FA and the carboxyl group of GlcUA within CS-E¹⁾. The raw material CS-E contains approximately 50% E-type structural disaccharide units (4 and 6-O-sulfations of GalNAc). After conjugation, FA exhibits strong fluorescence, whereas unconjugated FA shows only weak fluorescence (approximately 1/150 of that after conjugation)¹⁾, and therefore free FA has no significant effect on measurements.

Each vial contains 1 mg (FACS-E2(SqC)1) or 3 mg (FACS-E2(SqC)3) of lyophilized FACS-E. The excitation wavelength is 485—500 nm, and the emission wavelength is 510—525 nm. The attached Certificate of Analysis provides the actual measured values for each product specification.

Handling precautions:

- 1) The interior of the vial is depressurized. Please take care when opening the vial.
- 2) This product is hygroscopic. Store at -20 °C or below, protected from moisture.
- 3) After reconstitution, store the solution at -20 °C or below, protected from light. Preparing aliquots suitable for anticipated use is recommended.
- 4) Protect the product from strong light during handling. Exposure to intense light may cause photobleaching of FADS. The product may be used at room temperature if protected from strong light.
- 5) FACS-E is degraded by oxidising agents, reducing agents, strong acids, strong alkalis and heat.
- 6) **Fluorescence intensity varies depending on the pH of the solution and decreases under acidic conditions.** Please take note of the sample solution's pH when measuring fluorescence intensity.
- 7) This product is for research use only. It is not intended for use in diagnostic or therapeutic procedures.

Reference:

- 1) Ogamo A et al.: Carbohydr. Res., **105**, 69 (1982)

NOTICE: For R&D use only. Do not use for drug, household, cosmetically and others.

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