Fluoresceinamine-Labeled Sodium Hyaluronate (H1)

(The average molecular weight: $1200 \times 10^3 - 1600 \times 10^3$)

Product Code (Outside Japan): CSR-FAHA-H1

Product Code (Japan): FAHA-H1

Labeled Content: 3 mL/vial (1 mg/mL, in phosphate buffered saline (PBS) (–))

Appearance: yellow green solution

CAS number of Sodium Hyaluronate: 9067-32-7

Fluorescent Probe: Fluoresceinamine

CAS number of Fluorescent Probe: 3326-34-9

Overview: Hyaluronan (HA) is a glycosaminoglycan composed of repeating disaccharide units of N-acetyl-D-glucosamine (GlcNAc) and D-glucuronic acid (GlcUA). HA is abundant in synovial fluid, skin, the umbilical cord, and the vitreous body, where it exists as unbranched polysaccharide chains.

The Fluoresceinamine (FA)-Labeled Sodium Hyaluronate (FAHA) in this product is prepared by forming an amide bond between the fluorescent moiety FA and the carboxyl group of GlcUA within HA¹⁾. After conjugation, FA exhibits strong fluorescence, whereas unconjugated FA shows only weak fluorescence (approximately 1/150 of that after conjugation)¹⁾, and therefore free FA has no significant effect on measurements.

This product is sterilized by filtration and supplied as 3 mL per vial of a 1 mg/mL FAHA solution in PBS (–). The endotoxin level meets the product specifications. The excitation wavelength is 485–500 nm, and the emission wavelength is 510–525 nm. The average molecular weight is $1200 \times 10^3 - 1600 \times 10^3$. The attached Certificate of Analysis provides the actual measured values for each product specification.

Handling Precautions:

- 1) Protect the product from strong light as much as possible during handling. Exposure to strong light causes photobleaching of FAHA and degradation into lower molecular weight forms. The product may be used at room temperature if protected from strong light.
- 2) FAHA is degraded by oxidizing agents, reducing agents, strong acids, strong alkalis, and heat.
- 3) After thawing, gently and thoroughly mix the vial before use.
- 4) Store at $-20\text{ }^{\circ}\text{C}$ or below, protected from light. Preparing aliquots suitable for anticipated use is recommended.
- 5) **Fluorescence intensity varies depending on the pH of the solution and decreases under acidic conditions.** Please take note of the sample solution's pH when measuring fluorescence intensity.
- 6) This product is for research use only. It is not intended for use in diagnostic or therapeutic procedures.

Reference:

- 1) Ogamo A et al.: Carbohydr. Res., **105**, 69 (1982)

NOTICE: For R&D use only. Do not use for drug, household, cosmetically and others.

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