Fluoresceinamine-Labeled Sodium Heparan Sulfate (Porcine Kidney)

Product Code (Outside Japan): CSR-FAHS-P2(PgK)1, CSR-FAHS-P2(PgK)3

Product Code (Japan): FAHS-P2(PgK)1, FAHS-P2(PgK)3

Labeled Content: 1 mg/vial (FAHS-P2(PgK)1), 3 mg/vial (FAHS-P2(PgK)3)

Appearance: yellow—orange amorphous block or powder (lyophilizate)

Source of Sodium Heparan Sulfate: Porcine Kidney

Fluorescent probe: Fluoresceinamine

CAS number of fluorescent probe: 3326-34-9

Overview: Heparan sulfate (HS) is a sulfated glycosaminoglycan composed of repeating disaccharide units of N-acetyl-D-glucosamine (GlcNAc) or N-sulfo-D-glucosamine (GlcNS) and D-glucuronic acid (GlcIA) or L-iduronic acid (IdoUA). HS is abundant in the basement membrane, where it exists as unbranched polysaccharide chains covalently linked to the core proteins of proteoglycans. Tissues rich in the basement membrane, such as the kidney or blood vessels, are commonly used as raw materials for HS production.

The Fluoresceinamine (FA)-Labeled Sodium Heparan Sulfate (FAHS) in this product is prepared by forming an amide bond between the fluorescent moiety FA and the carboxyl group of GlcUA or IdoUA within HS¹⁾. After conjugation, FA exhibits strong fluorescence, whereas unconjugated FA shows only weak fluorescence (approximately 1/150 of that after conjugation)¹⁾, and therefore free FA has no significant effect on measurements.

Each vial contains 1 mg (FAHS-P2(PgK)1) or 3 mg (FAHS-P2(PgK)3) of lyophilized FAHS. The excitation wavelength is 485—500 nm, and the emission wavelength is 510—525 nm. The attached Certificate of Analysis provides the actual measured values for each product specification.

Handling precautions:

- 1) This product is hygroscopic. Store at -20 °C or below, protected from moisture.
- 2) After reconstitution, store the solution at -20 °C or below, protected from light. Preparing aliquots suitable for anticipated use is recommended.
- 3) Protect the product from strong light during handling. Exposure to intense light may cause photobleaching of FAHS. The product may be used at room temperature if protected from strong light.
- 4) FAHS is degraded by oxidizing agents, reducing agents, strong acids, strong alkalis and heat.
- 5) **Fluorescence intensity varies depending on the pH of the solution and decreases under acidic conditions.** Please take note of the sample solution's pH when measuring fluorescence intensity.
- 6) This product is for research use only. It is not intended for use in diagnostic or therapeutic procedures.

Reference :

- 1) Ogamo A et al.: Carbohydr. Res., **105**, 69 (1982)

NOTICE: For R&D use only. Do not use for drug, household, cosmetically and others.

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