ExoTrap™ Exosome Isolation Spin Column Kit for protein research

Cat. No. CSR-SHI-EXO-K010

Updated on Aug 8th, 2016

www.cosmobio.com

【1】 Background

Exosomes are 40 –100 nm membrane vesicles. They are found in many biological fluids including blood, urine, saliva, amniotic fluid, cell culture supernatants, and malignant ascites fluids. CD9 is a well known surface marker for exosomes. ExoTrap Spin Columns utilize immobilized CD9 antibody to capture exosomes from biological fluids in as little as 30 minutes. The contents of the captured exosomes are then readily processed for downstream analysis by western blotting, mass spectrometry, PCR, or sequencing, etc.

【2】 Component

Component	Quantity	Storage
ExoTrap™ Spin Column	10	4°C

Required but not provided:

- 1.7 mL tube
- PBS
- Elution Buffer

[3] Protocol for samples with low protein concentration like culture supernatant, urine, cerebrospinal fluid etc...

Note: Sample should be prepared at the discretion of the user. The following recommendations are only guidelines and may be altered to optimize or complement **the user's experimental design.**

Prepare Sample

1. Transfer 10 mL of culture supernatant to 15 mL tube.
2. Centrifuge the tube at 3,000 rpm for 15 minutes at 4°C .
3. Filter the supernatant by syringe and 0.22 µm disc filter and transfer to new tube.

Prepare Column

4. Open the top and bottom caps of ExoTrap™, set the column to 1.7 mL tube and centrifuge at 5,000 rpm for 1 minute at 20°C . Remove the medium.

Load Column with Sample

5. Apply 600 µL of sample to the column and centrifuge at 5,000 rpm for 1 minute at 20°C . Remove the flow through.
6. Repeat procedure 5 for 10 times. (6 mL of supernatant was used in total.)

Wash Sample-Loaded Column

7. Apply 600 µL PBS to the column and centrifuge at 5,000 rpm for 1 minute at 20°C and remove the flow through.
8. Repeat procedure 7 once again and wash the column.

Elute Exosomes

9. Transfer the column to new 1.7 mL tube.
10. Apply 50 µL of the desired elution buffer to the column filter (top of column).

For Western Blotting Analysis: use SDS sample buffer

For Mass spectroscopic Analysis: use a solution of 8M Urea, 50 mM ammonium bicarbonate


For RNA recover and miRNA Analysis: use a denaturing solution

11. Incubate at 37°C for 30 minutes.
12. Centrifuge at 7,500 rpm for 1 minute at 20°C and collect sample.

Isolation of exosomes from cell culture supernatant

Using ExoTrap™, exosome-derived protein was isolated by 50 µL SDS sample buffer and did western blotting analysis.

Sample: 293T cell culture supernatant, MKN7 cell culture supernatant
Primary Antibody: Anti CD9 (#CSR-SHI-EXO-M01)
(20 µL / lane was applied.)

【4】 Protocol for blood serum or blood plasma

Note: Sample should be prepared at the discretion of the user. The following recommendations are only guidelines and may be altered to optimize or complement the user's experimental design.

Prepare Sample

1. Dilute blood serum or blood plasma 6 times with PBS.
2. Centrifuge at 12,000 rpm for 15 minutes at 4°C .
3. Filter the supernatant by syringe and 0.22 μm disc filter or 0.22 μm spin filter tube and transfer to new tube.

Prepare Column

4. Open the top and bottom caps of ExoTrap™, set the column to 1.7 mL tube and centrifuge at 5,000 rpm for 1 minute at 20°C . Remove the medium.

Load Column with Sample


5. Apply 600 μL of sample to the column and centrifuge at 5,000 rpm for 1 minute at 20°C . Remove the flow through.
6. Repeat procedure 5 as needed. (In the example below, 200 μL versenized Human plasma was diluted 6 times and centrifuged twice.)

Wash Sample-Loaded Column

7. Apply 600 μL PBS to the column and centrifuge at 5,000 rpm for 1 minute at 20°C and remove the flow through.
8. Repeat procedure 7 once again and wash the column.

Elute Exosomes

9. Transfer the column to new 1.7 mL tube.
10. Apply 50 μL of the desired elution buffer to the column filter (top of column).
 - For Western Blotting Analysis: use SDS sample buffer
 - For Mass spectroscopic Analysis: use a solution of 8M Urea, 50 mM ammonium bicarbonate
 - For RNA recover and miRNA Analysis: use a denaturing solution
11. Incubate at 37°C for 30 minutes.
12. Centrifuge at 7,500 rpm for 1 minute at 20°C and collect sample.

Sample: Versenized Human plasma
primary antibody: Anti CD9 (#CSR-SHI-EXO-M01)
(20 μL / lane was applied.)

Isolation of exosomes from human plasma

Using ExoTrap™, exosome-derived protein was isolated by 50 μL SDS sample buffer and did western blotting analysis.

COSMO BIO Co., LTD.

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
TEL : (81)3-5632-9617
FAX : (81)3-5632-9618
e-mail : export@cosmobio.co.jp
URL : www.cosmobio.com

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ExoTrap™ Exosome Isolation Spin Column Kit for protein research

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