For Research use only. Not for clinical diagnosis



## Mount-Quick / Mount-Quick "Aqueous"

## Catalog No. DAI-DM-01 / DAI-DM-02

Updated on Oct 6th, 2020





Mount Quick

Mount-Quick "Aqueous"

Purpose of use	Sample mounting medium / mounter for permanent preservation
Features	Mount Quick
	ullet Rapid drying, highly transparent, and does not form bubbles. Highly soluble in xylene,
	does not degenerate or stain the specimen. Can withstand long-term storage after
	mounting.
	$\cdot$ Tube design provides easy control of amounts used and limits air exposure, preventing
	variability in concentration.
	Mount-Quick "Aqueous"
	<ul> <li>Maintains superior transparency when fully hardened.</li> </ul>
	Low fluorescence.
	<ul> <li>Mounting is simple with an easy-to-use dropper bottle.</li> </ul>
How to use	Mount Quick
	Squeeze out an amount approximately the same size as the sample section at one
	location on the cover glass or sample.
	Mount-Quick "Aqueous"
	Drop onto the cover glass or section from above.
	(Prevent bubble formation by inverting the container and touching the tip on the bottom
	of the drop.)
	* Before first use, cut open the container tip at a 45° angle with a box cutter or scissors.
Precautions	After use, always re-cap and store in a cool, dark place.



## **Application Data**

Use Mount Quick Aqueous when samples are not dehydrated and cleared (photograph as the color fades during long-term storage) and use Mount Quick Tube when samples are dehydrated and cleared (can be permanently preserved).

## Staining data chart for Mount Quick Tube and Mount Quick Aqueous

DM01 Mount Quick Tube (non-water soluble mounting medium for mounting after dehydration and clearing)		
Hematoxylin-eosin staining		
Azan stain		
Grimelius stain		
Silver stain		
Van Gieson's stain		
Periodic acid-Schiff (PAS) stain		
Giemsa stain		
Immunostaining (non-aqueous DAB)		
DM02 Mount Quick Aqueous (water soluble mounting medium for mounting after staining without dehydration)		
Lipid stains		
Immunohistochemical staining (AEC)		
Blood smear staining		
PAS, FITC, DAB staining		

Both Mount Quick Tube and Mount Quick Aqueous can be used with methyl green and BCIP/NBT staining, but are not suitable for long-term storage. After staining, we recommend photographing as soon as possible, as the color can fade in 24 hours after drying, dehydrating, clearing, and mounting. BCIP/NBT can stain red and black with immunostaining.



Cosmo Bio Co., Ltd.

[JAPAN] TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU. TOKYO 135-0016, JAPAN Phone: +81-3-5632-9610 FAX: +81-3-5632-9619 URL: https://www.cosmobio.co.jp/



[Outside Japan] 2792 Loker Ave West, Suite 101 Carlsbad, CA 92010, USA email: info@cosmobiousa.com URL: www.cosmobiousa.com Phone/FAX: (+1) 760-431-4600