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Mount Quick

Mount-Quick "Aqueous"

<b>Purpose of use</b>	Sample mounting medium / munter for permanent preservation	
<b>Features</b>	<p><b>Mount Quick</b></p> <ul style="list-style-type: none"> <li>• Rapid drying, highly transparent, and does not form bubbles. Highly soluble in xylene, does not degenerate or stain the specimen. Can withstand long-term storage after mounting.</li> <li>• Tube design provides easy control of amounts used and limits air exposure, preventing variability in concentration.</li> </ul>	<p><b>Specifications:</b></p> <p>Specific gravity: 0.89sp; Viscosity: 33cp; Refractive index: 1.41</p>
	<p><b>Mount-Quick "Aqueous"</b></p> <ul style="list-style-type: none"> <li>• Maintains superior transparency when fully hardened.</li> <li>• Low fluorescence.</li> <li>• Mounting is simple with an easy-to-use dropper bottle.</li> </ul>	<p><b>Specifications:</b></p> <p>Mobility: pale yellow, transparent, viscous aqueous solution; Specific gravity: 0.88sp; Viscosity: 33cp; Refractive index: 1.41</p>
<b>How to use</b>	<p><b>Mount Quick</b></p> <p>Squeeze out an amount approximately the same size as the sample section at one location on the cover glass or sample.</p>	
	<p><b>Mount-Quick "Aqueous"</b></p> <p>Drop onto the cover glass or section from above. (Prevent bubble formation by inverting the container and touching the tip on the bottom of the drop.)</p> <p>* Before first use, cut open the container tip at a 45° angle with a box cutter or scissors.</p>	
<b>Precautions</b>	After use, always re-cap and store in a cool, dark place.	

### Application Data

Use Mount Quick Aqueous when samples are not dehydrated and cleared (photograph as the color fades during long-term storage) and use Mount Quick Tube when samples are dehydrated and cleared (can be permanently preserved).

### Staining data chart for Mount Quick Tube and Mount Quick Aqueous

<b>DM01 Mount Quick Tube</b> <b>(non-water soluble mounting medium for mounting after dehydration and clearing)</b>	
Hematoxylin-eosin staining	○
Azan stain	○
Grimelius stain	○
Silver stain	○
Van Gieson's stain	○
Periodic acid-Schiff (PAS) stain	○
Giemsa stain	○
Immunostaining (non-aqueous DAB)	○
<b>DM02 Mount Quick Aqueous</b> <b>(water soluble mounting medium for mounting after staining without dehydration)</b>	
Lipid stains	○
Immunohistochemical staining (AEC)	○
Blood smear staining	○
PAS, FITC, DAB staining	○

Both Mount Quick Tube and Mount Quick Aqueous can be used with methyl green and BCIP/NBT staining, but are not suitable for long-term storage. After staining, we recommend photographing as soon as possible, as the color can fade in 24 hours after drying, dehydrating, clearing, and mounting.

BCIP/NBT can stain red and black with immunostaining.