For safety, this manual must be read carefully and the device be handled properly.

Cassette Electrophoresis Unit

Model DPE-1020 for Mini Gel

Cat. No. 303111

Version: June 12, 2019

www.cosmobiousa.com

Symbols and Conventions

Caution! This symbol indicates risk of harm or personal injury. Always consult the Instruction Manual before touching the area of the Cassette Electrophoresis Unit Model DPE-1020 that displays this symbol.

Caution! This symbol indicates risk of personal injury or harm by electrical shock. Always consult the Instruction Manual before touching the area of the Cassette Electrophoresis Unit Model DPE-1020 that displays this symbol.

Electric shock hazard

Cautions

This electrophoresis device must be operated together with an external DC power—supply designed specifically for electrophoresis applications under a high voltage. The output of this power supply must be isolated from external ground to insure that the DC voltage output floats with respect to ground. COSMO BIO is not responsible for any injury or damage caused by use of this device for purposes other than those for which it is intended or by modifications of the device.

The electric operating parameters for this device are as follows:

Voltage Limit 0-600 Vdc Current Limit 0-100 mA

Ambient Temperature Limit 0~35℃

Protective gloves and safety glasses must be worn when operating in a laboratory environment. This device is intended for in vitro use only and must be used with its Lid whenever it is operated. Gel Plates that are separately available (Please refer to Related Products) are made of glass. The Glass-Made-Plates are breakable if pressed improperly when attached or detached. Plates should be set again into the device when difficult to be assembled properly.

Limited Warranty

This device is warranted against manufacturing defects for a period of one year from the original purchase date. COSMO BIO will repair/exchange the defective parts under this warranty subject to the terms and conditions of repair of COSMO BIO provided that the device is operated according to the manual and labels.


www.cosmobiousa.com

I . Features

- 1. Unique Wedge system to firmly and easily fix the Gei Plate
- 2. Buffer filled in the Chambers to cool all the surface of the Gel Plate
- 3. Smaller quantity of buffer required
- 4. Easier-sample loading visible through the transparent body

II . Components

| Description | | Quantity |
|-------------|--|----------|
| 1 | Lid with power cables | 1 |
| 2 | Anode Buffer Chamber | 1 |
| 3 | Cathode Buffer Chamber Frame with electrodes | 1 |
| 4 | Gasket | 2 |
| (5) | Wedge | 2 |
| 6 | Acrylic Plate | 2 |
| 7 | Instruction Manual | 1 |

III . Specifications

| Lid | Acrylic |
|-------------------------------|------------------------|
| Anode Buffer Chamber | Acrylic |
| Cathode Buffer Chamber Frame | Acrylic |
| Gasket | Silicone |
| Wedge | Acrylic |
| Acrylic Plate | Acrylic |
| Overall size (LxWxH) | 8 × 18 × 18 cm |
| Weight (before buffer filled) | 0.5 Kg |
| Precast Gel Compatibility | MULTIGEL® II mini Gels |
| Electrical Limits | 600 V dc and 100 mA |
| Buffer Chamber Capacity | 500 mL |

Additional equipment and materials required (not included)

- Power supply
- Electrophoresis buffer (Please refer to Related Products)
- Sample loading solution (Please refer to Related Products)


Note: This device is not compatible with acetone or ethanol. Use of organic solvents voids the warranties.

www.cosmobiousa.com

IV . Unpacking


- 1. To remove the Lid, hold the device firmly and slide the Lid as indicated on the top.
- 2. Slide the Wedges towards the center of the Outer Chamber (Anode Buffer Chamber) to release the Inner Chamber (of Acrylic Plates and Cathode Buffer Chamber Frame).
- 3. Take out the Inner Chamber.


V . Assembly and Operation

Caution: Electrophoresis operates at a high voltage and glass-made Gel Plates are fragile by nature; some protective measure like wearing gloves should be taken when handling DPE-1020 not to injure yourself, hands, etc.

- 1. Ensure that the power supply is turned off and the Lid is not connected.
- 2. Ensure the Gaskets are set on the Inner Frame (Cathode-Buffer Chamber Frame). Otherwise place the Gaskets on the grooves of the Inner Frame and press the Gaskets so that they are fixed firmly into the grooves.
- 3. Pour the buffer (approximately 200 mL) in the Outer Chamber (Anode Buffer Chamber) up to around 4 cm from the bottom. Then place the Inner Frame into the Outer Chamber, matching respective colors of the plugs (red and black) to those as indicated on the Outer Chamber.

- 4. Take out a Gel Plate (Please refer to Related Products) from the bag and remove the comb gently. Insert the Gel Plate(s) placing the grooved side facing towards the inside of the Inner Frame. If only one piece of Gel Plate is placed, put the Acrylic Plate on the other side of the Inner Frame to fully build the Inner Chamber.
- 5. Stabilize the Inner Chamber by inserting the Wedges straight down into the corners between the Gel Plate(s) or Acrylic Plate and the Outer Chamber.
 - (Note: Place the grooved end of each Wedge facing outside, otherwise the Lid cannot be set place. Do not insert the top of the wedge below the top of the gel plate. Inserting the wedge too far can warp the gel plate, causing a gap between the glass plate and gel.)


www.cosmobiousa.com


6. Fill the Inner Chamber with running buffer up to 2 mm below the top edge of the Gel Plate(s).

After a few minutes, check for leakage (lowering of the level of the inside buffer). If leakage is observed, disassemble the Inner Chamber and repeat the above steps. (Note: Avoid a short circuit, do not overfill the Inner Chamber.)

7. Fill the Outer Chamber with running buffer up to 2.5 cm below the top edge of the Outer Chamber. (Note: Pour the buffer gently to avoid trapping air bubbles along the bottom of the Inner Chamber. Minor bubbles can be ignored.)

- 9. Ensure the power cables of the Lid are not connected to the power supply.
- 10. Place the Lid on the Outer Chambers and slide it so that the plugs can be connected with the jacks firmly.

- 11. Connect the power cables of the Lid to the power supply with the proper polarity. (Note: When attaching the Lid, the power supply must be off.)
- 12. Turn the power supply on to start electrophoresis.

Recommended conditions (for your reference)

SDS-PAGE: 30 mA DC constant for about 60 minutes a Gel Plate

Native-PAGE: 15 mA DC constant for about 100 minutes a Gel Plate

DNA-PAGE: 15 mA DC constant for about 100 minutes a Gel Plate

- 13. Turn the power supply off to stop electrophoresis when the tracking dye (Bromophenol blue or Xylene cyanol FF) reaches approximately 5 mm (around the center of Gasket) from the bottom edge of the Gel Plate(s).
- 14. Disconnect the power cables of the Lid from the power supply. Holding the device firmly, slide the Lid as indicated on the top for open.
- 15. Discard the buffers and release the Gel Plate (s) by sliding the Wedges towards the center of the Outer Chamber.
- * Staining or other procedures to follow the electrophoresis should be operated according to the individual instructions.
- * The device should be washed after each run with a gentle detergent and rinsed thoroughly with distilled or de-ionized water.

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