2D-Silver Stain Reagent II -For Electrophoresis-

Cat. No. DCB-423413

Last Updated: 2016/1/7

www.cosmobio.com

Background

Detection of proteins and nucleic acids migration by polyacrylamide gel electrophoresis and silver staining are now being highlighted as one of the most sensitive staining methods. However the original silver staining method is not necessarily easy to work and takes a long-time to complete. 2D-SILVER STAIN-II is an improved reagent kit developed to provide simple and fast staining.

[||] Features

- 1) Rapid staining results in a short time after electrophoresis
- 2) Staining with higher sensitivity

For Proteins: 50 to 100 times more sensitivity than CBB staining

For Nucleic Acids: 50 to 100 times more sensitivity than EB staining

- 3) Simple preparation and simple operation
- 4) No use of materials such as heavy metals which need to be controlled by regulations.
- 5) Staining process can be stopped any time during development for desired chromatic figures.

[|||] Applications

2D-SILVER STAIN II is applicable to the detection of proteins and nucleic acids subjected to polyacrylamide and SDS-polyacrylamide slab gel electrophoresis.

[IV] Components

Reagents	Main Ingredients	Volume
Fixing Reagent	Thiourea	100 ml x 1
Pretreatment reagent	Dithiothreitol, Glutar-aldehyde, Thiourea	100 ml x 1
Staining Solution A	Silver nitrate	100 ml x 1
Staining Solution B	Ammonium hydroxide, Sodium hydroxide	100 ml x 1
Concentrated developer	Citric acid, Folmaldehyde, Sodium thiosulfate	100 ml x 1
Stopper	Citric Acid	100 ml x 1

[V] Procedures

www.cosmobio.com

For Protein Staining, in addition to the reagents supplied with the kit, the following materials are required.

- 1) Methanol
- 2) Acetic Acid
- 3) Deionized water (less than 10⁻⁶ mho conductivity)
- 4) Measured cylinder
- 5) Flat-bottommed tray for developing gels
- 6) Measuring pipette

[V-1] Preparation of Reagents

Prepare the following solutiosn just before staining the gel. The volume of reagents use for preparation shown belowis based on teh size of 140 mm x 140 mm x 1.0 mm slab gel. Fro any different sizes of slab gels, the voulme of reagents should be prepared in proportion to the size of the gel to be used.

1	Fixing solution I Mix 100 ml of methanol and 20ml of acetic acid with 80ml of deionized water. Stir and use the mixture as fixing solution I.
2	Fixing solution II Mix 60 ml of methanol, 20 ml of acetic acid and 10ml of 1 Fixing Reagent with 110 ml of deionized water. Stir and use the mixture as fixing soultion II.
3	Pretreatment solution Mix 100 ml or methanol, 10 ml of 2 Pretreatment Reagent and 90ml of deionized water. Stir and use the mixture as pretreatment solution
4	Silver staining solution Mix 10ml of 3.) Staining Solution A and 10 mL of 4.) Staining Solution B and stir. Then add 180ml of deionized water and stir. Use the mixture as silver staining solution.
5	Developer Add 10 ml of 5.) Concentrated Developer to 190 ml of deionized water and stir. Use the mixutre as the developer

[V-2] Staining Procedure

[V-3] Notes

- As the silver staining solution produces explosive silver amide when left to stand, the solution must be treated with hydrochloric acid or sodium chloride to convert the silver into silver chloride (See Procedure 6.)
- Since Staining Solutions B among the components of 2D-Silver Stain-II, contains ammonium hydroxide and sodium hydroxide as an alkali substance respectively, take care to prevent them from contacting the eye and the skin. In case of eye contact rapidly wash the eye with flowing water, and consult a doctor. In case of contact with the skin and clothes, rapidly wash away with water, and consult a doctor if inflammation occurs.
- Staining Solution A and other components may also cause a stimulus to the mucous membrane and the skin. In case of contact with the eye, skin, and clothes, rapidly wash away with water, and consult a doctor if inflammation should be caused.
- Addition of the silver staining solution to the sample gel may color or make turbid the solution. Such coloration and tubidity do not affect the stain. (See Procedure 5.)
- Deionized water or distilled water of less than 10 -8 mho conductivity is recommended.
- The tray to be used for the staining should be clean and have a smooth surface.
- When a sample gel stained by the CBB staining is to be subjected to the silver staining, the sample gel should be washed with deionized water thoroughly. (3 hours or longer while changing deionized water several times prior to the silver staining.)
- Care should be taken not to damage the gel during operation.
- Care should be taken to see that the gel may not come out of the surface of solution.
- Use a pair of disposable gloves in staining operation not to touch the gel.
- If a proper shaker is not available, shake the tray by hand.
- The sensitiviey of a gel thicker than 1mm declines somewhat as compared with a thinner gel.
- igoplus Keep the reagents in a tightly closed container and refrigerated, and use within 6 months after unsealing.
- Use the reagents in a tightly closed container adn refrigerated, and use within 6 months after unsealing.
- Use the reagents (methanol, acetic acid) of superior quality

[V-4] Package Size and Storage

Enough to stain 10 slabs of gels (140 mm x 140 mm x 1.0 mm)

The kit should be stored in a dark place at 2-10 $^\circ$ C

The term of validity: Printed on the package.

[V-5] References

1. K. Ohsawa et. Al.: Silver Stain for Detecting 10-Femotogram Quantities of protein after Polyacrylamide Gel Electrophoresis, Anal. Biochem., 135, 409(1983)

2. S. Irie: A highly sensitive silver Staining for detection of proteins in polyacrylamide gels, biochemistry(Japan), 52, 411(1980)

3. B.R. Oakley et. al.: Visualization of proteins with silver "stain" for detecting proteins in polyacrylamide gels, Anal. Biochem., 105, 361(1980)

4. H.M. Poehling et. al.: Visualization of proteins with silver "stain", a critical analysis, electrophoresis, 2, 141 (1981)

5. U.K. Laemmli: Cleavage of structural proteins during the assembly of the head of bacteriophage T4., Nature, 227, 680(1970)

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