For Research Use Only

Multi Gel Dryer

- **Starter Set**

Cat. No. (Japan): 423512

Cat. No. (Outside Japan): DCB-423512

- **Reagent Set**

Cat. No. (Japan): 423505

Cat. No. (Outside Japan): DCB-423505

Created: March 18th, 2026

www.cosmobioussa.com

I. Introduction

The Gel Dryer was developed as an easy-to-use system for drying polyacrylamide gels into film after electrophoresis and staining. It's also ideal for preserving electrophoresis results using our Multi Gel® II Mini precast gels.

II. Features

1. Enables long-term storage of electrophoresis results.
2. No need to prepare reagents—simple operation.
3. Even high-concentration gels dry beautifully with minimal cracking.
4. No special equipment required.

III. Intended Use

For drying and preserving mini-slab polyacrylamide gels (PAG) and SDS-polyacrylamide gels (SDS-PAG) after electrophoresis and staining.

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IV. Packaging & Contents

Each set contains the following items.

	Starter Set (for 20 sheets) Cat. No.: 423512	Reagent Set (for 20 sheets) Cat. No.: 423505
Mounting Frames 2 sets	○	
Gel Drying Reagent 500 mL × 2	○	○
Pre-cut Cellophane 50 sheets	○	○


V. Storage Instructions

Keep tightly sealed at room temperature. Keep away from open flames (Gel Drying Reagent).


VI. Expiration Date

Labeled on the outer packaging

VII. Frame Components

<Assembly Diagram>

Each frame set includes the following parts:

① Clips 4 pieces


· Upper Frame 1 piece (Outer: approx. 14 × 14 cm, Inner: approx. 10.4 × 10.4 cm)

· Lower Plate 1 piece (approx. 14 × 14 cm)


* Pre-cut cellophane size: approx. 14 × 14 cm

VIII. Instructions for Use


- After staining (CBB, silver stain, etc.), immerse each gel in 50 mL of purified water and shake for 10 minutes. Replace the purified water and repeat this washing process three times in total.


- Pour off the purified water, then add 50 mL of Gel Drying Reagent per gel. (Insufficient reagent can cause the gel to crack.) Gently shake for 60 minutes.

- Soak two pre-cut cellophane sheets per gel in purified water until fully expanded. Place one sheet on the lower plate of the fixing frame.


- Carefully place the gel on top of the pre-cut cellophane sheet, making sure no air bubbles are trapped. Pour a small amount (a few mL) of the Gel Drying Reagent (the same one used in step 2) over the gel.


- 5 Place another precut cellophane sheet over the gel, covering the sides and making sure no air bubbles get trapped between the gel and the cellophane.


- 6 Position the top frame over the assembly.


- 7 Secure all four sides with clips.

- 8 Hold the frame at a slight angle and gently shake to drain off excess liquid.

- 9 Let the frame sit flat and undisturbed at room temperature for 1–2 days to dry completely.

- 10 Once the gel is fully dried, remove it from the frame, trim away any excess cellophane, and store it in a moisture-free environment. (Using an album or similar storage that is convenient.)

IX. Usage Notes

1. Always use the Gel Drying Reagent in a well-ventilated area. Do not ingest or inhale the vapors for extended periods. If it comes into contact with your eyes or mucous membranes, rinse thoroughly with water right away.
2. Keep the Gel Drying Reagent away from heat sources and areas above 40°C. Since the reagent is volatile, do not leave the cap off.
3. If the gel has any tears, those spots are more likely to crack. Handle the gel with care after electrophoresis. If you find any tears, trim them off with a knife before drying for best results.
4. Choose a container large enough for the gel to be fully immersed in the reagent and make sure it can be closed with a lid.
5. Avoid drying the gel with heat sources like dryers, as this can cause cloudiness in the gel or deform the frame.
6. While the gel is drying, try to keep it undisturbed and avoid touching or jarring it unnecessarily.
7. If you notice any cloudiness in the gel, leave it in the frame, gently rinse it with water to clear the haze, wipe away any excess moisture, and allow it to dry again.
8. When drying gels thicker than 1 mm, increase the amount of Gel Drying Reagent according to the volume of the gel and extend the processing time. Please also note that thicker gels are more prone to cracking.
9. Dried gels can absorb moisture, which may cause warping or cracking. Store dried gels in a dry place to prevent this.
10. Reusing Gel Drying Reagent can lead to cracks in the gel, so please avoid repeated use.
11. The drying condition of the gel may vary depending on gel concentration, temperature, and humidity. If you see cloudiness, shorten the treatment time with Gel Drying Reagent; if cracks appear, extend the treatment time as needed.

12. The gel may shrink slightly after treatment with the Gel Drying Reagent. Please keep this in mind when interpreting data.
13. The Gel Drying Reagent is classified as a hazardous material. Category 4, Class 2 Petroleum, Hazard Class III, Methanol Mixture, Water Soluble

X. Troubleshooting

Phenomenon	Cause & Solution
Gel cracks	<ul style="list-style-type: none"> * Insufficient rinsing <ol style="list-style-type: none"> 1. Make sure to rinse thoroughly. * Gel Drying Reagent not applied sufficiently <ol style="list-style-type: none"> 1. Shake well during treatment to ensure full coverage. 2. Extend the treatment time as needed. * Gel has tears <ol style="list-style-type: none"> 1. Take care to avoid creating tears in the gel. 2. If the gel already has tears, trim away those sections with a knife before drying. * Drying is taking too long <ol style="list-style-type: none"> 1. After setting the frame, remove any excess Gel Drying Reagent as much as possible. 2. If humidity is high, such as during summer, try to dry the gel in a well-ventilated area. * The gel was subjected to impact while drying. <ol style="list-style-type: none"> 1. During drying, avoid impacts and unnecessary handling.
The gel appears cloudy.	<ul style="list-style-type: none"> * Insufficient rinsing. <ol style="list-style-type: none"> 1. Make sure to rinse thoroughly. * Drying time was too short. <ol style="list-style-type: none"> 1. Avoid heating during drying. (If cloudiness occurs, please refer to item 7 in the "Usage Notes.")
Gel deformation or cracking during storage may occur.	<ul style="list-style-type: none"> * Insufficient drying. <ol style="list-style-type: none"> 1. Be sure to dry thoroughly. * The gel absorbed moisture. <ol style="list-style-type: none"> 1. After drying, store the gel away from humidity.

XI. References

- Juang R.H, et al.: Over-Drying Method for Polyacrylamide Gel Slab Packed in Cellophane Sandwich Anal. Biochem, 141, 348-350 (1984)
- Samal B.B.: Drying and Storage of Polyacrylamide Slab Gels; A Simple Procedure Anal. Biochem, 163, 42-44 (1987)

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