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TEL. 011-641-7507 FAX. 011-644-9209

Reagent for Measuring Potato Virus Y-N Strain

Product No. 309401

PVY-N ELISA HOKUDO*Hokudo***Reagent configuration**

96-well plate x 1

(Stored at 2-8°C)

Components	Packaging	Quantity	Remarks
Antibody solid-phase plate (Anti-PVY IgG solid-phase plate)	96 wells/plate	1	Dry plate
Enzyme-labeled antibody stock solution (ALP-labeled Anti-PVY IgG)	0.3 mL	1	100-fold concentrated solution
Concentrated homogenate solution (Phosphate buffer, Tween20)	100 mL	1	5-fold concentrated solution
Reaction buffer (PBS, BSA, Tween20)	26 mL	1	
Concentrated washing solution (Phosphate buffer, Tween20)	100 mL	2	5-fold concentrated solution
Enzyme substrate solution (PNPP solution)	26 mL	1	
Reaction stop solution (2 mol/L NaOH)	13 mL	1	
Positive control	1 mL	2	
Negative control	1 mL	2	
Plate seal		2	

*Containing 0.1% sodium azide: enzyme-labeled antibody stock solution, concentrated homogenate solution, reaction buffer, concentrated washing solution, positive/negative control

Operating method

■ Materials to be prepared

- Purified water, microplate reader (Dominant wavelength: around 405 nm), centrifuge (Capable of operating at 3000 rpm and 4°C)
- Micropipette (Variable type 200, 1000 µL) • microplate washer/reservoir
- Multichannel micropipettes (50-200 µL), disposable tubes (1.5 mL, 15 mL), paper towels

■ Preparation of reagents

1. Preparation of enzyme-labeled antibody solution: Dilute the enzyme-labeled antibody stock solution 100-fold using the reaction buffer, according to the required number of wells.
Example: 6 rows, 48 wells = 120 µL + 11.88 mL (Reaction buffer)
2. Preparation of the homogenate solution: To 100 mL of the concentrated homogenate solution, add 400 mL of purified water to obtain a 5-fold dilution.
3. Preparation of the washing solution: To 100 mL of the concentrated homogenate, add 400 mL of purified water to obtain a 5-fold dilution.


■ Measurement method **Caution: Allow each reagent to stand at room temperature before use.**

1. Preparation of samples for measurements

- 1) Weigh 1 g of each sample and cut into small pieces, with scissors in a mortar.
- 2) Add 5 mL of the diluted homogenate solution and mill it in a mortar.
- 3) Transfer the milled solution to a 15 mL centrifuge tube.
- 4) Afterwards, add 5 mL of the homogenate solution to the mortar and wash, and transfer the resulting washing solution to the same tube.
- 5) Centrifuge the mixture at 3000 rpm for 10 minutes at 4°C, and collect the supernatant for use as the 10-fold diluted sample.

2. Primary reaction

- 1) Remove the necessary number of antibody solid-phase plates from their aluminum bags. (Measurements in triplicate are recommended for the control.)
- 2) Add 200 µL each of the control and sample to wells using a micropipette.

- 3) Attach a plate seal, and allow it to react at room temperature for 1 hour.
3. Wash each well of the plate 5 times with washing solution, using the washing plate washer.
4. Secondary reaction
 - 1) After completely removing any washing solution that remains in the wells (by gently tapping the wells on a paper towel), transfer the necessary amount of the diluted enzyme-labeled antibody solution to the reservoir, and add 200 μ L of the solution to each well with a micropipette.
 - 2) Attach a plate seal, and allow it to react at room temperature for 1 hour.
5. Wash the wells in the same manner as described in 3. Cleaning.
6. Enzymatic reaction
 - 1) After completely removing any washing solution that remains in the wells, transfer the necessary amount of the enzyme substrate solution to the reservoir, and add 200 μ L of the solution to each well with a micropipette.
 - 2) Allow it to react at room temperature for 20 minutes.
 - *Be sure to use separate reservoirs for the enzyme substrate solution and the enzyme-labeled antibody solution.
 - *Do not return the enzyme substrate solution to the reagent bottle, once it has been transferred to the reservoir.
7. Stopping reaction After the stopping the reaction, add 50 μ L of the reaction stop solution transferred to the reservoir to each well using a micropipette, and mix gently.
8. Measure the absorbance Measure the absorbance at 405 nm using a microplate reader for absorbance measurement.
 - *Measure the absorbance immediately after stopping the reaction.
9. Evaluation of the results

Evaluate the absorbance of each sample by confirming that the absorbance of the negative control is not more than 0.2, and that of the positive control is not less than 1.0.

Precautions for operation

1. Read the instruction manual carefully before use, and make sure to understand the contents fully prior to making measurements.
2. Add the enzyme-labeled antibody solution or enzyme substrate solution promptly after washing.

Precautions for use

1. Do not expose the reagents to direct sunlight during storage and use.
2. Please strictly observe the specified reaction temperature and reaction time.
3. Allow the kit components to stand at room temperature (20-30°C) before use.
4. Wash the components well, as improper washing may result in variations in measurement results.
5. To prevent contamination between wells, take care to avoid foaming or adhering of the solutions to surrounding areas.
6. After unpacking the antibody solid-phase microplate, securely close the aluminum bag for storage.
7. Avoid microbial contamination of the samples and reagents, as well as contamination between reagents.
8. Since the enzyme substrate solution is easily oxidized by metal ions, always use new disposable tools, and use a disposable pipette when transferring the enzyme substrate solution from the reagent bottle to the reservoir. Do not return the enzyme substrate solution to the reagent bottle, once it has been transferred to the reservoir.
9. Handle the reagents with care, to prevent them from coming into contact with eyes, skin, etc. In case of contact with a reagent, wash the reagent off with copious water, and consult a specialist.

Storage

Store at 2-8°C.

Manufacturer

Expiration date

6 months after manufacture (indicated on the outer box and label)

Hokudo Co., Ltd.

10-4-28, 9 Jo-nishi, Hachiken, Nishi-Ku
Sapporo-shi, 063-0849
TEL.011-641-7507 FAX.011-644-9209