For general study

## Chitosan Nanofiber Coated Culture ware

● Culture plate

● Culture coverslip

FOR RESEARCH USE ONLY

This product is a cell culture ware that coats the Chitosan Nanofiber. Chitosan is a polysaccharide polymer produced by alkali processing (deacetylation) of chitin, which is natural material found in various living things such as shellfish (shrimp, crab), insects and mushrooms. Unique properties of Chitosan, such as biocompatibility and biodegradability, can be designed by controlling the molecular weight of Chitosan and the degree of deacetylation, conversion rate to the Chitosan.

One of distinguishing properties of Chitosan nanofiber is to become scaffolding material in cell culture. With Chitosan nanofiber, adhesiveness of cultured cell can be improved (compared to the culture in polystyrene plate or glass cover slip).

Chitosan nanofiber is hardly perishable and is an easy-to-use material for cell/tissue culture, compared to poly-L-lysine or collagen coated culture dish.

### ● Product information

CODE	WARE	KIND	QUANTITY	PACKAGING	STORED	EXPIRATION DATA
HSC13	cover slip	13mmφ	5	individual	room temperature	1 year
HSC15	cover slip	15mmφ	5	individual	room temperature	1 year
HSP12	plate	12well plate	1	individual	room temperature	1 year
HSP24	plate	24well plate	1	individual	room temperature	1 year
HSP96	plate	96well plate	1	individual	room temperature	1 year

We can produce 18,22mmφ cover slip and 48well plate according to your demand.( Entrusted production)


### ● Product specification

ITEM		
Deacetylation of Chitosan nanofiber	>90%	
Amount of fixed Chitosan nanofiber	0.05 ~ 0.1mg/cm <sup>2</sup>	For cover slip, Chitosan nanofiber is fixed on one surface.
Sterilization method	Ethanol	Ware is wet by ethanol.

### ● Stability

The expiration date is one year or more from the shipment day in the room temperature keeping. There is no problem in the quality though the Chitosan nanofiber might dry while keeping it and become white.

Manufactured by

Distributor

COSMO BIO Co., LTD.  
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>

e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

Phone : +81-3-5632-9617

FAX : +81-3-5632-9618


•how to use

**【cover slip】**

- Chitosan nanofiber is fixed on one surface of cover glass (cell culture side).
- The cover slip is doubly packaged by aluminum laminate packing (pack on the outside) and clear packing (pack on the inside). There are five clear packing in one aluminum laminate packing, and there is one cover slip in one clear packing.

- cover slip is packed by a clear packing so that Chitosan nanofiber-fixed side (cell culture side) turn to the side red seal is stuck. Please do not mistake the cell culture side.
- There is ethanol for disinfection in a clear packing for ethanol sterilization.
- When opening, at first, please open from the incision of the aluminum laminate packing and take out clear packings with tweezers etc. Next, the upper side of clear packing is cut with scissors etc., and please take out a cover slip with tweezers etc. and use it.

- The cover slip is wet in the ethanol. Please use it after leaving about 30 minutes in the room temperature in a clean bench. There is no problem in the quality though the Chitosan nanofiber(cell culture side) might dry while keepings it and become white. Please use it after wetting the culture solution etc enough.

**【multiwell plate】**

- The Chitosan nanofiber (cell culture part) is fixed at the bottom of well.
- Plate is singly packed by the vacuum packing. Please open the vacuum packing with scissors etc., and use it.

•References

- 1) Tsuneo Ohkuma ,et al.,Chitin and Chitosan Research12(2) 190-191(2006)  
「Preparation of Chitosan Nanofiber Mesh as a Scaffold for Regenerating Nerve Tissue」
- 2) Katsumi Konno ,et al.,Chitin and Chitosan Research13(2) 124-125(2007)  
「Manufacturing chitosan nanofiber by electrospinning process」
- 3) Katsuyoshi Sakai ,et al.,Chitin and Chitosan Research13(2) 126-127(2007)  
「Manufacturing cell culture substrates on chitosan nanofiber and manufacturing chitosan gel,chitosan sponge」
- 4) Kazuhiko Watabe ,et al.,Chitin and Chitosan Research13(2) 128-129(2007)  
「Neural tissue culture on chitosan nanofiber matrices」

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### Chitosan Nanofiber Coated Culture ware

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COSMO BIO CO., LTD.