Synthetic Zeolite Sheet-Type High-Performance Desiccant G

- Enables ultra-low humidity storage Humidity can be reduced to nearly zero.
- Also absorbs corrosive gases
- Convenient individually wrapped packets

Ideal for storing metals and pharmaceuticals, and perfect for situations where traditional desiccants like silica gel can't maintain sufficiently low humidity—give it a try for glove boxes, desiccators, and more.


This product features strong sheet-shaped synthetic zeolite 4A, with each piece individually sealed in an aluminum-laminated bag.

It not only removes moisture from the air but can also adsorb gases like H₂S, SO₂, and NH₃ for versatile use.

The main ingredient, synthetic zeolite 4A, is a type of crystalline hydrated aluminosilicate. Its chemical formula is Na₁₂AlO₂SiO₂12·27H₂O. The crystal surface features countless uniform micropores measuring 4Å, which are connected to the inner cavities. Sodium ions within the crystal structure also play a role in adsorption. This product offers both physical adsorption by van der Waals forces and chemical adsorption through ion exchange, making it especially effective in ultra-low humidity environments.

1. Compared to silica gel and other desiccants, this product absorbs moisture much faster and maintains high absorption even at very low humidity levels. This allows for storage in much drier conditions—humidity can be reduced to almost zero. Additionally, unlike silica gel, it does not release the absorbed moisture back into the environment.
2. Each unit is wrapped in PP film, preventing any powder from the desiccant from escaping.
3. Chemically stable, this desiccant retains its shape after use and is neither corrosive nor deliquescent.
4. Contains almost no heavy metals and uses no additives harmful to humans.
5. The shape is consistently uniform, offering excellent compressive strength and outstanding mechanical properties.

Isothermal Moisture Absorption Curve (25°C)


AZ-10G can absorb and hold up to about 2g of moisture, which is approximately 20% of its 10g weight. This is equivalent to the amount of saturated water vapor in 100L at 22°C (or about 80L at 25°C).

AZ-10G can adsorb the following substances:


In addition to hydrogen sulfide and sulfur dioxide shown below, it can adsorb molecules with an effective diameter of less than 4 nm.

Examples: CO₂, C₂H₅OH, etc.

Absorption of sulfur dioxide gas and ammonia (25°C)

Absorption of hydrogen sulfide (25°C)


Product Details

Catalog No.	Product Name	Packaging	Weight & Dimensions (H × W × D)
AZ10G-100	Synthetic Zeolite Desiccant AZ-10G	Individually packed in aluminum pouches	10g 60 × 40 × 4 (mm)

For storage, please use a sealed container such as a desiccator, or a dedicated high-barrier storage bag.

Humidity Changes in the Atmosphere When Using Oxygen Absorber (A-750HS) Together with a Desiccant

* Silica gel: Type A, 10g.

* Synthetic Zeolite Plate-Type High-Performance Desiccant, 10g

Synthetic Zeolite Desiccant

User Guide

Please read this manual carefully and use the product as directed.

About Synthetic Zeolite Desiccant

Synthetic zeolite 4A, the main ingredient, is a type of crystalline hydrated aluminosilicate. Its surface contains countless uniform pores measuring 4\AA , which connect to internal cavities. Sodium ions found within the crystal also play a role in adsorption. Its absorption properties are due to both physical adsorption driven by van der Waals forces and chemical adsorption involving ions. This desiccant is particularly effective in low humidity environments.

Features

1. Compared to silica gel and other desiccants, it absorbs moisture much faster and maintains high absorption even at low humidity levels, making it ideal for preserving items in very dry conditions. Additionally, unlike silica gel, it does not release absorbed moisture back into the environment.
2. Packed in PP film, you won't need to worry about powder leakage from the desiccant.
3. Chemically stable, with no change in shape or appearance after moisture absorption, and it is neither corrosive nor deliquescent.
4. Contains almost no heavy metals and uses no additives harmful to humans.
5. The product features a uniform shape, high compressive strength, and excellent mechanical properties.

Product Packaging Details and Precautions

- Synthetic Zeolite Sheet-Type Desiccant-10 is packed in perforated PP film and then individually wrapped in an aluminum pouch. Please use the perforated PP film as is—do not tear it.
- Sometimes the aluminum pouch may appear slightly puffed up. This is due to the release of air (oxygen and nitrogen) that was absorbed during manufacturing.

How to Use

1. Place your item in a dedicated storage bag or an airtight glass container.
2. Open the aluminum pouch, take out the Synthetic Zeolite Sheet-Type Desiccant-10, and put it into the storage bag or glass container.
3. Seal the storage bag using a heat sealer or a sealing clip. For glass containers, make sure they are tightly closed to keep out external air.

Precautions for Use

- If your item has protrusions that could damage the storage bag, wrap it first in a breathable container, or use paper or a poly bag (unsealed) before placing it in the storage bag.
- Once removed from the aluminum pouch, the desiccant will start working immediately. Please complete the process quickly.
- AirMake sure neither side of the air-permeable surface is blocked by the storage bag film or the items inside. Pay extra attention if using a vacuum sealer—if both sides are covered, the desiccant will not activate.

Using with Oxygen Absorbers (for creating an oxygen-free, dry environment)

If you need a dry and oxygen-free storage environment, try using the A-25AZS (Oxygen-Free & Dry Storage System).

Please note: If the moisture content of the oxygen absorber is depleted, its ability to absorb oxygen will be lost and it will no longer remove oxygen that enters through the storage bag.

Specifications

Ingredients: Synthetic Zeolite 4A $\text{Na}_{12}\text{AlO}_2 \cdot \text{SiO}_2 \cdot 27\text{H}_2\text{O}$ Weight: 10 Dimensions: $60 \times 40 \times 4$ mm

Moisture absorption: About 2g at 10% relative humidity

Other substances absorbed: H_2S , SO_2 , NH_3 and more

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