## **IST Fluolid-W**

### **Oligonucleotide Amine Labeling Kits**

IST Fluolid-W Orange 600 Oligonucleotide Amine Labeling Kit

IST Fluolid-W Yellow 540 Oligonucleotide Amine Labeling Kit

IST Fluolid-W Green 520 Oligonucleotide Amine Labeling Kit

#### **Contents**

IST Fluolid-W succinimidyl ester in DMSO	60 $\mu$ l	× 3
0.2 M Sodium bicarbonate buffer (pH8.5)	120 $\mu$ l	× 3

#### **Application**

Labeling of amine-modified oligonucleotides (5 – 10 nmol amino group)

#### **Characteristics**

- 1. high quantum yield in the solid state**
- 2. no photobleach**
- 3. high stability for light, heat and pH**
- 4. labeling rate is higher than that of traditional dye**


#### **Attention**

Product is research use only.

The fluorescense unit of "IST Fluolid" has high stability under various conditions.

To keep the activity of the succinimidyl ester group, please store at 2-8 °C.

Contents are subject to change without notice.

#### **Manufacture**

International Science Technology CO., LTD.  
Kasuga Lab.

#### **Distributor**

COSMO BIO CO., LTD.  
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

<http://www.cosmobio.co.jp>

e-mail : [export@cosmobio.co.jp](mailto:export@cosmobio.co.jp)

Phone : +81-3-5632-9617

FAX : +81-3-5632-9618

## ***IST Fluolid-W Oligonucleotide Amine Labeling Kits***

**IST Fluolid-W Orange 600 Oligonucleotide Amine Labeling Kit, cat# ISU-IST001**

**IST Fluolid-W Yellow 540 Oligonucleotide Amine Labeling Kit, cat# ISU-IST002**

**IST Fluolid-W Green 520 Oligonucleotide Amine Labeling Kit, cat# ISU-IST003**

---

### **Contents**

IST Fluolid-W succinimidyl ester in DMSO	60 $\mu$ l $\times$ 3
0.1 M sodium borate buffer (pH8.5)	120 $\mu$ l $\times$ 3

---

### **Attention**

Product is research use only.

The fluorescence unit of "IST Fluolid" has high stability under various conditions.

To keep the activity of the succinimidyl ester group, please store at 2-8 °C.

Contents are subject to change without notice.

---

### **Protocol**

The kits have been optimized for labeling 5 to 10 nmol 5'-amine modified oligonucleotide (18–25 base).

Slightly shorter or longer oligonucleotides may be labeled by the same procedure. However, adjustments of the protocol may be necessary for greatly shorter or longer oligonucleotides.

It may be usable under the low DMSO concentration. In the case of the precipitations, they typically have no effect on the reaction. However, adjustments of the protocol may be necessary for yield loss.

Take care of yield loss, when the conditions of the protocol are changed. In this instance, please investigate the conditions.

### **(1) Preparation of Amine-Modified Oligonucleotide Solution**

Dissolve the 5-10 nmol of amine-modified oligonucleotide in 100  $\mu$ l of 0.1M borate buffer (pH8.5).

### **(2) Reaction Mixture (Labeling Reaction)**

Mix the 50  $\mu$ l IST Fluolid-W succinimidyl ester/DMSO solution with the amine-modified oligonucleotide solution prepared in Section (1). The reaction mixture is stirred slowly for 1-2 hours at room temperature. In the case of the precipitations, they have no effect on the reaction.

### **(3) Purification of Labeled Oligonucleotide**

The labeled oligonucleotide can be purified from the reaction mixture by preparative gel electrophoresis or reverse-phase HPLC.

One example of the purification is shown below.

#### **Purification by HPLC**

Labeled oligonucleotides can be purified by reverse-phase HPLC using a standard analytical C8 or C18 column. Load the labeled oligonucleotides solution onto the column in 0.1M triethylammonium acetate (TEAA) buffer (pH7.0) and a linear 5-95% acetonitrile gradient over 30 minutes.

- Peaks may vary depending on the length of the oligonucleotide and the purity of the sample. We recommend to run a small portion of the sample on the column first, to ascertain the identity of the peaks before running of preparative HPLC.

- To determine the identity of the peaks, monitor the absorbance at both 260 nm and at the absorbance maxima ( $\lambda_{\max}$ ) for the dye. For instruments with only one detector, two small samples should be run, each monitored at a different wavelength. Unlabeled oligonucleotide will show an absorbance at 260 nm only. Both the free dye and the labeled oligonucleotide will have absorbance at both 260 nm ( $A_{260}$ ) and at the absorbance maximum of the dye ( $A_{\max}$ ); however, the labeled oligonucleotide will have a higher  $A_{260}:A_{\max}$  ratio.

### Characteristics

Fluorescent Dye	$\lambda_{\text{abs}}^*/\text{nm}^{\dagger}$	$\lambda_{\text{em}}^*/\text{nm}^{\dagger}$	$\epsilon/\text{M}^{-1}\text{cm}^{-1}^{\dagger}$	Molecular Weight
Fluolid-W Orange 600	440	602	15,000	570
Fluolid-W Yellow 540	410	541	12,000	540
Fluolid-W Green 520	395	522	10,000	510

\* Absorbance and fluorescence emission maxima.

† Extinction coefficient at  $\lambda_{\text{abs}}$ .

If there is an unknown point, please contact us in the following.


Figure 1. Protocol for Oligonucleotide Amine Labeling with Oligonucleotide Amine Labeling Kits.

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