
Code No.KAL-KH007

For research use only

## Anti Red Sea Bream Vitellogenin Monoclonal Antibody (Clone No. 5A4)

Package Size 100  $\mu$ g (100  $\mu$ L/vial)  
 Format Mouse monoclonal antibody 1.0mg/mL  
 Buffer 10mM sodium phosphate , 0.15M sodium chloride , 1%BSA , pH7.4  
 Storage Store below -20°C  
 Once thawed, store at 4°C. Repeated freeze-thaw cycles should be avoided.  
 Clone No. 5A4  
 Subclass IgG1  
 Immunogen source Red sea bream vitellogenin purified from plasma of 17  $\beta$  -estradiol induced fish  
 Working dilution Immunoblotting 0.1~1  $\mu$ g/ml  
 ELISA 0.25~2  $\mu$ g/ml (Fixed quantity system by ELISA has not been established.)  
 Specificity < cross reactivity >

	Red sea bream	Carp	Mummichog	Medaka
ELISA	(+)	(-)	(-)	(-)
Immunoblotting	(+)	(+)	(-)	(-)

1 2 3 4

←Vitellogenin

Sample

1. Red sea bream
2. Carp
3. Mummichog
4. Medaka

### ELISA

Antigen: Red sea bream mummichog and medaka vitellogenins were purified from male fish induced with 17  $\beta$  -estradiol and carp vitellogenin was purified from female. Each antigen were diluted to 1  $\mu$ g/mL and added to microtiter plate.

Primary antibodies: Purified monoclonal antibody (5A4)


### IMMUNOBLOTTING

Sample: Plasma(dilution 1:40) from female fish (carp) or male fish induced with 17  $\beta$  -estradiol (red sea bream, mummichog and medaka).

Primary antibody: 5A4 diluted 1:1000

## Anti Red Sea Bream Vitellogenin Monoclonal Antibody (Clone No. 5A4)

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