

# Yatalase<sup>TM</sup>

from Corynebacterium sp. OZ-21

Cat#OZK-OZ-10EX

Package Size: 2 g

Yatalase is used to lyse cell walls of filamentous fungi. The product is prepared from culture supernatants of *Corynebacterium* sp. OZ-21, and consists mainly of chitinase, chitobiase and β-1,3-glucanase.

#### Feature

Has excellent thermostability and can be stored at room temperature.

Efficiently digests native chitin.

Has the revitalization of Chitinase, Chitobiase, Chitosanase, and β-1,3-Glucanarse Can be used alone to prepare protoplasts from filamentous fungi.

#### Source

Corynebacterium sp. OZ-21 Form: Lyophilized powder(containing lactose)

## Storage

Lyophilized preparation is stable for at least 1 year at 4°C. Maintain at -20°C for long term strage. Avoid freeze/thaw cycles.

## Specific activities

Chitinase activity: Approximately 50 units/g powder Chitobiase activity: Approximately 500 units/g powder

Lytic activity against cell walls: Approximately 10,000 units/g powder

#### **Properties**

Optimum pH: pH 5~8

Optimum temperature for enzyme activity: 30~50

### **UNIT DEFINITION**

#### Chitinase activity:

One unit of chitinase activity is determined as the amount required to release  $1\mu$ mol of N-acetylglucosamine from chitin in 1 minute.

#### Chitobiase activity:

One unit of chitobiase activity is defined as the amount of enzyme required to release  $1\mu$ mol of p-nitrophenol from p-nitrophenyl-N-acetyl- $\theta$ -D-glucosaminide in 1 minute.

#### Lytic activity:

One unit of enzyme activity is defined as the amount required to cause a 1% decrease in absorbance in 1 hour.

Table 1 Enzyme activities found in Yatalase

Enzyme activity	Activity (units/g power)
Chitinase	50
Chitobiase	500
Chitosanase	19
-1,3-glucanase	300
Protease	31
Cell-well lysing activit	y 10,000

Table 2 Yields of protoplasts

Species	Yield of protoplasts (× 10g-wet cells)
Aspergillus oryzae	70
Aspergillus kawachii	2
Aspergillus terreus	21
Penicillium citrinum	20
Penicillium lanosum	50
Trichoderma koningi	; 100
Monascus sp.	14
Mucor hiemalis	8
Rhizopus nigricans	21
	11
Pleurotus ostreatus	20
Coprinus cinereus	2
Lentinus edodes	



# **PROTOCOLS**

## Conditions for protoplast preparation

Strains	Culture media	Culture conditions	Conditions for preparing protoplasts
* Aspergillus oryzae	Czapek-Dox + 0.5% casamino acid (pH5.6)	30°C, 20 hours aeration (Rotary, 140 rpm)	2% Yatalase <sup>TM</sup> solution 0.6 M (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> 50 mM maleate buffer (pH5.5) 30°C, Reciprocal shaking (at 60-70 rpm), 2-3 hours
* Aspergillus kawachii	Czapek-Dox + 0.5% casamino acid (pH5.6)	30°C, 20 hours Stationary culture	2% Yatalase <sup>TM</sup> solution 0.6 M (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> 50 mM maleate buffer (pH5.5) 30°C, Reciprocal shaking (at 60-70 rpm), 2-3 hours
* Aspergillus terreus * Penicillium citrinum * Penicillium lanosum * Tricoderma koningii	Dextrin-peptone (pH5.5)	30°C, 20 hours aeration (Rotary, 140 rpm)	2% Yatalase™ solution 0.6 M (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> 50 mM maleate buffer (pH5.5) 30°C, Reciprocal shaking (at 60-70 rpm), 2-3 hours
* Mucor hiemalis * Rhizopus nigricans	2% Malt extract	30°C, 12 hours Stationary culture (Sporangiospores in germ)	2% Yatalase™ solution 0.5 M MgSO <sub>4</sub> 50 mM maleate buffer (pH5.5) 30°C, Reciprocal shaking (at 60-70 rpm), 4 hours
* Pleurotus ostreatus * Corpinus cinereus * Lentinus edodes	OSG or MYG medium (pH5.5)	25-30°C, 3-4 days Stationary culture	2% Yatalase™ solution 0.6 M MgSO <sub>4</sub> 50 mM maleate buffer (pH5.5) 30°C, Reciprocal shaking (at 60-70 rpm), 2-3 hours
* <i>Monascus</i> sp.	Dextrin-peptone (pH5.5)	25°C, 20 hours aeration (Rotary at 140 rpm)	2% Yatalase <sup>™</sup> solution 0.6 M (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub> 50 mM maleate buffer (pH5.5) 30°C, Reciprocal shaking (at 60-70 rpm), 2-3 hours



# COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME, KOTO-KU, TOKYO 135-0016 JAPAN [International/Outside Japan] [国内連絡先]

TEL: (81)3-5632-9617 / FAX: (81)3-5632-9618

TEL: 03-5632-9610 / FAX: 03-5632-9619

e-mail: export@cosmobio.co.jp / URL: www.cosmobio.com e-mail: support@cosmobio.co.jp / URL: www.cosmobio.co.jp