Labiase

(from *Streptomyces fulvissimus* TU-6)

Labiase, produced by a submerged culture of *Streptomyces fulvissimus* TU-6, is a new enzyme preparation that lyses effectively cell walls of numerous lactic acid bacterium. This preparation is used as a tool for studies of cell walls structure of lactic acid bacterium and the preparation of plasmid DNA, intracellular enzyme and protoplast from lactic acid bacterium and so on.

SPECIFICATIONS:

Appearance.....	Lyophilized powder (containing lactose)
Activity.....	≥ 10 units/g (β -N-Acetyl-D-glucosaminidase)
Other activities contained	
Lysozyme.....	$\geq 2 \times 10^4$ units/g
Optimum pH and temperature.....	See Fig. 1,2
pH and thermal stability	See Fig. 1,2

ASSAY FOR β -N-ACETYL-D-GLUCOSAMINIDASE ACTIVITY:

Unit Definition

One unit of enzyme catalyzes the release of 1 μ mole of *p*-nitrophenol from *p*-nitrophenyl N-acetyl- β -D-glucosaminide per minute at 37°C, pH4.0.

Method

Reaction mixture

Substrate.....	5mM <i>p</i> -nitrophenyl N-acetyl- β -D-glucosaminide in 0.2M
	Citrate-phosphate buffer, pH 4.0
Enzyme.....	suitably diluted enzyme
	Total volume 1.00 ml

Cat#OZ-30EX

Procedure

After incubation for 10 minutes at 37°C, terminate by 2ml of 1M Na₂CO₃. Read A₄₀₅.

Calculation

$$\text{Enzyme unit} = \frac{A_{405} \times 3}{17.7 \times 10^3 \times 10} \times 10^3$$

The molar extinction coefficient of free *p*-nitrophenol under these conditions is 17,700.

ASSAY FOR LYtic ACTIVITY TOWARD LACTIC ACID BACTERIUM:

Method

Reaction mixture

Buffer.....	50mM Citrate-sodium citrate buffer, pH 4.0	1ml
Substrate...	Lactic acid bacterium suspension(A ₆₆₀ =5.0)	1ml
Enzyme.....	500 μ g protein/ml solution	1ml
Distilled water		2ml

Total volume 5ml

Procedure

(Before adding the enzyme solution preincubation is carried out for 5 minutes at 37°C.) After incubation for 2 hours at 37°C with gentle shaking, A₆₆₀ of the mixture is determined. As a reference, 1 ml of distilled water is used instead of enzyme solution.

Calculation

$$\text{Percentage decrease in A}_{660} = \frac{(d0 - dt) - (D0 - Dt)}{d0} \times 100$$

d(0 or t):A₆₆₀ of reaction mixture after 0 or 2 hours
 D(0 or t):A₆₆₀ of reference mixture after 0 or 2 hours

Result

See Table 1.

The extent of lysis of lactic acid bacterium cells by **Labiase** varies with lactic acid bacterium strain, growth stage of lactic acid bacterium, or cultural condition.

STORAGE:

Lyophilized preparation is stable for at least 1 year when stored at 4°C.


Fig.1 pH-stability and activity

Stability:18 hours treatment at 25°C

Activity:37°C


Fig.2 Thermal-stability and activity

Stability:10 minutes treatment

Activity:pH4.0

Table 1. Lysis of lactic acid bacterium cells by Labiase

Strains	Lysis (%)
<i>Leuconostoc mesenteroides</i> IFO 3832	92.8
<i>Leuconostoc lactis</i> IFO 12455	99.6
<i>Streptococcus faecalis</i> IFO 3971	96.5
<i>Lactobacillus acidophilus</i> IFO 13951	97.9
<i>Lactobacillus plantarum</i> IAM 1041	96.8
<i>Lactobacillus casei</i> subsp. <i>casei</i> IFO 3533	59.4
<i>Lactobacillus casei</i> subsp. <i>ramnosus</i> IAM 1118	91.8
<i>Lactobacillus brevis</i> IFO 3345	95.6
<i>Lactobacillus fermentum</i> IFO 3071	70.4
<i>Lactobacillus</i> sp. (<i>buchneri</i>) IFO 3961	89.4
<i>Lactobacillus fructivorans</i> IFO 13954	96.9
<i>Lactobacillus delbrueckii</i> subsp. <i>bulgaricus</i> JCM 1002	97.7
<i>Lactobacillus delbrueckii</i> subsp. <i>delbrueckii</i> JCM 1012	92.0
<i>Lactobacillus helveticus</i> JCM 1120	98.7
<i>Lactobacillus helveticus</i> JCM 1005	84.5
<i>Lactobacillus kefir</i> JCM 5818	93.0
<i>Lactobacillus sanfrancisco</i> JCM 5668	82.1
<i>Bifidobacterium bifidum</i> JCM 1255	62.9
<i>Streptococcus salivaris</i> subsp. <i>thermophilus</i> IFO 13957	94.0
<i>Streptococcus mutans</i> JCM 5705	73.9

Please note that the above data is an example. The values vary widely depending on the condition.

Manufacturer: **ÖZEKI** OZEKI CORPORATION

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