FITC-labeled Gelatin-zymography Kit (Cosmo Bio-type)

Catalog No. AK83 (Japan)
PMC-AK83 (Outside Japan)

January 16, 2026

Introduction

Matrix metalloproteinases (MMPs) belong to the family of metalloproteinases, which consists of at least 20 members, and known to be involved in the metabolism of extracellular matrix proteins. MMPs are widely detected by zymography. The FITC-labeled Gelatin-Zymography Kit (Cat.No.PMC-AK83-COS) provides an easy system of the electrophoresis for zymography. This product is used for detecting ProMMP-2, MMP-2 and ProMMP-9 in blood, body fluid, secretion, cell lysate, cell culture medium, and other samples.

Components

Component	Quantity	Storage
Precast gel for 13 samples	5 pieces	4°C , Do not freeze
Electrophoresis Buffer (× 10)	100 mL × 2	4°C
Washing Buffer (× 10)	100 mL × 1	4°C
Reaction Buffer (× 10)	25 mL × 1	4°C
Sample Preparation Buffer	5 mL × 1	4°C
MMP markers (ProMMP-2, MMP-2, ProMMP-9)	0.2 mL	4°C

The size of the glass plate gel is 100 mm (W) × 100 mm (H) and the thickness is 1mm. One kit contains reagents for 60 samples.

Additional Materials Required

- Distilled water

Features

- Precast gel and reagents for zymography are included in this kit.
- You can check the enzyme reaction continuously.

Preparation of reagents

- Electrophoresis Buffer

Dilute Electrophoresis Buffer (× 10) with distilled water to a 1x concentrated solution.

Note: Diluted buffer should be stored at 4°C.

- Washing Buffer

Dilute Washing Buffer (× 10) with distilled water to a 1x concentrated solution.

Note: Diluted buffer should be stored at 4°C.

- Reaction Buffer

Dilute Reaction Buffer (× 10) with distilled water to a 1x concentrated solution.

Note: Diluted buffer should be stored at 4°C.

Preparation of samples

Mix Sample Preparation Buffer and a samples in the ratio of 1 to 1 and leave it in room temperature for 15 minutes.

Note: Do not heat the samples.

Protocol

1. Load 100-150 mL of Electrophoresis Buffer to the lower chamber (anode). Refer to the instruction manual of your electrophoresis tank/chamber to know appropriate volume of the buffer.
2. Take out the comb on precast gel carefully and set the gel to electrophoresis chamber. The side of sample holes should be set on the upper side. If the sample holes are disturbed, fix them up with needle etc.
3. Load around 100 mL of Electrophoresis Buffer to the upper chamber (cathode).
4. Apply samples and MMP markers to gel plate. (Fig.1: 10 μ L of MMP markers applied and enzymatic reacted for 24 hrs at 37°C).
5. Run electrophoresis at 15 mA constant current. (If you use 2 gels, set the current at 30 mA)
6. *MMP markers can be used without sample preparation buffer.
7. After the run is completed, turn off the electrophoresis chamber and take out the gel plate from the chamber.
8. Remove the upper glass plate and peel off the gel from the lower glass plate (the biggest glass plate) carefully with a spatula.
9. Put the gel in a tray with 200 mL of Washing Buffer. Incubate with shaking at Room Temperature.
10. Put the gel in a container with 50 mL of Reaction Buffer and seal up the container. Incubate the gel at 37°C for 20-40 hrs. (Lower enzyme concentration needs longer reaction time.)
11. After enzymatic reaction, remove the gel and take a photo by gel imaging system. (fluorescent wavelength: 535 nm)

COSMO BIO Co., LTD.

[JAPAN]

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,
KOTO-KU. TOKYO 135-0016, JAPAN
Phone: +81-3-5632-9610
FAX: +81-3-5632-9619
URL: <https://www.cosmobio.co.jp/>

COSMO BIO USA

[Outside Japan]

2792 Loker Ave West, Suite 101
Carlsbad, CA 92010, USA
e-mail: support@cosmobioua.com
Phone/FAX: (+1) 760-431-4600
URL: www.cosmobioua.com

Electrophoresis patterns


Fig 1. Electrophoresis pattern of MMP marker


Fig 2. Electrophoresis pattern for human tear samples

- ① Molecular-weight marker
- ② MMP marker
- ③~⑧ Human tear samples (6 lanes)

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