For research use only. Not for clinical diagnosis.

Arginase Activity Assay Kit

(Arginase Activity Assay Kit, Catalog No. AK89N)

August 9, 2024

Arginase is an enzyme that hydrolyzes L-Arginine to L-Ornithine and urea. It is known to be strongly expressed when macrophages are activated in M2 of the M1 (inflammatory)/M2 (non-inflammatory) known as one of the functional classifications of macrophages and is frequently used as a marker of M2 macrophage activation. This product extracts intracellular arginase and evaluates arginase activity by detecting urea produced by arginase from L-arginine. Compared to conventional assay methods, the stability and sensitivity of the reagent have been improved, and arginase activity can be easily measured.

This product is based on a patent owned by Professor Yasuhiko Tabata of Kyoto University and TOKYO OHKA KOGYO CO., LTD.

I-1. Kit Components

Contents	Volume	Number of bottles	Storage Temperature	Precautions for handling
Enzyme activation solution	2.5mL	1	4°C	When handling, wear protective equipment such as glasses and gloves, and take sufficient care to avoid contact with the human body.
Substrate solution	5mL	1		
Urea standard solution (160μg/mL)	1mL	1		
Urea Detection Solution A (light-shielded container)	2.2mL	1		Please check the label for hazardous indications. When handling, wear protective equipment such as glasses, mask, and gloves, and take sufficient care to avoid contact with the human body. Urea Detection Solution B reacts with ultraviolet rays, so be careful not to expose it to natural or fluorescent light for a long time.
Urea Detection Solution B (light-shielded container)	4mL	1		
Detection solution mixing vessel	-	2	Room temperature	Since the mixed urea detection solution reacts with ultraviolet light, it should be treated in the light-shielded container provided when mixed and heated.
Plate seal	-	2		Use the plate seal provided since the detection solution containing high concentrations of sulfuric acid is heated during the coloring reaction.

This product can measure 96 samples on a 96-well plate.

I-2. Preparations required] *Not included in the product.

- Gloves, mask, safety glasses
- Aluminum foil
- Microplate reader (measurement wavelength: 540 nm)
- Micropipettes, tips, microtubes with lids
- 96-well plate (transparent, flat bottom)
- Reservoir
- Thermostatic bath capable of heating up to 95°C and hot plate capable of heating up to 95°C (a heat block or thermal cycler the size of a 96-well plate can be substituted)
- 37°C thermostatic incubator or CO₂ incubator
- **Concentrated sulfuric acid (purity: 98%)**
- Cell extract: 0.1% Triton[®] X-100 aqueous solution plus cOmplete[™] protease inhibitor cocktail (part number: 11697498001, Roche Diagnostics, 1 tablet/50 mL)
- UV lamp for coloration (UV lamp attached to clean bench or UV gel imaging equipment can be substituted if available; if no UV lamp is available, natural or fluorescent light will also progress coloration, but it will take time)

II. Preparation for measurement Preparation before measurement

1. Urea Detection Solution A

Prepare 22mL of urea detection solution A by adding 9.9mL of concentrated sulfuric acid and 9.9mL of ultrapure water to the urea detection solution A bottle.

Note: Urea Detection Solution A after addition of sulfuric acid can be stored at 4°C for 1 year.

Note: When handling concentrated sulfuric acid, use appropriate protective equipment in a ventilated area. Follow the facility's regulations for handling and disposal.

2. Urea Detection Solution (Preparation for use)

Mix Urea Detection Solution A and Urea Detection Solution B at ratio 5:1 (A:B = 5:1) in the attached container for mixing detection solutions.

(180 µL of detection solution is required per well of a 96-well plate, so adjust the volume according to the number of wells to be used.)

Note: Since the mixed solution contains high concentrations of sulfuric acid and ethanol and the solution will be heated, be sure to check that the lid is tightly closed to prevent volatilization or leakage. Do not stir while heating. Use gloves, masks, protective glasses, etc. when handling.

After heating the solution in a thermostatic bath at 95°C for 30 minutes, allow it to cool sufficiently before stirring gently to use it as a urea detection solution.

(If there are insufficient mixing containers when multiple measurements are desired, prepare Urea Detection Solution A or B using a light-shielded container.)

3. Urea standard solution

Prepare a dilution series of Urea standard solution by 2-fold step dilution with ultrapure water.

Concentration: 0 (ultrapure water only), 2.5, 5, 10, 20, 40, 80, 160 µg/mL

Preparation volume: 100 µL each

Storage period of diluted urea standard solution: 4°C, 1 month

III. Sample Preparation

Cells are to be lysed by adding 40 µL of the above cell extract solution to a 96-well plate (bottom area: 0.32 cm²/well), excluding the supernatant, and the cells are to be used as the sample.

The volume of the extraction solution should be adjusted according to the bottom area of the plate.

IV. Procedure

1. Add 20 µL of enzyme activation solution and 40 µL of substrate solution to 20 µL of sample and urea standard solution and mix.

*Duplicate assay is recommended.

2. Cover the plate with the attached plate seal or lid and incubate in a 37°C incubator for 1 hour to allow the enzyme-substrate reaction to take place.

3. After the enzyme-substrate reaction, add urea detection solution prepared in advance at 180 µL/well.

*Please be careful to avoid exposure to light as much as possible when adding the urea detection solution, for example by using a reservoir to add it quickly.

4. Seal the plate tightly with the attached plate sealer and shield the 96-well plate from light with aluminum foil or the like.

Note: Since the urea detection solution contains high concentrations of sulfuric acid and ethanol and will be heated, the plate seal should be tightly adhered to prevent volatilization or leakage. Use gloves, masks, protective glasses, etc. when handling.

5. Heat the product at 95°C for 2 hours under shaded light using a hot plate, etc.

6. After heating, remove the aluminum foil and accelerate the reaction under a UV lamp until sufficient color is observed. 15 minutes to 1 hour under a UV lamp is sufficient. Under natural light or fluorescent light, it takes about 1 to 2 hours. Refer to the example in Figure 1 for the degree of coloration.

7. Measure absorbance (wavelength: 540 nm)

Urea concentration

160 μ g/mL

80 μ g/mL

40 μ g/mL

20 μ g/mL

10 μ g/mL

5 μ g/mL

2.5 μ g/mL

0 μ g/mL


Fig.1: Example of coloration of urea standard solution

V. Calculation of standard curve and arginase activity

The absorbance of the urea standard solution at each concentration is averaged, and a standard curve is calculated from the average value by fitting a four-parameter logistic curve fit (4-PL) or a quadratic regression curve (either function will fit).


Fig.2: Example of a 4PL curve (software: Prism8)
Excel)


Fig.3: Example of a regression curve (software:

Arginase activity (units/mL) is defined as the conversion of 1 μmol of L-Arginine to 1 μmol of L-Ornithine and urea per minute at 37°C, which can be obtained using the following formula

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