**For research use only. Not for clinical diagnosis.**

# Bone Marrow Monocyte Precursor Cell

Catalog No. BMMC (Japan)  
PMC-BMMC (Outside Japan)

April 3, 2026

## I. Product Overview

Monocytes are produced by the differentiation of monocyte precursor cell in bone marrow. In tissue, monocytes mature into macrophages and are involved in immunity, recovering tissues, etc. Bone Marrow Monocyte Precursor Cell (Code: BMMC) contains frozen monocyte precursor cells derived from rat bone marrow, and the cells can be differentiated into mature monocytes using Wash Medium (Code: BMMW) and Culture Medium (Code: BMMG).

## II. Product Warranty

**Cosmo Bio warrants that product (cells) will be viable until the expiration date, and is valid only if the product is stored and cultured according to the information indicated in this product data sheet. Cosmo Bio has optimized the cell culture media formulation which is ideal for the product. While other, unspecified cell culture media may also produce satisfactory results, a change in cell culture media or the absence of an additive(s) from the recommended cell culture media may affect recovery, growth and/or function of the product. If an alternative cell culture medium formulation is used to culture the product, the Cosmo Bio warranty for cell viability is no longer valid.**

## III. Product Components

Product Name	Size	Quantity	Storage Conditions	Expiration date
Bone Marrow Monocyte Precursor Cell (Catalog No. BMMC)	2 x 10 <sup>6</sup> cells/vial	1	Liquid Nitrogen	6 months

\*Shipping: dry ice

## IV. Dedicated Media (Serum Containing Media) (sold separately)

Product Name	Size	Quantity	Storage Conditions	Expiration date
Monocyte Wash Medium (Catalog No. BMMW)	50 mL	1	-20 °C (Store at 4 °C after thawing)	- Written on the bottle (stored at -20 °C) - Use within 3 months after thawing. (stored at 4 °C)
Monocyte Culture Medium (Catalog No. BMMG)	25 mL	1		- Written on the bottle (stored at -20 °C) Use within 2 weeks after thawing. (stored at 4 °C)

\*Shipping: dry ice

### Culture Medium components:

Catalog No. BMMW : MEM-α, FBS, antibiotic, etc.

Catalog No. BMMG : MEM-α, FBS, antibiotic, 100 ng/mL M-CSF, etc.

**For research use only. Not for clinical diagnosis.**

※ The culture medium can be refrozen only once. It is highly recommended to aliquot out to avoid freeze-thaw cycles.

It should be noted that repeated freeze-thaw of these products will damage or cause quality loss of these products.

## V. Instructions For Use

※ This product [is not capable of being subcultured].

## VI. Materials required but not provided

- Variable volume pipettes
- Monocyte Wash Medium (Catalog No. BMMW)
- Monocyte Culture Medium (Catalog No. BMMG)
- HydroCell (CellSeed Inc.) Culture plate, 12-well
- 15 mL centrifuge tube.

## VII. Precautions

- Read the instructions carefully before beginning the culture.
- This kit is for research use only, not for human or diagnostic use.
- Always wear gloves and lab coat when handling the cell culture.


## VIII. Protocols

Prepare a collagen type I coated plate.

1. Thaw the Monocyte Wash Medium (Catalog No. BMMW) and Monocyte Culture Medium (Catalog No. BMMG) at 4 °C.
2. Quickly place Monocyte Precursor Cell vial in a 37 °C water bath until its content is completely thawed. (for 2 minutes)
3. Transfer thawed cells into a 15 mL centrifuge tube containing 10 mL of Wash Medium and centrifuge at 4°C at 430 x g for 5 minutes to pellet the cells.
4. Remove the supernatant, and re-suspend cells in 10 mL of Wash Medium and centrifuge at 4 °C at 430 x g for 5 minutes.
5. Remove the supernatant, and re-suspend the cell pellet in a total volume of 5.5 mL of Culture Medium.
6. Transfer the cell suspension to culture plate (1 mL/well/12well-plate) and incubate at 37 °C in a 5 % CO<sub>2</sub>, humidified incubator.
  - ※ Low cell binding culture plates such as HydroCell (CellSeed Inc.), are recommended to culture monocytes as the cells are capable of strongly adhering to normal culture plates.
  - ※ The size of cells is about several micrometers in the first cell culture and may reach 10 micrometers in 3-4 days of culture (Figure 1). When cells are cultured for a period of time they are able to form several tens of colonies.

**For research use only. Not for clinical diagnosis.**

**Figure 1**

## **IX. References**

- (1) Sunao Takesita, Keisuke Kaji, Akira Kudo. Identification and Characterization of the New Osteoclast Progenitor with Macrophage Phenotypes Being Able to Differentiate into Mature Osteoclasts. JOURNAL OF BONE AND MINERAL RESEARCH Volume 15, Number 8(2000) .1477-1488.

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