Product manual

For research use only. Not for clinical diagnosis.

Catalog No. CMC12C

Cardiomyocyte (Mouse)

May. 20, 2019

Principle

Cardiomyocyte is one of the cell groups that compose heart. Cardiomyocyte is known as beating involuntary striated muscle cells. The beat that is the feature of the cardiomyocyte can be widely used for a pharmacologic and electrophysiology assays, because it responds to a variety of stimuli including hormones, drugs and electricity. Cells are isolated from mouse embryos' heart, and cardiomyocytes are enriched by removal of non-cardiomyocytes. The involuntary beat of the cells can be seen in culture using Cardiomyocyte (Catalog No. CMC12C) .

Warranty

Cosmo Bio warrants that product (cells) will be viable until the expiration date, and is valid only if the product is stored and cultured according to the information indicated in this product data sheet. Cosmo Bio has optimized the cell culture media formulation which is ideal for the product. While other, unspecified cell culture media may also produce satisfactory results, a change in cell culture media or the absence of an additive(s) from the recommended cell culture media may affect recovery, growth and/or function of the product. If an alternative cell culture medium formulation is used to culture the product, the Cosmo Bio warranty for cell viability is no longer valid.

Components/Storage

Product Name	Size	Quantity	Storage Conditions	Expiration date
Cardiomyocyte (Mouse) (Catalog No. CMC12C)	2×10 ⁶ cells/vial (0.5mL/vial)	1	Liquid Nitrogen	6 months

*Shipping: dry ice

Optimized Culture Medium (selling separately)

Product Name	Size	Quantity	Storage Conditions	Expiration date
Cardiomyocyte Culture Medium (Catalog No. CMCM)	500 mL	1	-20°C (Store at 4°C after thawing)	- Written on the bottle (stored at -20°C)
Fibronectin coat solution (Catalog No. SFN01)	12 mL	1		- Please use up early after thawing (stored at 4°C)

Culture Medium components:

Catalog No. CMCM : D-MEM/F-12, FBS, antibiotic, etc.

Catalog No. SFN01 : bovine, plasma derived, concentration 20 µg / mL

*Shipping: dry ice

Materials required but not provided

- Variable volume pipettes
- Cardiomyocyte Culture Medium (Catalog No. CMCM)
- Fibronectin coat solution (Catalog No. SFN01)
- 24-well culture plate (flat bottom)
- 15 ml centrifuge tube.

Precautions

- Read the instructions carefully before beginning the culture.
- This kit is for research use only, not for human or diagnostic use.
- Always wear gloves and lab coat when handling the cell culture.
- Upon receiving, cryovials are immediately transferred to and stored in liquid nitrogen until ready to use. It is recommended to use cell vial immediately upon receiving.

Protocols

< Prepare a fibronectin coated Culture 24-well plate >

【Things to prepare】

- 24 well plate for culture
- Fibronectin coat solution (Catalog No. SFN01)
- Sterilized ultra pure water

1. Thaw fibronectin solution in a refrigerator (4 ° C) in advance.
2. Add 0.5 ml of fibronectin solution to each well.
(Add 0.25 ml/cm² of solution to each well in other types of culture plates.)
3. Leave the plate at 37°C in 5% CO₂ humidified incubation overnight.
4. Aspirate remaining material and rinse plates with sterile ultra pure water twice. Plates are ready for use.

<1-1. Cultured with the 24well-plate>

A seeding density of 2.0×10⁵ cells/well is recommended.

1. Thaw the Cardiomyocyte Culture Medium (cat# CMCM) at 4 ° C .
2. 4.5 mL of Cardiomyocyte Culture Medium into a 15 ml centrifuge tube and warm the to 37 ° C.
3. Quickly thaw the Cardiomyocyte vial for 1min 15sec in a 37°C water bath.
4. Transfer the vial contents of thawed cells into a 15 ml centrifuge tube containing 4.5 ml of Culture Medium.
5. Gently mix gently the cell suspension by slow pipetting up and down, and adjust cell density to 2.0×10⁶ cells/5 ml solution in the tube.
6. Transfer 0.5 ml of cell suspension to each well of fibronectin-coated 24-well plate.
7. Incubate the plate at 37°C in a 5 % CO₂ humidified incubator.
8. The next day, gently add 0.5 ml of fresh pre-warmed culture medium to each well.
9. Exchange the medium with fresh and pre-warmed culture medium every day or every other day until the culture reaches 80-100% of cell confluent.

* Culture reaches its confluent within 3-5 days, and cardiomyocytes begin beating.

* Do NOT use cold culture medium. Please pre-warm culture medium at 37°C before use to ensure the viability of cardiomyocytes.

* If you incubate at a volume of 0.5 mL or more, the beating rate may drop or the beating may stop.


Fig. 1 Morphology of cardiomyocyte

<1-2. Culture adding 5-Bromo-2'-deoxyuridine>

It is recommended to supplement 5-Bromo-2'-deoxyuridine (BrdU) to culture in order to increase the purity of cardiomyocytes (BrdU prevents non-cardiomyocytes from developing). Preparation of BrdU: Dissolve BrdU (Sigma Cat. No. B-9285 or the similar product) in culture medium to a concentration of 10mM and filter sterilization. Fresh BrdU solution should be used for cardiomyocyte culture.

10. Prepare cell suspension as described above in steps 1-5(<1-1. Cultured with the 24well-plate>).
11. Transfer 0.5 ml of cell suspension to each well of fibronectin-coated 24-well plate and add 5 µl of 10mM BrdU to each well (0.1 mM BrdU-final concentration).
12. Incubate the plate at 37°C in a 5 % CO₂ humidified incubator.
13. The next day, exchange the medium with fresh 0.5 ml of the medium containing 0.1 mM BrdU.
14. On day 3 of the culture, replace the medium with fresh 0.5 ml of the medium without BrdU.

For research use only. Not for clinical diagnosis.

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COSMO BIO Co., LTD.