# Osteoblast from cranial bone

(Rat Osteoblast from cranial bone)

May 14, 2019

## Principle

Bone metabolism is composed of balanced bone formation of osteoblasts and bone resorption of osteoclasts. Osteoblast from cranial bone (OBC02C) contains frozen osteoblasts isolated from rat calvariae. Osteoblast from cranial bone (OBC02C) is a useful to study osteoblast and osteogenesis.

## Warranty

Cosmo Bio warrants that product (cells) will be viable until the expiration date, and is valid only if the product is stored and cultured according to the information indicated in this product data sheet. Cosmo Bio has optimized the cell culture media formulation which is ideal for the product. While other, unspecified cell culture media may also produce satisfactory results, a change in cell culture media or the absence of an additive(s) from the recommended cell culture media may affect recovery, growth and/or function of the product. If an alternative cell culture medium formulation is used to culture the product, the Cosmo Bio warranty for cell viability is no longer valid.

## Precautions

- Read the instructions carefully before beginning the culture.
- Remove the cryovial from the dry ice packaging and immediately place into liquid nitrogen storage until use.

## Components

| Product Name              | Size                              | Quantity | Storage Conditions             | Expiration date |
|---------------------------|-----------------------------------|----------|--------------------------------|-----------------|
| Osteoblast, cryopreserved | 1.0x10 <sup>6</sup><br>cells/vial | 1        | Liquid Nitrogen<br>vapor phase | 6 months        |

\*Shipping: dry ice

### Optimized Culture Medium (selling separately)

| Catalog No. | Product Name              | Size   | Quantity | Storage Conditions | Expiration date         |
|-------------|---------------------------|--------|----------|--------------------|-------------------------|
| OBCM        | Osteoblast Culture Medium | 500 mL | 1        | -20°C<br>Freezer   | - Written on the bottle |
|             |                           |        |          | 4°C                | 3month after thawing    |

Culture Medium components:  $\alpha$ -MEM, FBS, antibiotics, etc.

### Materials required but not provided

- ✧ Variable volume pipettes
- ✧ Culture vessels
- ✧ 0.05% Trypsin/EDTA
- ✧ HBSS or PBS(-)

### Protocol

- A) Cultured with the 25 cm<sup>2</sup> flask
- 1) Carefully remove the cryovial from liquid nitrogen and thaw cells in a water bath at 37°C for 120 seconds.
  - 2) Transfer thawed cells into a 15 ml centrifuge tube containing 10 ml of culture medium and Transfer 1mL of culture medium in the same conical tube back to the cryovial and pour the contents back to 15mL conical tube.
  - 3) Centrifuge for 5 minutes at 4°C at 600 x g for 5 minutes.
  - 4) Remove the supernatant, and re-suspend the cell pellet in approximately 5 ml of culture medium.
  - 5) Transfer the cell suspension to 25 cm<sup>2</sup> flask and incubate the flask at 37°C under 5% CO<sub>2</sub> and 100% humidity.
  - 6) The next day, Replace the medium with fresh pre-warmed culture medium.  
※Approximately 2-3 days of culture, cells become confluent. For subculture, please refer to the protocol below. Subculture of the cells can be performed up to passage 2.

B) Subculture

- 1) When the cells reach 70 - 90% of confluent, they should be subcultured.
  - 2) Aspirate the medium. Rinse the dish with 10mL of HBSS or PBS (-). Repeat twice.
  - 3) Remove HBSS or PBS (-) and then add 3 ml of trypsin/EDTA solution into flask (25 cm<sub>2</sub> flask).
  - 4) Gently rock the flask to make sure that the cells are covered by trypsin/EDTA solution and then immediately remove trypsin/EDTA solution.
  - 5) Incubate the flask in a 37°C incubator until cells are completely rounded up (monitored with inverted microscope). Approximately it takes 2 to 3 minutes.
  - 6) Add culture medium to the flask and transfer detached cells to centrifuge tube, and then centrifuge the centrifuge tube at 4°C at 600 x g for 5 minutes.
  - 7) After removing the supernatant, re-suspend cells in culture medium and centrifuge .for 5 minutes at 4°C at 600 x g.
  - 8) Remove the supernatant, and re-suspend cells in culture medium. Count cells and plate cells in a new plate or flask (Adjust cell density to the desired experiment).
- ※Approximately 2-3 days of culture, cells become confluent when seeding density is 30,000 cells/ cm<sup>2</sup>.

COSMO BIO Co., LTD.

[JAPAN]

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: 81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>

COSMO BIO USA

[Outside Japan]

2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobioussa.com](mailto:info@cosmobioussa.com)  
Phone/FAX: (+1) 760-431-4600  
URL: [www.cosmobioussa.com](http://www.cosmobioussa.com)

**For research use only. Not for clinical diagnosis.**

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