**For research use only. Not for clinical diagnosis.**

# Osteoclasts V-2, Rat

Catalog No. OSC12C

April 15, 2020

## Principle

Bone metabolism is composed of balanced osteogenesis and bone resorption. Research studies have shown that bone marrow cells can be differentiated into osteoclasts using M-CSF (Macrophage Colony Stimulating Factor) and RANKL (Receptor Activator of NF kappa B Ligand). #OSC12C is useful to evaluate osteoclast formation and activation.

## Warranty

Cosmo Bio warrants that product (cells) will be viable until the expiration date, and is valid only if the product is stored and cultured according to the information indicated in this product data sheet. Cosmo Bio has optimized the cell culture media formulation which is ideal for the product. While other, unspecified cell culture media may also produce satisfactory results, a change in cell culture media or the absence of an additive(s) from the recommended cell culture media may affect recovery, growth and/or function of the product. If an alternative cell culture medium formulation is used to culture the product, the Cosmo Bio warranty for cell viability is no longer valid.

## Precautions

- Read the instructions carefully before beginning the culture.
- Remove the cryovial from the dry ice packaging and immediately place into liquid nitrogen storage until use.

## Components

Product Name	Size	Quantity	Storage Conditions	Expiration date
Osteoclast precursor, Rat, cryopreserved	2.0x10 <sup>6</sup> cells/vial	2 vials	Liquid Nitrogen vapor phase	6 months

**\*Shipping: dry ice**

## Optimized Medium (selling separately)

Catalog No.	Product Name	Size	Quantity	Storage Conditions	Expiration date
OSCMW	Osteoclast Wash Medium	100 mL	1	-20°C Freezer	- Written on the bottle
OSCMR	Rat Osteoclast Culture Medium (w/ RANKL and M-CSF)	50 mL	1	-20°C Freezer	- Written on the bottle

**Culture Medium components: α-MEM, FBS, antibiotics, etc.**

## Materials required but not provided

- ✧ Osteoclast Wash Medium (#OSCMW)
- ✧ Osteoclast Culture Medium (#OSCMR)
- ✧ Variable volume pipettes
- ✧ Culture vessels
- ✧ 15 ml centrifuge tube

## Protocol

- A) Cultured with the 96-well culture plate
- 1) Carefully remove the cryovial from liquid nitrogen and thaw cells in a water bath at 37°C with gently shaking.
  - 2) Transfer thawed cells into a 15 ml centrifuge tube containing 10 ml of Osteoclast Wash Medium (#OSCMW) and Transfer 1mL of culture medium in the same conical tube back to the cryovial and pour the contents back to 15mL conical tube.
  - 3) Centrifuge for 5 minutes at 4°C at 200 x g for 5 minutes.
  - 4) Remove the supernatant, and re-suspend the cell pellet in 2.5 to 5.0 ml of Rat Osteoclast Culture Medium (#OSCMR).
  - 5) Dispense 100  $\mu$  L of the cell suspension to each well of 96-well culture plate, and incubate the flask at 37°C under 5% CO<sub>2</sub> and 100% humidity.
  - 6) Three days later, Replace the medium with fresh culture medium.
  - 7) Change the medium every other day.
- ※Osteoclasts will begin to fuse and form osteoclasts after 4 or 5 days of incubation

For study of osteogenesis, add the reagent to the medium at various stages of osteogenesis. To visualize the osteogenesis, the optional TRAP Staining Kit (#AK04F) is useful for staining of tartrate-resistant acid phosphatase.


Fig.1 Osteoclast induced of differentiation with M-CSF/RANKL


Fig. 2 Pit on the slice of ivory.  
(HE staining)


Fig. 3 Pit on the slice of ivory.  
(SEM picture)

## References

1. Takeshita et al.(2000) Identification and Characterization of the New Osteoclast Progenitor with acrophage Phenotypes Being Able to Differentiate into Mature Osteoclasts. J. Bone Miner. Res. 15, 1477-1488.
2. Scheven et al. (1998) A sequential culture approach to study osteoclast differentiation from nonadherent porcine bone marrow cells. In Vitro Cell Dev Biol Anim 34, 568–577.
3. Martha et al. (1995) Enhanced Expression of  $\alpha$ V Integrin Subunit and Osteopontin during Differentiation of HL-60 Cells along the Monocytic Pathway. Exp. Cell Res. 216, 335-341.
4. Itonaga et al.(1999) 1,25-Dihydroxyvitamin D3 and Prostaglandin E2 Act Directly on Circulating Human Osteoclast Precursors. Biochem. Biophys. Res. Commun. 242, 703-709.

**COSMO BIO Co., LTD.**

**【JAPAN】**

TOYO EKIMAE BLDG. 2-20, TOYO 2-CHOME,  
KOTO-KU. TOKYO 135-0016, JAPAN  
Phone: +81-3-5632-9610  
FAX: +81-3-5632-9619  
URL: <https://www.cosmobio.co.jp/>

**COSMO BIO USA**

**【Outside Japan】**

2792 Loker Ave West, Suite 101  
Carlsbad, CA 92010, USA  
email: [info@cosmobiousa.com](mailto:info@cosmobiousa.com)  
Phone/FAX: (+1) 760-431-4600  
URL: [www.cosmobiousa.com](http://www.cosmobiousa.com)

**For research use only. Not for clinical diagnosis.**

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- (1) Takeshita *et al.*(2000) Identification and Characterization of the New Osteoclast Progenitor with Macrophage Phenotypes Being Able to Differentiate into Mature Osteoclasts. *J. Bone Miner. Res.* 15, 1477-1488.
- (2) Scheven *et al.* (1998) A sequential culture approach to study osteoclast differentiation from nonadherent porcine bone marrow cells. *In Vitro Cell Dev Biol Anim* 34, 568-577.
- (3) Martha *et al.* (1995) Enhanced Expression of  $\alpha$ v Integrin Subunit and Osteopontin during Differentiation of HL-60 Cells along the Monocytic Pathway. *Exp. Cell Res.* 216, 335-341.
- (4) Itonaga *et al.*(1999) 1,25-Dihydroxyvitamin D<sub>3</sub> and Prostaglandin E<sub>2</sub> Act Directly on Circulating Human Osteoclast Precursors. *Biochem. Biophys. Res. Commun.* 242, 703-709.

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COSMO BIO CO., LTD.