Rat Subcutaneous White Adipocyte

Principle

White adipocytes, a type of subcutaneous adipocytes, play an important role in energy storage. White Adipocyte Culture Kit, V-1 (Rat) contains preadipocytes isolated from rat subcutaneous adipose tissues and culture medium that induces differentiation of precursor cells into mature adipocytes. The kit provides a convenient system for studying the mechanism of adipogenesis as well as for analyzing effects of drugs on metabolic syndrome such as obesity, diabetes and hypertension.

Warranty

Cosmo Bio warrants that product (cells) will be viable until the expiration date, and is valid only if the product is stored and cultured according to the information indicated in this product data sheet. Cosmo Bio has optimized the cell culture media formulation which is ideal for the product. While other, unspecified cell culture media may also produce satisfactory results, a change in cell culture media or the absence of an additive(s) from the recommended cell culture media may affect recovery, growth and/or function of the product. If an alternative cell culture medium formulation is used to culture the product, the Cosmo Bio warranty for cell viability is no longer valid.

Precautions

- Read the instructions carefully before beginning the culture.
- Remove the cryovial from the dry ice packaging and immediately place into liquid nitrogen storage until use.

Components

Product Name	Size	Quantity	Storage Conditions	Expiration date
Rat Subcutaneous White Adipocyte, Cryopreserved	3 x 10 ⁶ cells/vial	1 vial	Liquid nitrogen storage	6months

***Shipping: dry ice**

Optimized Culture Medium (selling separately)

Catalog No.	Product name	Size	Storage Conditions	Expiration date
SACMR	Subcutaneous Adipocyte Culture Medium	250mL	-20°C Freezer	- Written on the bottle

Materials required but not provided

- ✧ Variable volume pipettes
- ✧ 15 ml centrifuge tube
- ✧ 24-well culture plate
- ✧ Subcutaneous Adipocyte Culture Medium (Catalog No. SACMR)

Protocols

1. Thaw the Subcutaneous Adipocyte Culture Medium (Catalog No. SACMR) at 4°C
2. Thaw vials of frozen cells by warming them in 37°C warm water for 2 minutes.
3. Transfer thawed cell solution to a 15 mL centrifuge tube containing 10 mL of Subcutaneous Adipocyte Culture Medium (Catalog No. SACMR). Mix gently and centrifuge at 4°C at 200 x g for 5 minutes.
4. Remove the supernatant and re-suspend cells in 10 mL of Subcutaneous Adipocyte Culture Medium (Catalog No. SACMR) and centrifuge at 4°C at 200 x g for 5 minutes.
5. Remove the supernatant and re-suspend the cell pellet in 12.5 mL of Subcutaneous Adipocyte Culture Medium (Catalog No. SACMR)
6. Dispense 0.5 mL of cell suspension to each well of 24-well plate and incubate the plate at 37°C in a 5% CO₂ humidified incubator.
7. Gently add 0.5 ml of Subcutaneous Adipocyte Culture Medium (Catalog No. SACMR) to each well of 24-well plate on the following day.
8. Carefully change the medium 1 ml per well on the second day after seeding.
9. Change the medium every other day. Be careful to not disturb the cell layer.
 - i. Approximately 3-4 days of culture, preadipocyte culture becomes confluent.
 - ii. Approximately 5 days of culture, cells turn to mature adipocytes.
 - iii. Approximately 8 days of culture, cells become hypertrophic and start detaching from the bottom of the well.

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COSMO BIO CO., LTD.