Attention (for increase in binding capacity)

"Standard Protocol" is suited for purification of "high-binding-capacity" antibodies such as human or mouse loG. "Saturation Protocol". described below, is recommended for such as rat IgG or "low-binding-capacity" antibodies or first trial of the item.

Saturation Protocol Step 1. Equilibration

Note: If IgG-binding is too weak, use 1.5 M Glycine-3.0 M NaCl (pH 9.0) as binding buffer, after confirmation that antibody is not ·Same as "Standard Protocol". inactivated.

Step 2. Sample Apply

•Plug into outlet of a column tightly and add prepared sample to the column.

• Close a screw cap and shake the column vigorously for 1-2 hrs to avoid sinking gel.

• Put off the outlet plug, set the column into a 2mL-tube and centrifuge at 2,000 x g for 5 seconds.

Step 3. Wash

•Put off the screw cap, add 0.6 mL of Binding Buffer, plug into column-outlet and shake for 5 min. • Put off the outlet plug, set the column into a 2 mL-tube and centrifuge at 2,000 x g for 5 seconds. •Repeat these steps for more 2 times.

Step 4. IaG-Elution

- Plug into outlet of a column tightly and add 0.2 mL of Elution Buffer.
- Close a screw cap and shake the column vigorously for 5 min to avoid the gel sinking.
- Put off the outlet plug, set the column into a 1.5 mL-tube containing Neutralization Buffer. Eluate is collected into the tube after centrifugation at 2,000 x g for 5 seconds.
- Set the column into another 1.5 mL-tube containing Neutralization Buffer and repeat the same procedure. (If pH of Elution Buffer is higher than 3, repeat once again and collect 3rd eluate.)

Note: • In some species of antibody, binding to the gel may be weak.

In some molecular species of Rat IgG2a, binding to the gel may be weak (EX : about 1mg/mL gel) ·In mouse IgM, there are 2 type of molecular species. "High-binding" type can be purified with this protocol, but "low-binding" type is difficult to be purified.

Order Information

Product Name	Contents	Code No.
Ab-Capcher ExTra	2 mL	P-003-2
	10 mL	P-003-10
	100 mL	P-003-100
Ab-Rapid SPiN EX (Spin column)	0.1 mL gel/column x 10 (20-empty 2 mL-tubes included)	P-014-10
Buffer Kit	1 kit (Bind. 200mL, Elut. 30mL, Neutr. 1mL)	P-011

There are cases that prices will be changed without notice. For research use only.

COSMO BIO CO., LTD. Inspiration for Life Science TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME KOTO-KU, TOKYO 135-0016, JAPAN TEL:+81-3-5632-9617 FAX:+81-3-5632-9618

ProteNova Co.Ltd. Takamatsu Lab. 2217-44 Hayashi-cho, Takamatsu Kagawa 761-0301 Japan TEL +81-87-897-2073 FAX +81-87-816-2073 URL http://protenova.com




Ab-Rapid SPiN EX

Users Manual P-014-10 Ver.1.2

Ab-Rapid SPiN EX Specifications

This column includes endotoxin-tested Ab-Capcher ExTra.

•Gel volume: •Gel matrix:	0.1 mL (25% gel slurry, 0.4mL) Highly crosslinked-agarose (Rapid Run)	reverse side for closing
•Column volume: •Particle size: •Ligand:	0.8 mL 35 μm Alkali-resistant Protein A derivative (Protein A-R28)	Screw cap
 Binding Capacity: Form: Storage: Accessories: 	approx. 7.5 mg human IgG /column 20% ethanol 4-8°C 20 empty 2 mL-tubes	← Lower filter ← Snap-off plug



How to use Snap-off plug

Snap off the outlet plug and use its

Materials

- bench-top centrifuge (1,000 2,000 × g)
- •1.5mL micro-centrifuge tube
- Buffers
 - Binding Buffer: PBS Elution Buffer: 0.1 M Glycine-HCl, pH 2.5 - pH 3.0 Neutralization Buffer: 1 M Tris
- * Buffer Kit (Set of buffers needed for antibody purification) is available. (See Order Information)

Sample preparation (example)

- ◆ Ascites :3 x dilution with Binding Buffer.
- \bullet Serum: precipitation with 50%-saturated (NH₄)₂SO₄ or 5x dil. with Binding Buffer
- Cultured medium: Adjust pH to neutral.

Recommended pre-treatments of sample before applying to the column.

- Centrifugation; $10,000 \times q$, 10 min
- Filtration; 0.45µm filter

(Please use low-protein-adsorption types)

* If there are insolubles in the sample, make sure to do pre-treatments.

Preparation for 50% ammonium sulfate precipitation

- 1. Prepare saturated ammonium sulfate.
- Add equal volume of saturated ammonium sulfate gradually to serum and mix.
- 2. Stand on ice for more than 1hr.
- 3. After centrifugation at 4°C, remove the supernatant.
- Wash precipitate with 50%-saturated ammonium sulfate.
- 4. Resolve the precipitate with small volume of Binding Buffer. The precipitate contains antibody.
- 5. Exchange to Binding buffer with dialysis or desalting column.

"Standard Protocol"	Required Time : 10) min	ſ
Preparation of tubes for Elution Add Neutralization Buffer to 1.5mL r Elution Buffer(pH 2.5) · · · 1st tube, Elution Buffer(pH 2.8) · · · 1st tube,	5 μL; 2nd tube, 9 μL		Neutralization Buffer
Step 1. Equilibration •If gel attached to top of column, mix •Put off an outlet plug and set a colu •Remove preservative solution by ca (If lid prevent centrifuge, cut off lids •Put off a screw cap, add 0.6 mL of for 5 seconds. If buffer remains in the	imn into a 2mL-tube (included). entrifugation at 2,000 x g for 5 sec of tubes before centrifugation.) Binding Buffer and centrifuge at 2	2,000 x g	1 min
Step 2. Sample Apply •Plug into outlet of the column tightly •Close the screw cap and incubate f 30-60 seconds. •Put off the outlet plug, set the colur centrifuge at 2,000 x g for 5 seconds	or 4 min with mixing every nn into a 2 mL-tube and	column. Sample→	4 min
Step 3. Wash •Put off the screw cap, add 0.6 mL suspended, and centrifuge at 2 •Repeat this step more 2 times in b (If non-specific proteins should be referred)	2,000 x g for 5 seconds. ench-top centrifuge.	S Was	1 min
Step 4. IgG-Elution Plug into outlet of the column tightl Close the screw cap, mix by tappir Put off the outlet plug, set the colur Neutralization Buffer and collect elu at 2,000 x g for 5 seconds. Repeat the same steps, collect 2nd including Neutralization Buffer. (80% of purified IgG is collected in If higher concentration of IgG is nee 2nd eluate is also available.)	g and leave to stand for 1 min. nn into a 1.5 mL-tube including ate in the tube by centrifugation eluate in another 1.5 mL-tube 1st eluate and 20% of it is in 2nd o	Elution → Buffer 0.2 mL Ig eluate. Ne	1 min × 2
For increase in binding cap "Saturation Protocol"	acity) Required T		e Next Page °s
 Incubation time of sample and ge In Washing and Elution (Step 3 a 			teps.