

A solution of human IgG solution (30mg/mL, 8 mL) was applied to the column and IgG was purified according to "Standard Protocol". Total 191.8 mg of IgG was purified with 4-times elution. 1st fraction eluted with 2 bed volumes (5 mL) contained 37% of total purified IgG. When 2nd fraction eluted with 1 bed volume (2.5 mL) was added to 1st fraction, it was increased to 79.9%. When 3rd fraction eluted with 1 bed volume (2.5 mL) was added to the mixture of 1st and 2nd fractions, 97.7% of total purified IgG was recovered.

Order Information

Product Name	Contents	Code No.		
Ab-Capcher ExTra	2 mL	P-003-2		
	10 mL	P-003-10		
	100 mL	P-003-100		
Ab-Rapid SPiN EX (Spin column)	0.1 mL gel/column x 10 (20 empty 2-mL tubes included)	P-014-10		
	2.5 mL gel x 2, empty column x 2	P-014-5-1		
Buffer Kit	1 kit (Bind. 200mL, Elut. 30mL, Neutr. 1mL)	P-011		
There are cases that prices will be				

changed without notice. For research use only.



TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME KOTO-KU, TOKYO 135-0016, JAPAN TEL : +81-3-5632-9617 FAX : +81-3-5632-9618 ProteNova Co.Ltd. Takamatsu Lab. 2217-44 Hayashi-cho, Takamatsu Kagawa 761-0301 Japan TEL +81-87-897-2073 FAX +81-87-816-2073 URL http://protenova.com

Ab-Rapid SPiN EX-5mL

P-014-5-1 Ver.1.2



Ab-Rapid SPiN EX 5mL specifications

•Gel volume:	Ab-Capcher ExTra, 2.5mL x 2 (50% gel slurry, 5mL x 2)
 Spin column: 	2 (Max. 22 mL / column)
 Top cap 	2
 Bottom cap: 	2
 Gel matrix: 	Highly closslinked agarose
 Particle size: 	approx. 35µm
Ligand:	Alkali-resistant Protein A derivative (Protein A-R28)
 Binding capacity : (saturated) 	approx. 190 mg human IgG / column
• Form :	20% Ethanol
 Storage: 	4-8°C

How to use of Snap-Off plug Before use of the column, snap-off the plug and use the bottom cap. Top cap Bottom cap Lower filter Plug

Materials

- Centrifuge (swing or angle type, < 3,000 rpm, the column can be used at 3,000 × g)
 50mL conical tube
- Buffers

Binding Buffer: PBS Elution Buffer: 0.1 M Glycine-HCl, pH 2.5 - pH 3.0 Neutralization Buffer: 1 M Tris

* Buffer Kit (Set of buffers needed for antibody purification) is on sale. (See Order Information)

Sample preparation (example)

- Ascites :3 x dilution with Binding Buffer.
- \bullet Serum : precipitation with 50%-saturated (NH₄)₂SO₄ or 5 x dilution with Binding Buffe
- Cultured medium: Adjust pH to neutral.

Recommended pre-treatments of sample before applying to column.

- Centrifugation; 10,000 \times g, 10 min
- Filtration; 0.45µm filter

(Please use low-protein-adsorption types)

* If there are insolubles in the sample, make sure to do pre-treatments.

Preparation for 50% ammonium sulfate precipitation

1. Prepare saturated ammonium sulfate.

Add equal volume of saturated ammonium sulfate gradually to serum and mix.

2. Stand on ice for more than 1hr.

3. After centrifugation at 4°C, remove the supernatant.

Wash precipitate with 50%-saturated ammonium sulfate.

4. Resolve the precipitate with small volume of Binding Buffer. The precipitate contains antibody.

5. Exchange to Binding buffer with dialysis or desalting column.

olug	Protocol Requir	red Time : 2 hrs	
off the	Preparation of column and recovery tubes for Elution •Mix a bottle of Ab-Capcher ExTra 2.5mL and add all •Add Neutralization Buffer to 50 mL conical tubes. •••1st tube, 100μL; 2nd or later tubes, 50 μL (In case	l of content to an empty column. se of Elution Buffer at pH 2.8)	
D	Step 1. Equilibration •Snap off the outlet plug of the column and set it into •Remove preservative solution by centrifugation at 3, •Close the bottom cap, add 7.5 mL of Binding Buffer, into the 50 mL conical tube and centrifuge at 3,000 x Repeat this step once more. If buffer remains in the column, centrifuge for longer times	a 50mL conical tube. 000 x rpm for 30 seconds. agitate the column, set it rpm for 30 seconds. ne.	
000 × g)	Step 2. Sample Apply •Close the bottom cap tightly and add prepared samp •Close the top cap and incubate for 1~2 hrs with mixi •Put off the bottom cap, set the column into a 50 mL centrifuge at 3,000 x rpm for 30 seconds.	ole to the column. ng for the gel not to sink . conical tube and Sample→ a	Flow Through
	 Step 3. Wash Put off the top cap, add 10 mL of Binding Buffer, closshake for 5 min. Put off the bottom cap and centrifuge at 3,000 x rpm Repeat this step more 2 times. (If non-specific proteins should be reduced, repeat to be addressed on the state of the sta	ose the bottom cap and n for 30 seconds. n for 5 times)	Wash
ing Buffer	r Step 4. IgG-Elution · Close the bottom cap tightly and add 5 mL of Elution · Close the top cap and shake for 5 min. · Put off the bottom cap, set the column into a 50 mL including Neutralization Buffer and collect eluate in the centrifugation at 3,000 x rpm for 30 seconds. · For 2nd and 3rd elution, add 2.5 mL of Elution Buffer in another 50 mL conical tube including Neutralization	n Buffer. conical tube Elut. buffer – tube by r to the column, shake 5 min, on Buffer, respectively.	IgG eluate Neutral. buffer
	Approx. 37% of purified IgG is collected in 1st eluate and 43 If higher concentration of IgG is needed, use mixture of 1st a If higher amount of IgG is needed, use mixture of 1st to 3rd	% of it is in 2nd eluate. and 2nd eluate (approx. 80% recovery) I eluate (approx. 97% recovery).).
	Note • Spin column (included) is available up to 3,000 x g • Centrifuge conditions: Swing rotor (r = 10 cm), at 3 (standard protocol) Angle rotor (r = 8.1 cm), at 3 • Some of angle rotors, if top of the column touch the • If sample volume is over the tube capacity, dischar	3,000 rpm (1,000 x g), for 30 sec. 3,000 rpm (815 x g), for 30 sec. e lid of centrifuge, cut short the co rge of the sample cannot complete	onical tube. e at 1st

centrifugation. After 1st centrifugation, discard the fluid in the tube and centrifuge again.