## **Technical Information**

<u>Protein A-R28</u> is an alkali-tolerant IgG-binding protein derived from protein A, which is developed with ProteNova's patent technology. Protein A-R28 strongly binds to various species and subclasses of IgG, compared with Protein A and G. The coupling to resin (Ab-Capcher) provides an alkali-washable unique affinity medium with high binding capacity for immunoglobulin, which is useful for purification of human, rabbit, mouse and rat IgGs. Ab-Capcher is also useful for immuno-precipitation experiments.

Table 1.	Binding	properties	of Ab-	-Capcher	(Ab-Rapid	PuRe)
----------	---------	------------	--------	----------	-----------	-------

Species	Sub class	Ab-Cap	Protein A	Protein G
Mouse	IgG1 IgG2a	++++	+	++
Rat	IgG1 IgG2a	++++ +++	- -	+ +
Goat	lgGs	++++	-	+
Chicken	lgY	-	-	-
Human	lgG	+++++	++++	++
Rabbit	lgG	+++++	++++	++

Note: In some species of antibody, binding to the gel may be weak.

#### **Order Information**

Product Name	Contents	Code No.
Ab-Capcher ExTra	2 mL	P-003-2
	10 mL	P-003-10
	100 mL	P-003-100
Ab-Rapid PuRe EX 2	Column x 2, 2.5 mL syringe x 1	P-015-2
Ab-Rapid PuRe EX 10	Column x 10	P-015-10
Ab-Rapid SPiN EX (Spin column)	0.1 mL gel/column x 10 (20 empty 2 mL-tubes included)	P-014-10

There are cases that prices will be changed without notice. For research use only.



# COSMO BIO CO., LTD.

Inspiration for Life Science

TOYO EKIMAE BLDG. 2-20, TOYO 2CHOME KOTO-KU, TOKYO 135-0016, JAPAN

TEL: +81-3-5632-9617 FAX: +81-3-5632-9618 ProteNova Co.Ltd. Takamatsu Lab. 2217-44 Hayashi-cho, Takamatsu Kagawa 761-0301 Japan TEL +81-87-897-2073 FAX +81-87-816-2073 URL http://protenova.com



# **Ab-Rapid PuRe EX**

Users Manual P-015-2 Ver.1.4

In some molecular species of Rat IgG2a, binding to the gel may be weak (EX : about 1mg/mL gel)

<sup>•</sup>In mouse IgM, there are 2 type of molecular species. "High-binding" type can be purified with this protocol, but "low-binding" type is difficult to be purified.

## Ab-Rapid PuRe Specifications

"Endotoxin-tested " Ab-Capcher ExTra is prepacked.

•Gel volume: 0.5 mL

•Gel matrix: Highly-crosslinked-agarose

(Rapid Run)

•Column volume: 0.8 mL •Particle size: 35 µm

Ligand: Alkali-resistant Protein A derivative

(Protein A-R28)

Binding Capacity: >45 mg human lgG /column

•Form: 20% ethanol •Storage: 4-8°C

## Materials

2.5 mL syringe

Microcentrifuge tube

D. # - "

Buffers

Binding Buffer: PBS

Elution Buffer: 0.1 M Glycine-HCl, pH 2.8

Nuetralization Buffer: 1 M Tris

\* Buffer Kit (PN-011) is also available from ProteNova.

Buffer kit contains Binding Buffer, Elution Buffer and Nuetralization Buffer.

## Sample preparation (example)

- ◆ Ascites : 3 x dilution with Binding Buffer.
- ◆Serum: precipitation with 50%-saturated (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> or 5x dil. with Binding Buffer
- ◆Cultured medium: Adjust pH to neutral.

Recommended pre-treatments of sample before applying to the column.

- Centrifugation; 10,000 × g, 10 min
- Filtration; 0.45µm filter

(Please use low-protein-adsorption types)

\* If there are insolubles in the sample, make sure to do pre-treatments.

#### Preparation for 50% ammonium sulfate precipitation

- 1. Prepare saturated ammonium sulfate.
- Add equal volume of saturated ammonium sulfate gradually to serum and mix.
- 2. Stand on ice for more than 1hr.
- 3. After centrifugation at 4°C, remove the supernatant.
- Wash precipitate with 50%-saturated ammonium sulfate.
- 4. Resolve the precipitate with small volume of Binding Buffer. The precipitate contains antibody.
- 5. Exchange to Binding buffer with dialysis or desalting column.

## Protocol for IgG Purification

Column pressure is higher and flow rate is slower than those of Ab-Rapid PuRe.

#### Preparation

Luar lock adaptor

Ab-Capcher ExTra gel

Upper filter

Lower filter

Plua

\* If air bubble is present in the space between gel and column


using 10 mL syringe (flow-rate approx. 5mL/min). It is important to

Before use, pass through 10 mL of Binding Buffer into the column by

press the gel bed. Repeat this procedure until the bubble disappears.

Add Neutralization Buffer to a microcentrifuge tube. (1/30 volume of eluate; 30-35µL to 1 mL of eluate)

Remove a cap on the top of column and fit a luar lock adaptor.




#### Step 1. Equilibration of Column

Fill a syringe with 2.5 mL of Binding Buffer and connect to the top of column.

After removal of a seal and a plug from the bottom of column, pass through 2.5 mL of Binding Buffer at a flow rate of 2.5 mL/min.


Then, remove the syringe.



#### Step 2. Sample Apply to Column

Apply the sample solution to column using the syringe at a flow rate of 2 mL/min.

Then, remove the syringe.




### Step 3. Washing of Column

Wash the column with 2.5 mL of Binding Buffer with the syringe at a flow rate of 2.5 ml/min.

\* If serum was directly applied to the column, washing with more than 5 mL is recommended.


Then, remove the syringe.



### Step 4. Elution of IgG

Elute with 2 mL of Elution Buffer using the syringe at a flow rate of 2 mL/min.

For this, discard first 0.3-0.5 mL of eluate and collect the following 1.0 - 1.5 mL of eluate as IgG fraction into a microcentrifuge tube, in which Neutralization Buffer (1/30 volume of eluate) is pre-added, and mix.



Time: within 10 min

#### \* Storage and reuse of column

- •Ab-Rapid PuRe EX column is alkali-washable. After 10 times washing ({5 min x 3} x 10 times), it is confirmed that remained binding capacity is 95% with 0.5N NaOH and 88% with 1.0N NaOH.
- •When the column is reused, wash the column with 2.5 mL of 0.1-0.5 N NaOH by syringe after elution of IgG. Using immediately after washing, equilibrate with 2.5 mL of Binding Buffer twice. Then, apply the sample.
- •For storage of column, add 2.5 mL of 20% EtOH, tightly close a cap and a plug, and store at 4-8 °C.