


# Ab-Rapid SPiN™ IgG purification from various animal sera

## Operation flow chart

Various serums 0.3mL  
↓  
2-fold dilution with PBS  
↓  
Add to Ab-Rapid SPiN™  
↓  
Reaction (1 hour)  
↓  
Wash (PBS), 5 minutes x 3 times  
↓  
0.1 M Glycine-HCl (pH 2.5), 5 minutes  
mixing, elution  
(Add 1 M Tris to neutralize)

## Comparison of amount of IgG in serum



## Electrophoresis result

### SDS-PAGE


Gel: 12.5% (Tris-Glycine)

Stain: CBB

Lane 1: Serum

Lane 2: Flow through

Lane 3: Eluate



IgG purified from the serum of 5 different animal species using Ab-Rapid SPiN™ was subjected to yield and purity analysis. The amount of purified IgG was 20 mg/mL gel or more from the sera of all animals but not from FBS. Given the IgG content of serum is typically approximately 10 mg/mL, each purification was performed on 30 mg/mL gel amount of serum. Thus, IgG could be purified from various animal sera with sufficient recovery rate. Since FBS originally has a low IgG content, it is speculated that the yield was low in this purification.

Protenova Co., Ltd.

〒769-2604

1488 Nishimura, Higashikagawa City, Kagawa Prefecture

TEL 0879-49-0702 / FAX 0879-49-0703

Home page <http://protenova.com>