

Ab-Carrier™

Species-specific IgG transfection efficiency into HeLa cells

Experimental protocol

Seeding HeLa cells in a 12-well plate

 $(1.0 \times 10^5 \text{ cell/well})$; medium volume MEM (+ 10% FBS) 1 mL/well)

↓ 37° C in the presence of 5% CO2 for 24 hours

Ab-Carrier TM 1 μ L was added to 20 μ L of FITC-labeled IgG (0.1 mg/mL) and mixed well.

↓ room temperature for 20 minutes

Add 21 $\mu\text{L/well}$ of reaction solution to HeLa cells after 24-hour culture

 \downarrow 37° C in the presence of 5% CO2 for 4 hours

Remove the medium and wash with PBS 1 mL/well $\, imes\,$ 2 times

 \downarrow

Fluorescent microscope observation (EVOS digital inverted microscope; LMS)

1

Add 0.25% Trypsin-EDTA 100 µL/well

 \downarrow Incubation at 37° C in the presence of 5% CO2 for 2 minutes

Add 1 mL of MEM (+ 10% FBS)

1

FACS (CYTOMICS FC 500; Beckman)

Table 1. Antibody transfer activity (evaluated by FACS)

Upper row: Ratio of cells into which antibody was introduced (%)

Lower row: Amount of antibody introduced per cell (Relative fluorescence intensity)

Goat polyclonal IgG	93.2%
	14.0
Rabbit polyclonal IgG	92.6%
	17.9
Mouse monoclonal IgG1	81.1%
	3.6
Mouse monoclonal IgG2a	78.1%
	4.8




Figure 1. Fluorescence microscope image of HeLa cells 4 hours after introduction of FITC-labeled antibody

- (A) Goat polyclonal IgG
- (B) Rabbit polyclonal IgG
- (C) Mouse monoclonal IgG1
- (D) Mouse monoclonal IgG2a

Green: FITC labeled antibody

Scale bar: 200 µm

HeLa cells (1.0 x 10⁵ cells/mL) 24 hours after seeding in a 12-well plate were treated with a mixture of antibody transfection reagent Ab-CarrierTM and FITC-labeled IgG antibodies (goat, rabbit, mouse (IgG1, IgG2a). The cells were cultured for 4 hours and then observed with a fluorescence microscope (phase difference). Fluorescence was observed in most cells (Fig. 1). The proportion of cells into which the antibody was introduced (antibody introduction rate) and the amount of antibody introduced per cell (relative fluorescence intensity) were evaluated by FACS. Goat/rabbit IgG showed a high antibody transfer rate of 90% or more, mouse IgG1 of 80% or more, and mouse IgG2a of 75% or more.

Protenova Co., Ltd. 7769-2604

1488 Nishimura, Higashikagawa City, Kagawa Prefecture TEL 0879-49-0702 / FAX 0879-49-0703

Home page http://protenova.com