

Ab-Capcher ExTra™ and Protein L Purification of human IgG Fab fragment

Reaction protocol

Human IgG (10 mg/mL)

IgG papain digestion

- ↓ Add 1/10 volume of papain (1 mg/mL) and react at 37° C for 1 hour
- ↓ Add 1M iodoacetamide (final 15mM) to stop the reaction

Ab-Capcher ExTra™ packed column (equilibrated with PBS)

Undigested IgG and Fc fragment removal


- ↓ Add papain digest to the column
- ↓ Collect the raw material (FT) fraction

Protenova Protein L-immobilized gel packing column (equilibrated with PBS)

Recovery of Fab fragment

- ↓ Ab-Capcher ExTra™ FT fraction added to the column
- ↓ FT fraction collection
- ↓ Wash the column with PBS
- ↓ Elute with 0.1M glycine-HCl pH 2.5
- ↓ Neutralize the elution fraction with 1M Tris

Purification completed



SDS-PAGE

Gel: 15% a.a. (Tris-Glycine)

Stain: CBB

Lane 1: MW marker

Lane 2: human IgG

Lane 3: human IgG after papain digestion

Lane 4: Ab-Capcher ExTra™ FT fraction

Lane 5: Protein L elution fraction

Fab fragments were purified from human IgG using Ab-Capcher ExTra™ and Protenova Protein L-immobilized gel. Papain digestion reduced human IgG into Fab and Fc fragments (lane 3). By adding this papain digest to an Ab-Capcher ExTra™ packed column, undigested IgG and Fc fragment were adsorbed to the column, and the FT fraction containing Fab was collected (lane 4). This FT fraction was added to a Protenova Protein L packed column, and the eluted fraction was collected. SDS-PAGE confirmed that high purity Fab fragments were obtained in the protein L column elution fraction (lane 5). Thus, it was found that the combination of Ab-Capcher ExTra™ and Protenova Protein L-immobilized gel was effective for purification of Fab fragments.

Protenova Co., Ltd.

〒769-2604

1488 Nishimura, Higashikagawa City, Kagawa Prefecture

TEL 0879-49-0702 / FAX 0879-49-0703

Home page <http://protenova.com>