

Blend Taq-Plus-

Code: TYB-BTQ-201 (250 units)

Storage: Store at -20°C

Concentration: Blend Taq-Plus- 2.5 units/ul
Anti-Taq DNA polymerase antibody 0.5 ug/ul

Unit Definition: One unit of enzyme is defined as the amount of enzyme that will incorporate 10 nmoles of dNTPs into acid insoluble materials in 30 minutes at 75°C.

Assay Condition: 25 mM TAPS (pH 9.3)
50 mM KCl
10 mM MgCl₂
200 uM each of dATP, dGTP, dTTP
100 uM [α -³²P]-dCTP
20 ug activated salmon sperm DNA per 50 ul reaction

Storage Buffer: 20 mM Tris-HCl (pH 8.0)
100 mM KCl
0.1 mM EDTA
1 mM DTT
0.5 % Tween-20
0.5 % Nonidet P-40
50 % Glycerol

Materials Provided: 10× Buffer for Blend Taq (Mg²⁺ concentration: 20mM)
dNTPs : 2mM dATP, dGTP, dCTP, dTTP each

Quality Control

1. PCR Assay: The 17.5kbp fragment of human β globin gene could be amplified using human genomic DNA.

Purchase of this product is accompanied by a limited license to use it in the Polymerase Chain Reaction(PCR) process for Research Field in conjunction with a thermal cycler whose use in the automated performance of the PCR process is covered by the up-front license fee, either by payment to Perkin-Elmer or as purchased, i. e., an authorized thermal cycler.

PCR conditions

Reaction mixture for General PCR Amplification

10× Buffer for Blend Taq	5ul
Blend Taq-Plus- (2.5 U/ul)	0.5ul (1.25U)
Template	1-50ng (plasmid) 10-1000ng (Genome)
Primer	10-50pmoles (final conc.: 0.2-1.0uM)
2mM dNTPs	5ul (final conc.: 0.2mM)
Total volume	50ul

Suggested Cycling parameters

Segment	Target			Number of Cycles
	<1kb	1kb – 6kb	>6kb	
1	94°C 2min			1
2 Denaturation	94°C 30sec (see notes #1)			30
3 Annealing	Primer Tm-5°C (55-60°C) 30sec		68°C	
4 Extension	72°C 1min	72°C 1min/kb PCR target	1min/kb PCR target (see notes #2, #3)	

Notes

1. The denaturation time and temperature may require optimization. Typically denaturation time will range 25 - 30sec at 94°C or 5 - 10sec at 98°C.
2. Toyobo suggests >70°C melting temperature of primer (Tm) for long target (>6kb).
3. A dNTP concentration range of 0.3 - 0.4mM total is recommended for longer target (>10kb).
4. Blend Taq produces amplification products that are ready to clone directly into TA cloning vectors without any special optimization. To get sufficient numbers of transformants and white colony, longer ligation time will be efficient, recommended 2 hours – overnight for the amplification products of >1kb.
5. Blend Taq-Plus- is designed to be suitable for Hot Start PCR. Non-specific amplification due to mispriming and/or formation of primer dimer before thermal cycling can be prevented.

For research use only; not for use as a diagnostic.

Distributor



COSMO BIO CO., LTD.
Inspiration for Life Science

TOYO 2CHOME, KOTO-KU, TOKYO, 135-0016, JAPAN

http://www.cosmobio.co.jp

e-mail : export@cosmobio.co.jp

Phone : +81-3-5632-9617

FAX : +81-3-5632-9618