

www.toyobo.co.jp/e/bio

Instruction manual THUNDERBIRD Probe qPCR Mix 0910

A4250K

THUNDERBIRD™ Probe qPCR Mix

QPS-101T 1 mL x 1 QPS-101 1.67 mL x 3

Store at -20°C, protected from light

Contents

- [1] Introduction
- [2] Components
- [3] **Primer/Probe design**
- [4] Template DNA
- [5] Protocol
 - 1. Standard reaction set up
 - 2. Cycling conditions
 - 2-1. Real-time PCR conditions using Applied Biosystems 7900HT
 - 2-2. Real-time PCR conditions using Roche LightCycler 1.1
- [6] Related Protocol

cDNA synthesis

- [7] **Troubleshooting**
- [8] Related products

CAUTION

All reagents in this kit are intended for research purposes only. Do not use for diagnosis or clinical purposes. Please observe general laboratory precautions and observe safety procedures while using this kit.

-LightCycler™ is a trademark of Idaho Technology, Inc. and Roche Molecular Systems, Inc.

-TaqMan[®] is a registered trademark of Roche Molecular Systems, Inc. -SYBR[®] is a registered trademark of Roche Molecular Systems, Inc.

-SYBR is a registered trademark of Roche Molecular Systems, Inc.

JAPAN TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



[1] Introduction

Description

THUNDERBIRDTM Probe qPCR Mix is a highly efficient 2x Master Mix for real-time PCR using TaqMan[®] probes. The master mix contains all required components, except for ROX reference dye, probe and primers (50x ROX reference dye is individually supplied with this kit). The master mix facilitates reaction setup, and improves the reproducibility of experiments.

This product is an improved version of Realtime PCR Master Mix (Code No. QPK-101). In particular, reaction specificity and PCR efficiency is enhanced.

Features

-High specificity

The specificity for the detection of low-copy targets is improved.

-Homogeneous amplification

The dispersion of PCR efficiency between targets is reduced by a new PCR enhancer*. (*Patent pending)

-Broad dynamic range

High specificity and effective amplification enable the detection of a broad dynamic range.

-Compatibility for various real-time cyclers.

The reagent is applicable to most real-time cyclers (i.e. Block type and glass capillary type). Because the 50x ROX reference dye is individually supplied with this kit, the kit can be applied to real-time cyclers that require a passive reference dye.

-Hot start PCR

The master mix contains anti-Taq DNA polymerase antibodies for hot start technology. The antibodies are easily inactivated in the first denaturation step, thereby activating the DNA polymerase.

About the fluorescent probe detection system

The TaqMan[®] probe system utilizes fluorescence emission from the probes. The probes hybridize to the target amplicons and then emit fluorescence upon degradation by the 5'-3' exonuclease activity of Taq DNA polymerase. This type of detection system can achieve higher specificity in real-time PCR assays than the SYBR[®] Green I detection system.

[2] **Components** This kit includes the following components for 40 reactions (QPS-101T) and 200 reactions (QPS-101), with 50 µl per reaction. All reagents should be stored at -20°C.

<QPS-101T> THUNDERBIRD™ Probe qPCR Mix 1 ml x 1 50x ROX reference dye 50 µl x 1

JAPAN TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



Notes:

-THUNDERBIRD[™] Probe qPCR Mix can be stored, protected from light, at 2-8°C for up to 3 months. For longer storage, this reagent should be kept at -20°C and protected from light. No negative effect was detected by 10 freeze-thaw cycles of THUNDERBRID[™] Probe qPCR Mix. This reagent does not contain the ROX reference dye.

-50x ROX reference dye can be stored, protected from light, at 2-8°C or -20°C. For real-time cyclers that require a passive reference dye, this reagent must be added to the reaction mixture at a concentration of 1x or 0.1x. The master mix solution with the ROX reference dye can be stored, protected from light, at 2-8°C for up to 3 months. For longer storage, this reagent should be kept at -20°C and protected from light. The pre-mixed reagents can be prepared according to the following ratios. [5] Table 1 shows the optimal concentration of the ROX dye.

1x solution

THUNDERBIRDTM Probe qPCR Mix : 50x ROX reference dye = 1.67 ml : 66.8 µl THUNDERBIRDTM Probe qPCR Mix : 50x ROX reference dye = 1 ml : 40 µl

0.1x solution

THUNDERBIRDTM Probe qPCR Mix : 50x ROX reference dye = $1.67 \text{ ml} : 6.7 \mu \text{l}$ THUNDERBIRDTM Probe qPCR Mix : 50x ROX reference dye = $1 \text{ ml} : 4 \mu \text{l}$

For real-time cyclers that do not require a passive reference dye, THUNDERBIRDTM Probe qPCR Mix without the ROX reference dye can be used.

[3] Primer/Probe design

1. Primer conditions

Highly sensitive and quantitative data depend on primer design. The primer should be designed according to the following suggestions;

-Primer length: 20-30 mer

-GC content of primer: 40-60%

- -Target length: \leq 200 bp (optimally, 80-150 bp)
- -Melting temperature (Tm) of primers: 60-65°C
- -Purification grade of primers: Cartridge (OPC) grade or HPLC grade

Notes:

-Longer targets (>200 bp) reduce efficiency and specificity of amplification.

-Tm of the primers can be flexible, because the Tm value depends on the calculation formula.

JAPAN

TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



2. Fluorescent probe

The probes should be designed according to the guidelines of each probe system. Because insufficiently purified probes may inhibit the reaction, HPLC-grade probes should be used.

[4] **Template DNA** The following DNA samples can be used as templates.

1. cDNA

Non-purified cDNA, generated by reverse transcription reactions, can be used directly for real-time PCR using THUNDERBIRDTM Probe qPCR Mix. Up to 10% of the volume of a cDNA solution can be used for a real-time PCR reaction. However, excess volume of the cDNA may inhibit the PCR. Up to 20% (v/v) of the cDNA solution from ReverTra Ace[®] qPCR RT Kit (Code No. FSQ-101) can be used for real-time PCR (see [6]).

2. Genomic DNA, Viral RNA

Genomic DNA and viral RNA can be used at up to 200 ng in 50 μl reactions.

3. Plasmid DNA

Although super-coiled plasmids can be used, linearized plasmid DNA produces more accurate assays. The copy number of the plasmid DNA can be calculated by the following formula.

Copy number of 1µg of plasmid DNA = 9.1×10^{11} / Size of plasmid DNA (kb)

JAPAN

TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



[5] Protocol

1. Reaction mixture setup

	Reaction vo	lume	Final	
Reagent	50 µl	20 µl	Concentration	
DW	X μl	Xμl		
THUNDERBIRD™ Probe qPCR Mix	25 µl	10 µl	1x	
Forward Primer	15 pmol	6 pmol	$0.3 \ \mu M^{*1}$	
Reverse Primer	15 pmol	6 pmol	$0.3 \ \mu M^{*1}$	
TaqMan [®] Probe	10 pmol	4 pmol	$0.2 \; \mu M^{*_1}$	
50x ROX reference dye	1µl / 0.1 µl	0.4µl / 0.04	$\mu l = 1x / 0.1x^{*}$	
DNA solution	Υµl	Yμl		
Total	50 µl	20 µl		

Notes:

^{*1} Primer / probe concentration should be determined according to the manufacturer's instructions.

Higher primer concentration tends to improve the amplification efficiency, and lower primer concentration tends to reduce the non-specific amplification. The primer concentration should be set between $0.2-0.6 \ \mu M$.

*2 50x ROX reference dye must be added when using real-time cyclers that require a passive reference dye, according to Table 1. Table 1 shows the optimum concentration of the ROX reference dye. This dye is not necessary for real-time cyclers that do not require a passive reference dye.

Table 1 Recommended ROX dye concentration

Real-time cycler	Optimal dye concentration
	(dilution ratio)
Applied Biosystems 7000, 7300, 7700, 7900HT etc.	1x (50:1)
Applied Biosystems 7500, 7500Fast,	0.1x (500:1)
Stratagene cyclers (Optional) etc.	
Roche' cyclers, Bio-Rad cyclers, BioFlux cyclers etc.	Not required

Notes:

The ROX dye in Realtime PCR Master Mix (Code No. QPK-101) corresponds to 1x concentration.

JAPAN

TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp **CHINA** TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



2. PCR cycling conditions

The following table shows the recommended thermal conditions using primers designed according to the recommended primer and probe conditions described in [3]. Almost all targets can also be amplified using the ongoing conditions with other real-time PCR reagents.

<3-step cycle>	Temperature	Time	Ramp	
Pre-denaturation:	95°C	$20-60 \text{ sec}^{*1}$	Maximum	
Denaturation:	95°C	$1-15 \text{ sec}^{*2}$	Maximum	4 0 avalas
Extension:	60°C*3	30-60 sec ^{*4}	Maximum	40 cycles
	(data collection	should be set at t	he extension step)	

^{*1} Due to the anti-Taq antibody hot start PCR system, the pre-denaturation can be completed within 60 sec. The pre-denaturation time should be determined according to the recommendations of each real-time cycler. If the optimal pre-denaturation time cannot be determined, the time should be set at 60 sec.

Table 2	The recommended	pre-denaturation time for each real-time cycler
---------	-----------------	---

Real-time cycler	Pre-denaturation time
High speed cycler (e.g. Applied Biosystems 7500Fast)	20 sec
Capillary cycler (e.g. Roche LightCycler [™] 1.x, 2.0)	30 sec
General real-time cyclers (Applied Biosystems 7700, 7500,	60 sec
7900HT (normal block), Stratagene cyclers, BioFlux	
cyclers	

^{*2} The following table shows the optimal denaturation times for each real-time cycler. If the optimal denaturation time cannot be determined, the time should be set at 15 sec.

Real-time cycler	denaturation time
High speed cycler (e.g. Applied Biosystems 7500Fast)	3 sec
Capillary cycler (e.g. Roche LightCycler [™] 1.x, 2.0)	5 sec
General real-time cyclers (Applied Biosystems 7700, 7500,	15 sec
7900HT (normal block), Stratagene cyclers, BioFlux	
cyclers	

 Table 3
 The recommended denaturation time for each real-time cycler

*3 Insufficient amplification may be improved by decreasing the extension temperature, and non-specific amplification (e.g. abnormal shapes of the amplification curve at low template concentrations) may be reduced by increasing the extension temperature. The extension temperature should be set at 56-64°C.

^{*4} If the target size is smaller than 300 bp, the extension time can be set at 30 sec on almost all real-time cyclers. Instability of the amplification curve or variation of data from each well may be improved by setting the extension time at 45-60 sec. Some real-time cyclers or software need over 30 sec for the extension step. In these cases, the time should be set according to each instruction manual (e.g. Applied Biosystems 7000/73000: \geq 31 sec; Applied Biosystems 7500: \geq 35 sec.).

JAPAN	CHINA
TOYOBO CO., LTD.	TOYOBO Bio
Tel(81)-6-6348-3888	Tel(86)-21-58
www.toyobo.co.jp/e/bio	
tech_osaka@toyobo.jp	

TOYOBO Bio-Technology, CO., LTD. [el(86)-21-58794900.4140



2-1. Real-time PCR conditions using Applied Biosystems 7900HT

(Normal block type, software version 2.2.2)


The following is an example of a TaqMan[®] assay using Applied Biosystems 7900HT.

(1) The cycling parameters should be set according to the following "Thermal Cycler Protocol" window under the "Instrument" tab.

Setup Instrument		
Thermal Cycler Protocol		?
Mode: C Fast @ Standa	rd 🔘 9600 Emulation	Sample Volume (µL): 25
	nent Ramp Rate Data Collection	
Stage 1	Stage 2 Repeats 40	
95.0	95.0	
1:00	0.15	
	60.0	
	0:45	
1		
Add Cycle Add H	Add Step Delete Step	Add Dissociation Stage

Notes:

- Appropriate sample volumes should be set.
- ≥ 45 sec is necessary for the extension step.
- (2) Click the "Data collection" tab.



- (3) Insert the PCR plate
- (4) Start the program

JAPAN
TOYOBO CO., LTD.
Tel(81)-6-6348-3888
www.toyobo.co.jp/e/bio
tech_osaka@toyobo.jp

CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



2-2. Real-time PCR conditions using Roche LightCycler 1.1

(Software version 3.5)

The following is an example of a TaqMan[®] probe assay using Roche LightCycler 1.1.

(1) The cycling parameters should be set according to the following window. Analysis and Acquisition mode of the initial denaturation step must be set at "None".

Cycle Program Data Cycles 1	Analysis Mode None Quantification Melting Curves
Temperature Targets	
	e (°C) me (hrs:min:sec) ─Secondary Target Temperature (°C) ─Step Size (°C) ─Step Delay (cycles) ─Acquisition Mode
(ns) 95 ‡ 30 ‡ 20.00 ‡ 0	₿0.0 ₿0 ₿ NONE ₿ (De)

(2) The cycling parameters should be set according to the following window. Analysis mode of the PCR step must be set at "Quantification". Acquisition modes of 95°C and 60°C must be set at "None" and "Single", respectively.

Cycle Program	ı Data	Analysis	Mode None
Cycles 40 🖡			Quantification Melting Curves
Temperature ⁻	Targets		
	get Temperatura Incubation Tir Tem	mè (hrs:min:s perature Tra —Secondar	sec) Instition Rate (°C / s) Y Target Temperature (°C) p Size (°C) C Step Delay (cycles) Acquisition Mode
(Ins) 95 算 5	20.00	₿0.0	
(hs) 60 ₿ 30	20.00	0.0	SINGLE

JAPAN

TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140

7



(3) The cycling parameters should be set according to the following window. Analysis and Acquisition mode of the cooling step must be set at "Non".

Cycle Program Data	Analysis Mode <mark> None</mark>
Cycles 1	Quantification Melting Curves
Temperature Targets	
Target Tempera	
	Time (hrs:min:sec)
	emperature Transition Rate (*C / s)
	Secondary Target Temperature (°C)
	Step Delay (cycles)
	Acquisition Mode
(ins) 40 \$ 30 \$ 20.00 \$	
Ins	

(4) Insert the capillaries in the carousel, and start the cycling program.

JAPAN

TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



[6] Related Protocol 1. cDNA synthesis

cDNA synthesized by various cDNA synthesis reagents can be used with THUNDERBIRD[™] Probe qPCR Mix. However, cDNA synthesized by a reagent specialized for real-time PCR can increase sensitivity.

ReverTra Ace[®] qPCR RT Kit (Code No. FSQ-101) is a cDNA synthesis kit suitable for real-time PCR. Here, the protocol with ReverTra Ace[®] qPCR RT Kit is described. However, for the detailed protocol, please refer to the instruction manual of the kit.

(1) Denaturation of RNA

Incubate the RNA solution at 65°C for 5 min and then chill on ice.

Notes:

-This step can be omitted. But this step may increase the efficiency of the reverse transcription of RNA, which forms secondary structures.

-Do not add 5x RT Buffer and/or enzyme solution at this step.

(2) Preparation of the reaction solution

Reagent	Volume (amount)
Nuclease-free Water	X μl
5x RT Buffer	2 µl
Primer Mix	0.5 µl
Enzyme Mix	0.5 µl
RNA solution	0.5 pg-1 μg
Total	10 µl

(3) Reverse transcription reaction

-Incubate at 37°C for 15 min. <Reverse transcription>

-Heat at 98°C for 2 min. <Inactivation of the reverse transcriptase>

-Store at 4°C or -20°C.*

*This solution can be used in the real-time PCR reaction directly or after dilution.

Notes:

The above temperature conditions are optimized for ReverTra Ace® qPCR RT Kit.

JAPAN

TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



[7] Troubleshooting

Symptom	Cause	Solution
Loss of linearity in the high cDNA/DNA concentration region.	Inhibition by the components in the cDNA/DNA solution.	 -DNA: The DNA sample may contain PCR inhibitors. The DNA samples should be repurified. -cDNA: The components in the cDNA synthesis reagent may inhibit the PCR reaction. The cDNA sample should be used after dilution.
Lost of linearity or	The template DNA is insufficient. Adsorption of the DNA to the tube wall.	When the DNA/cDNA copy number is lower than 10 copies per reaction, the linearity of the reaction tends to be lost. The template concentration should be increased. The diluted DNA templates tend to be absorbed onto the tube wall. Dilution should be performed
lower signal in the low DNA/cDNA concentration region.	Competition with primer dimer formation.	just prior to experiments. In the probe assay, primer dimers are not detected. However, dimer formation may reduce the amplification efficiency of the target, especially for reactions at low template concentration. The reaction conditions should be optimized or the primer sequences should be changed.
Loss of linearity of the amplification carves.	Competition with non-specific amplification.	In the probe assay, non-specific amplification is not detected. However, non-specific amplification may reduce the amplification efficiency of the target. The reaction conditions should be optimized or the primer sequences should be changed.
The PCR efficiency is	Inappropriate cycling conditions.	Optimize the cycling conditions according to [5].
lower than 90% (slope:	Degradation of the primers.	Fresh primer solution should be prepared.
<-3.6)	The calculation of the PCR	The Ct value on the linear region should be used
The DOD officiant in	efficiency is inappropriate. The calculation of the PCR	to calculate PCR efficiency.
The PCR efficiency is higher than 110%.	efficiency is inappropriate.	The Ct value on the linear region should be used to calculate PCR efficiency.
Reproducibility is not good.	Poor purification of the template DNA.	Low-purity DNA may contain PCR inhibitors. Re-purify the DNA samples.
	Absorption of the template DNA to the tube wall.	Diluted DNA templates tend to be absorbed onto the tube wall. Dilution of the template DNA/cDNA should be performed just prior to experiments.
	Plasmid DNA or PCR product is used as a template.	In general, plasmid DNA or PCR product is used at low concentration. Diluted DNA templates tend to be absorbed onto the tube wall. Dilution of the template DNA/cDNA should be performed just prior to experiments. Dilution with a carrier nucleic acid solution (Yeast RNA) is also effective in improving linearity.
	Inappropriate thermal conditions.	Optimize the thermal conditions according to [5].
	Low purity of the primers or probes.	Different lots of primers or probes may show different results. When the lot is changed, prior testing of the primer or probe should be performed.

JAPAN

TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp CHINA

TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140

10



Symptom	Cause	Solution	
Amplification from the	Contamination or carry over of the	Change the contaminated reagent.	
non-template control	PCR products.		
(NTC).	Inappropriate settings of	In multiplex experiments, inappropriate setting of	
	fluorescence measurement, such as	fluorescence measurement may cause the	
	in the case of multiplex PCR.	detection of noise by the cross talk of fluorescent dyes. Settings should be reconfirmed.	
Low amplification	Excessive amount of ROX reference	Excessive amount of ROX reference dye may	
curve signal /	dye.	cause low signal. 50x ROX reference dye should	
Unstable amplification		be used according to [5] Table 1.	
curve signal.	Inappropriate settings of	Settings should be confirmed according to the	
	fluorescence measurement.	instruction manual of each detector.	
	Low purity of fluorescent probes.	Low purity of the probe may increase the base	
		line. HPLC grade probes should be used.	
	Excessive intensity of the quencher	Certain quenchers (e.g. TAMRA) may cause a	
	Dye.	higher baseline because of its fluorescence. Use	
		of a non-fluorescent quencher may improve the	
		high baseline.	
	Degradation of the probe.	Store the probes according to the manufacture's	
		recommendations.	
	Insufficient fluorescence	Certain detection systems require a longer time to	
	measurement time.	detect the fluorescent signal. Longer extension	
		(measurement) time (45-60 sec) may improve the	
		unstable signal.	
	Insufficient reaction volume.	Low reaction volume may cause an unstable	
		signal. Increase the reaction volume.	

JAPAN

TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140

11



[8] Related products

Product name	Package	Code No.
High efficiency real-time PCR master mix	200	QPS-201
THUNDERBIRD [™] SYBR [®] qPCR Mix	reactions	
High efficiency cDNA synthesis kit for real-time PCR	200	FSQ-101
ReverTra Ace [®] qPCR RT Kit	reactions	
One-step Real-time PCR master mix for probe assay	0.5 mL x 2	QRT-101T
RNA-directTM Realtime PCR Master Mix	0.5 mL x 5	QRT-101
One-step Real-time PCR master mix for SYBR® Green assay	0.5 mL x 2	QRT-201T
RNA-direct TM SYBR [®] Realtime PCR Master Mix	0.5 mL x 5	QRT-201

JAPAN

TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140



NOTICE TO PURCHASER: LIMITED LICENSE

A license to perform the patented 5' Nuclease Process for research is obtained by the purchase of (i) both Authorized 5' Nuclease Core Kit and Licensed Probe, (ii) a Licensed 5' Nuclease Kit, or (iii) license rights from Applied Biosystems.

This product is an Authorized 5' Nuclease Core Kit. Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,079,352, 5,789,224, 5,618,711, 6,127,155, 5,677,152, 5,773,258, 5,407,800, 5,322,770, 5,310,652, 5,210,015, 5,487,972, and claims outside the US corresponding to US Patent No. 4,889,818. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. Separate purchase of a Licensed Probe would convey rights under the applicable claims of US Patents Nos. 5,538,848, 5,723,591, 5,876,930, 6,030,787, 6,258,569, 5,804,375 (claims 1-12 only), and 6,214,979, and corresponding claims outside the United States. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained from the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

JAPAN

TOYOBO CO., LTD. Tel(81)-6-6348-3888 www.toyobo.co.jp/e/bio tech_osaka@toyobo.jp CHINA TOYOBO Bio-Technology, CO., LTD. Tel(86)-21-58794900.4140