


### Example 5 Polarized cell culture using permeable collagen membranes

[Reference : Yoshizato K, et al. *J Cell Sci.* 91(Pt 4): 491-499]



When cultured using collagen membranes, it was observed that the glucose consumption and lactic acid excretion in the medium on the basal side were significantly higher than on the apical side. On the other hand, a decrease in pH on the basal side was observed, but when a membrane transporter inhibitor was added to the medium, the difference in pH decreased. These results indicate that culture on collagen membranes can maintain cell polarity and maintain a functional membrane transport mechanism. (Reference 5)

### Reference

- Murakami Y, et al. Reconstitution of the embryonic kidney identifies a donor cell contribution to the renal vasculature upon transplantation. (2019) *Sci Rep.* Feb 4;9(1):1172.
- Matsumoto A, et al. Induction of enamel matrix protein expression in an ameloblast cell line co-cultured with a mesenchymal cell line in vitro. (2011) *In Vitro Cell Dev Biol Anim.* 47(1):39-44.
- Ueno T, et al. Characteristic Gene Expression Profiles of Human Fibroblasts and Breast Cancer Cells in a Newly Developed Bilateral Coculture System. (2015) *Biomed Res Int.* 2015;2015:960840.
- Hatou S, et al. Functional corneal endothelium derived from corneal stroma stem cells of neural crest origin by retinoic acid and Wnt/ $\beta$ -Catenin signaling. (2013) *Stem Cells Dev.* 22(5):828-39.
- Yoshizato K, et al. Functionally polarized layers formed by epidermal cells on a permeable transparent collagen film. (1988) *J Cell Sci.* 91(Pt 4):491-499.

Cat. No.	Description	Quantity	Storage
KOU-CM-6	Atelocollagen Permeable Membrane (6 well)	12 pcs/box	4°C
KOU-CM-24	Atelocollagen Permeable Membrane (24 well)	24 pcs/box	4°C
KOU-CLF-01	Atelocollagen Permeable Membrane Sheet	1 sheet/box	4°C

World distributor




2792 Loker Avenue West, Suit 101, Carlsbad, CA 92010 USA  
Phone : +1 760-431-4600 / Fax : +1 760-431-4604

Please don't hesitate to contact us with any questions related to product selection or use.

<https://www.cosmobioussa.com/>

# Atelocollagen, Permeable membrane



**Culture on membrane**

**Air liquid Interface culture**

**Co-culture**

**Cell sheet transplantation**

# Atelocollagen, Permeable membrane

CM-6, CM-24, CLF-01


## Product summary

Atelocollagen, Permeable membranes are permeable collagen membranes that are made from atelocollagen and are permeable to substances. By seeding cells on both sides of the membrane, it is possible to observe cell-cell interactions without cells mixing. In addition, when epithelial cells are cultured on this membrane, it is possible to culture them while maintaining the polarity of the cells in the same way as in vivo, such as absorption of nutrients and excretion of metabolic wastes through the collagen membrane.

## Applications


Culture on membrane, Air liquid Interface culture, Co-culture, Cell sheet transplantation

## How to use



## Example 1 On-membrane culture of renal organoids using permeable collagen membranes


[Reference : *Sci Rep.* 2019 Feb 4;9(1):1172. Created by modifying figures 3A and 3C-G. ©Murakami Y., et al. 2019 (Licensed under CC BY 4.0) <https://creativecommons.org/licenses/by/4.0/>]



When kidney organoids reconstituted using a low-adhesion plate were cultured at the air-liquid interface on a collagen membrane, extensive ureteric bud molecules, surrounding vascular networks, and the formation of numerous glomeruli were observed. The above results show that, although there are certain limitations, it is possible to create renal organoids with vascular structures using this reconstruction method and air-liquid interface culture on collagen membranes. (Reference 1)

## Example 2 Induction of ameloblast differentiation using the permeable collagen membrane


(Taniguchi A, et al., Natl. Inst. for Material Science, Japan.)



For SW culture, rat odontogenic epithelium cell line (HAT-7) and bovine dental follicle cell line (BCPb8) were seeded on opposite sides of a collagen membrane. As a result, expression of Ameloblastin and Amelogenin (important proteins for enamel formation) were markedly increased compared to SC culture or monolayer culture. These experiments revealed that the permeability of collagen membrane is advantageous for remodeling cell-cell interactions in a manner similar to in vivo. (Reference 2)

## Example 3 Co-culture of Cancer Cells and Fibroblasts Using a Permeable Collagen Membrane


[Biomed Res Int. 2015;2015:960840. Created by modifying figure 1 and Tables 5-6. ©Ueno T., et al. 2015 (Licensed under CC BY 3.0)]



Using a collagen membrane, three types of breast cancer cell lines—HCC1937, MCF-7, and SK-BR-3—were co-cultured with fibroblasts, and gene expression changes were evaluated via microarray analysis. As a result, among the top 100 genes that were upregulated or downregulated in HCC1937 cells, 77% of the upregulated genes showed downregulation in fibroblasts, and 82% of the downregulated genes showed upregulation, indicating a reciprocal effect through co-culture. Furthermore, while all breast cancer cells exhibited reprogramming toward aerobic metabolism, the fibroblasts appeared to be reprogrammed toward a Warburg-like phenotype. These findings suggest that the collagen membrane is effective for evaluating intercellular interactions under non-contact co-culture conditions (Reference 3).

## Example 4 Induction of corneal endothelial differentiation and transplantation with the permeable collagen membrane

(Hato S, et al., Keio University, Japan.) (Licensed under CC BY 3.0)



Multipotent cornea-derived precursors isolated from mouse and human corneal stroma were cultured on collagen membrane and differentiated to functional corneal endothelium. Transplanting TECE attached collagen membrane to donor cornea where corneal endothelium was stripped suppressed edematous and notably increased the transparency of cornea compared to control groups. These results show that collagen membrane, with its high transparency and biocompatibility, is suitable for transplantation, such as corneal transplantation. (Reference 4)