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For serum-free cell culture

Kyokuto FM-1 Cryopreservation Medium

This cryopreservation medium for serum-free cell cultures was introduced by Ohno and colleagues in 1988. Traditionally, cell freezing has relied on media containing dimethyl sulfoxide, glycerin, and serum. However, as serum-free cell culture has become more common, the demand for serum-free freezing media has grown. Recently, it has been found that using methylcellulose as a serum substitute can improve cell recovery rates. This medium is formulated based on that approach.

This medium can also be used to freeze cells cultured with regular serum. The procedure is the same as for serum-free cultures.

Composition	Per 1 L
Eagle MEM Medium (for autoclave sterilization)	9.4g
Methylcellulose	2.0g
Dimethyl sulfoxide (DMSO)	200.0mL
L-Glutamine	0.3g
Sodium bicarbonate	1.0g
Purified water	800.0mL

Composition of Eagle MEM Medium		(mg/L)	(mg/L)
L-Arginine hydrochloride	126.0	Calcium Pantothenate	1.0
L-Cysteine hydrochloride (monohydrate)	31.4	Pyridoxal Hydrochloride	1.0
L-Tyrosine	36.0	Riboflavin	0.1
L-Histidine hydrochloride (monohydrate)	42.0	Thiamine Hydrochloride	1.0
L-Isoleucine	52.0	Biotin	0.02
L-Leucine	52.0	Sodium Chloride	6,800.0
L-Lysine Hydrochloride	73.0	Potassium Chloride	400.0
L-Methionine	15.0	Calcium Chloride (anhydrous)	200.0
L-Phenylalanine	32.0	Magnesium Sulfate (anhydrous)	93.5
L-Threonine	48.0	Sodium Dihydrogen Phosphate (anhydrous)	115.0
L-Tryptophan	10.0	Glucose (anhydrous)	1,000.0
L-Valine	46.0	Succinic Acid	75.0
Choline Bitartrate	1.8	Sodium Succinate (hexahydrate)	100.0
Folic Acid	1.0	Kanamycin	60.0
Inositol	2.0	Phenol Red	6.0
Nicotinamide	1.0		

Operation Method

FM-1-Based Cell Freezing Procedure

Please perform all procedures under sterile conditions.

This medium contains cryoprotectant at double the usual concentration.

1. Base Medium

For each tube of medium (5 mL), add 5 mL of cell suspension (prepared in your usual culture medium) for a 1:1 ratio. Alternatively, add an equal volume of this medium to the cell suspension.

2. Transfer the mixture from step 1 into a sterile serum tube or ampoule, seal tightly, and place it in a freezing container made of Styrofoam,

-80 °C or lower for storage. You may also freeze using a programmable freezer.

Prepare the cell density to be
10 cells/mL

Please prepare the cell suspension at a
concentration of at least 10^6 cells/mL.

Note: You may notice methylcellulose precipitates in this medium, but its performance remains unaffected.

Thawing and Re-culturing Frozen Cells

Quickly immerse the ampule stored at -80°C or below into a 37°C water bath and agitate gently by hand to rapidly thaw.

After thawing, open the ampule and dilute the contents 5- to 10-fold with standard culture medium, then centrifuge at about 1000 r.p.m. for 3–5 minutes.

Discard the supernatant, add the required amount of standard culture medium, transfer to a culture vessel, and start incubation. If needed, wash the cells again with fresh medium. For improved results, replace the medium after overnight culture.

Storage Method

Store in a cool, dark place (2–10°C).

Packaging Unit

5mL × 5
bottles

Expiration Period

Valid for 1 year

References

1. Tadao Ohno: Cell Engineering, 7(2):171-172, (1988)
2. T. Ohno, K. Kurita, S. Abe, N. Eimori and Y. Ikawa: A simple freezing medium for serum-free cultured cells. Cytotechnology 1: 257-260 (1988)