

**Blockmaster™ PA1080** [Code No. : J-PA1080]

**PRODUCT DESCRIPTION**

**Blockmaster™ PA Series** is a synthetic, virus and animal protein free blocking reagent that reduces protein and cell adsorption to solid substrates.

**Blockmaster™ PA Series** consists of our proprietary aqueous polymer with a hydrophobic unit for physical adsorption to solid surface.

**Blockmaster™ PA Series** can be used either as a substitute for BSA or along with BSA.

**Features**

- Virus and animal protein free
- Quality can be controlled (Specified molecular weight range)
- Water soluble
- Adsorbs physically to solid substrates e.g. polystyrene, glass, polydimethylsiloxane, etc.
- Prevents protein and cell adsorption to solid substrates

**Example Applications**

Blocking reagent for microfluidics, storage vessels for proteins, cell culture plates, immunoassay

**SPECIFICATIONS**

Package volume	100 mL
Solid content in solution	1 wt%
Solvent	Water containing 0.1% ProClin950 as a preservative
Appearance	Colorless or slightly yellow, transparent
Expiration date	Printed on the label

**STORAGE**

**Blockmaster™ PA Series** is stable when stored at 2-8 °C. Do not freeze the vial.

**DISPOSAL**

Observe all federal, state and local laws when considering most appropriate disposal method.

**IMPORTANT NOTICE**

This product is for research use only and not intended for therapeutic or diagnostic use.

**RECOMMENDED PROTOCOL**

**Protocol to prevent protein adsorption to substrates**

1. Add 200  $\mu$ L of **Blockmaster™ PA1080** into the wells of 96 wells plate.
  2. Incubate for 30min at room temperature.
  3. Remove **Blockmaster™ PA1080** and wash the wells with 350  $\mu$ L water 3 times.
  4. Add 100  $\mu$ L protein solution into the wells.
- Blockmaster™ PA 1080** coating prevents protein adsorption to the wells.

**Reagent and equipment requirement**

- Substrate ; 96 wells polystyrene microplate, 96wells glass microplate
- Blocking reagent ; **Blockmaster™ PA1080**
- Protein sample ; Anti-mouse IgG conjugated HRP at 200 ng/mL concentration in PBS

**Protocol to prevent cell adsorption to substrates**


1. Add 200  $\mu$ L **Blockmaster™ PA1080** into the 96 wells plate.
  2. Incubate for 30 min at room temperature.
  3. Remove **Blockmaster™ PA1080** and wash the wells with 200  $\mu$ L PBS 3 times.
  4. Add 100  $\mu$ L of cell suspension into the wells.
- Blockmaster™ PA 1080** coating prevents cell adsorption to the wells.

**Reagent and equipment requirement**


- Substrate ; 96 wells polystyrene microplate, 96wells glass microplate
- Blocking reagent ; **Blockmaster™ PA1080**
- Cell sample ; HT29 cell at  $2.5 \times 10^4$  cells/mL concentration in PBS

**REFERENCE**

**Protein adsorption to the non-treated polystyrene microplate after incubation for 1 hr**



**Cell adsorption to the non-treated polystyrene microplate after incubation for 15 hr**



**CONTACT INFORMATION**

Medical & Biological Laboratories Co., Ltd.  
 IVD Materials Department, Sales & Marketing Division  
 E-mail: [IVDMaterial@mbi.co.jp](mailto:IVDMaterial@mbi.co.jp)  
 URL: <https://www.mblbio.com/e/>