8486
Lot 008FA~

Printed November 29, 2022
Version 6.0

For Research Use Only. Not for use in diagnostic procedures.

Autophagy Watch

CODE No. 8486

For Autophagy Flux Assay and LC3 Immunocytochemistry

Description

Autophagy has been suggested to have an important role in development, differentiation, neurodegeneration, infection and cancer. Autophagy is a process of intracellular bulk degradation in which cytoplasmic components including organelles are sequestered within double-membrane vesicles. As shown in Fig. 1, a flat vesicle called an isolate membrane appears in the cytoplasm (1). Subsequently, the membrane extends while taking in the cytoplasm (2), its edges fusing to compose an AP (Autophagosome) (3). Mitochondria and other large organelles are also contained within the AP. When the AP fuses with a lysosome (4), its contents are degraded (5). Autophagosomes deliver the contents to the lysosome/vacuole for degradation. This process aids in cell survival during starvation periods by eliminating nutrients being needed for unnecessary cellular processes. .

Detecting LC3 protein by using Western blotting application is widely used for monitoring autophagic activity. LC3 is most actively studied and frequently used as a mammalian autophagy marker. Shortly after translation (proLC3), LC3 is processed at the C-terminus into LC3-I. Upon induction of autophagy, LC3-I is conjugated to the substrate phosphatidylethanolamine (PE). The PE-LC3-I conjugate is referred to as LC3-II. Despite having a higher molecular weight than LC3-I, LC3-II is more hydrophobic and shows a higher mobility in sodium dodecyl sulfate-polyacrylamide gel electrophoresis (SDS-PAGE). Although the amount of LC3-II is correlated with the number of autophagosomes, it can be difficult to interpret LC3 immunoblot analysis.

Autophagy Watch is a combination of antibodies and lysosomal inhibitors that can detect autophagy flux in Western blotting, and visualize autophagosomes in Immunocytochemistry. *Autophagy Watch* is ready to use for the rapid and quantitative monitoring of autophagic activity.



Fig. 1. Autophagy pathway and lysosomal inhibitors

Features of Autophagy Watch

1. Anti-LC3 antibody and lysosomal inhibitors that detect autophagy flux in Western blotting are included.

***About autophagy flux assay**

Autophagy flux can be measured by comparing the amount of LC3-II in the presence and absence of lysosomal inhibitors (Fig. 2). If the enhancement of autophagy flux is occurring, the amount of LC3-II will be higher in the presence of the inhibitor.

(Reference) Mizushima, N. and Yoshimori, T. *Autophagy* **3**, 542-545 (2007)

Fig. 2. Interpretation of autophagy flux

2. This kit includes HRP conjugated antibodies for Western blotting application of the autophagy flux assay [Anti-LC3 mAb HRP-Direct (code no. M186-7Y), Anti- α -Tubulin pAb HRP-Direct (code no. PM054-7Y)]. These antibodies are time saving and easy to use.
3. Inhibitors included with the kit are diluted 1,000-fold.
4. This kit includes Anti-LC3 antibody for Immunocytochemistry [Anti-LC3 mAb (code no. M152-3YY)]. This antibody is used to visualize autophagosomes.

Kit Components

Code No.	Products	Clone	Isotype	Quantity	Application	Crossreactivity
M186-7Y	Anti-LC3 mAb-HRP-Direct	8E10	Mo IgG2ak	100 μ L	WB	Hu, Mo, Rat, Hm
M152-3YY	Anti-LC3 mAb	4E12	Mo IgG1k	50 μ L, 2 mg/mL	WB (weak), IC, IP, FCM, Immuno-EM	Hu, Mo, Rat, Hm
PM054-7Y	Anti- α -Tubulin pAb-HRP-Direct	Polyclonal	Rab IgG (aff.)	100 μ L	WB	Hu, Mo, Rat, Hm, Chi
PM036-PNY	Positive control for anti-LC3 antibody			100 μ L (20 tests)		
	Chloroquine solution (x1000)			100 μ L		
	Bafilomycin A1 solution (x1000)			100 μ L		
	Cell lysis buffer (x5)			1 mL x2		

WB: Western Blotting, IC: Immunocytochemistry, IP: Immunoprecipitation, FCM: Flow Cytometry, Immuno-EM: Immuno-electron microscopy, Hu: Human, Mo: Mouse, Rab: Rabbit, Hm: Hamster, Chi: Chicken

Formulation

M152-3YY, PM036-PNY, PM054-7Y: PBS containing 50% Glycerol (pH 7.2). No preservative is contained.

M186-7Y: PBS containing 50% Glycerol (pH 7.2)/ Preservative /Stabilizer

PM036-PNY: human LC3B overexpressed 293T 1×10^5 cells were suspended with 100 μ L of Laemmli sample buffer. This product contains 5% 2-mercaptoethanol. Handle with care.

Chloroquine solution ($\times 1000$): 50 mM Chloroquine/PBS

Bafilomycin A1 solution ($\times 1000$): 100 μ M Bafilomycin A1/DMSO

Cell lysis buffer ($\times 5$): 0.25 M Tris-HCl (pH 7.5), 0.75 M NaCl, 0.25% Nonidet P-40

Storage

This kit is stable for one year from the date of purchase when stored at -20°C .

Example of protocols using this kit

1. Autophagy flux assay by Western blotting

The following protocol is an example of LC3 autophagy flux assay detected by Western blotting under starvation conditions of mammalian cells (MEF cells).

Materials required but not provided

- Nutrient medium [RPMI-1640 (Sigma; code no. R8787), 10% FCS, 100 U/mL Penicillin, 50 μ g/mL Streptomycin]
- Starvation medium [HBSS- Hanks' Balanced Salt solution (Sigma; code no. H9269)]
- 0.25% trypsin/1 mM EDTA solution (Wako; code no. 201-16945)
- PBS
- Laemmli sample buffer
- Protein concentration assay kit (*e.g.*, Bio-Rad; code no. 500-001JA)
- Skimmed milk
- PBS-T (0.05% Tween-20/PBS)
- Chemiluminescent substrate (GE Healthcare; code no. RPN2106)
- Culture dish (90 mm) (Thermo Scientific; code no. 172931)
- Cell scraper (Sumitomo Bakelite; code no. MS-93170)
- PVDF membrane (Merck Millipore; code no. IPVH00010)

*The product name/manufacturer name is reference information.

For details on each product, please check with the respective manufacturer.

Method

a) Preparation of loading samples

Preparation of reagents

Cell lysis buffer

Dilute Cell lysis buffer ($\times 5$) with 5 times its volume with distilled water. Store at 4°C.

(e.g., dilute 1 mL of Cell lysis buffer ($\times 5$) with 4 mL of distilled water.)

1) (First day)

Maintain MEF cells in a 90 mm culture dish in nutrient medium (10% FCS/RPMI-1640) until 70~80% confluence at 37°C with 5% CO₂.

2) (Second day)

Detach MEF cells from the dish using trypsin/EDTA solution and count the cells. Seed 5×10^5 cells per dish with 10 mL nutrient medium. (Prepare eight dishes.)

3) Incubate cells at 37°C with 5% CO₂ overnight.

4) (Third day)

Prepare four types of medium in 50 mL tubes as shown in Table 1.

Table 1. Preparation of culture medium for Autophagy flux assay

	Medium		Inhibitor	
a	10%FCS/RPMI-1640	20 mL	(-)	
b	HBSS	20 mL	(-)	
c	HBSS	20 mL	Chloroquine solution	20 μ L
d	HBSS	20 mL	Bafilomycin A1 solution	20 μ L

- 5) To wash the cells, remove medium supernatant from the culture dishes, add 10 mL of PBS gently, and remove PBS.
- 6) Add 10 mL of medium prepared in step 4 (a, b, c, d in Table 1) to each two dishes, respectively.
- 7) Incubate at 37°C with 5% CO₂ for 3 hours. (Autophagy induction)
- 8) Wash the cells like step 5.
- 9) Add 5 mL PBS and detach the cells from the dish using cell scraper.
- 10) Place cell suspension of a, b, c and d into 15 mL tubes and centrifuge at $700 \times g$ at 4°C for 5 minutes.
- 11) Carefully remove and discard the supernatant. *Cell pellet can be stored at -80°C.
- 12) Add 100 μ L of ice-cold Cell lysis buffer to each cell pellet and vortex for 5 seconds.
- 13) Transfer the samples to 1.5 mL sampling tubes and centrifuge at $15,000 \times g$ at 4°C for 5 minutes.
- 14) Transfer supernatants from step 13) into new 1.5 mL sampling tubes and measure total protein concentration. (e.g., Bradford method).
(Protein concentration of the sample will be approximately 1~2 mg/mL. If concentration is low, increase the starting culture dishes.)
- 15) Add Laemmli sample buffer until final protein concentration 0.5 mg/mL and boil the samples for 5 minutes. *Prepared samples can be stored at -20°C.

b) Western blotting

Reagents Preparation

- Membrane blocking buffer (10% skimmed milk/PBS)
Dissolve 5 g skimmed milk in 50 mL PBS.
- Antibody dilution buffer (1% skimmed milk/PBS)
Dissolve 0.5 g skimmed milk in 50 mL PBS.
- Antibody reaction reagent
Dilute Anti-LC3 mAb-HRP-Direct (code no. M186-7Y) (1:1,000) with 1% skimmed milk/PBS.
(e.g., dilute 5 µL antibody with 5 mL 1% skimmed milk/PBS.)
*In the case of detecting loading control, dilute Anti- α -Tubulin pAb-HRP-Direct (code no. PM054-7Y) (1:1,000) with 1% skimmed milk/PBS.

- 1) Load 20 µL of the sample per lane in 15% SDS-polyacrylamide gel.
*In the case of using positive control, load 5 µL of the PM036-PNY per lane.
- 2) Run the gel electrophoresis at 15 mA for 2 hours.
- 3) Blot the protein to a PVDF membrane at 1 mA/cm² for 1 hour in a semi-dry transfer system. (See the manufacturer's manual for precise transfer procedure.)
- 4) To reduce nonspecific binding, soak the membrane in 50 mL of 10% skimmed milk/PBS overnight at 4°C.
- 5) Wash the membrane with PBS-T (0.05% Tween-20/PBS) for 5 minutes. Repeat 3 times.
- 6) Incubate the membrane with antibody diluted with 1% skimmed milk for 1 hour at room temperature. *Do not allow the membrane to dry during incubation.
- 7) Wash the membrane as described in step 5.
- 8) Detect LC3 band using chemiluminescence substrate.

To detect loading control using PM054-7Y (Anti- α -Tubulin pAb-HRP-Direct), soak the used membrane in 50 mL 10% skimmed milk/PBS for 1 hour at room temperature or overnight at 4°C and repeat procedures step 5 to step 8.

Results

Fig. 3 shows the results of autophagy flux assay by Western blotting using the above protocol. Starvation (HBSS culture) induces increasing the amount of LC3-II (Lane 1 vs. Lane 2).

Additionally, the amount of LC3-II increases in the presence of Chloroquine or Bafilomycin A1 (Lane 2 vs. Lane 3, 4). These results suggest the accumulation of LC3-II under starvation condition indicates the enhancement of the autophagic flux.


Fig. 3. Results of autophagy flux assay (Cells: MEF)

*MEF Atg5^{-/-} was kindly provided by Dr. Noboru Mizushima. (Department of Biochemistry and Molecular Biology, Graduate School and Faculty of Medicine, The University of Tokyo)

2. Immunocytochemistry

The following protocol is the example of observing autophagosomes by LC3 Immunocytochemistry under starvation conditions of mammalian cells (MEF cells).

Materials required but not provided

- Nutrient medium
[RPMI-1640 (Sigma; code no. R8787), 10% FCS, Penicillin 100 U/mL, Streptomycin 50 µg/mL]
- Starvation medium [HBSS- Hanks' Balanced Salt solution (Sigma; code no. H9269)]
- 0.25% trypsin/1 mM EDTA solution (Wako; code no. 201-16945)
- PBS
- 4% Paraformaldehyde (PFA)/PB Solution (Wako; code no.163-20145)
- Digitonin (*e.g.* Sigma; code no. D141-100MG)
- DAPI solution (*e.g.* Dojindo; code no. 340-07971 diluting 1000-fold in PBS)
- Mounting medium (Life Technologies; code no. P36930)
- Secondary antibody [Anti-Mouse IgG Alexa Fluor® 488 conjugate (Life Technologies; code no. A11001)]
- 8 wells chamber slide (Thermo Scientific; code no. 154941)

*The product name/manufacturer name is reference information.

For details on each product, please check with the respective manufacturer.

Preparation of reagents

- Cell permeabilization solution
Prepare 100 µg/mL Digitonin/PBS just before use.
(*e.g.*, dissolve 1 mg Digitonin in 10 mL PBS.)
- Primary antibody reaction solution
Dilute anti-LC3 mAb (code no. M152-3YY) (1:50) with PBS.
(*e.g.*, dilute 20 µL antibody in 1 mL PBS.)
- Secondary antibody reaction solution
Dilute anti-Mouse IgG Alexa Fluor® 488 conjugate (1:500) with PBS.
(*e.g.*, dilute 2 µL antibody in 1 mL PBS.)

Method

1) (First day)

Maintain MEF cells in a 90 mm culture dish in nutrient medium (10% FCS/RPMI-1640) until 70-80% confluence at 37°C with 5% CO₂.

2) (Second day)

Detach MEF cells from the dish using trypsin/EDTA solution and count the cells.

Seed 5×10³ cells per wells with 200 µL nutrient medium.

3) Incubate at 37°C with 5% CO₂ overnight.

4) (Third day)

Prepare four types of medium in 1.5 mL sampling tubes as shown in Table 2.


Table 2. Preparation of culture medium for Immunocytochemistry.

	Medium		Inhibitor	
a	10%FCS/RPMI-1640	1 mL	(-)	
b	HBSS	1 mL	(-)	
c	HBSS	1 mL	Chloroquine solution	1 µL
d	HBSS	1 mL	Bafilomycin A1 solution	1 µL

- 5) Remove the culture supernatant by careful aspiration. Do not allow the slides to dry during steps 5-17.
- 6) To wash the cells on slides, add 200 μ L HBSS in each well and remove HBSS.
- 7) Add 200 μ L medium prepared in step 4 (a, b, c, d in Table 2) to each well, respectively.
- 8) Incubate at 37°C with 5% CO₂ for 3 hours.
- 9) Remove supernatant and wash gently the cells with PBS. Repeat 2 times.
- 10) Fix the cells adding 200 μ L of 4% PFA/PB solution to each well for 10 minutes at room temperature.
- 11) Wash the cells as in the step 9.
- 12) Permeabilize the cells with 200 μ L of 100 μ g/mL Digitonin/PBS for 10 minutes at room temperature.
- 13) Wash the cells as described in the step 9.
- 12) Add 200 μ L of primary antibody reaction reagent for 1 hour at room temperature.
- 14) Wash the cells as described in the step 9.
- 15) Add 200 μ L of secondary antibody reaction reagent and incubate for 30 minutes at room temperature. Keep out light.
- 16) Add a drop of DAPI solution and incubate for 5 minutes at room temperature.
- 17) Wash the cells as described in the step 9.
- 18) Add mounting medium onto the slide, then put a cover slip on it and image the slide with fluorescence microscope.

Results

Fig. 4 shows the results of LC3 Immunocytochemistry using the above protocol. Under starvation (HBSS culture) condition, the accumulation of autophagosomes is detected. In the presence of an inhibitor, the number of autophagosomes is further increased.

Microscope: BZ-9000 Generation II (Keyence)
 Primary antibody: Anti-LC3 mAb (code no. M152-3YY)
 Secondary antibody: Anti-Mouse IgG Alexa Fluor® 488 conjugate

Fig. 4. Results of LC3 Immunocytochemistry (Cells: MEF)

Troubleshooting

1. It is difficult to separate the bands of LC3-I and LC3-II using Western blotting.
LC3 are detected as two bands at 16 kDa and 14 kDa in various cell lysates, respectively. We recommend using 15% SDS-polyacrylamide gel to obtain optimal performance.
2. LC3-I and LC3-II band cannot be detected.
Consider increasing the volume of loading sample or increasing the concentration of anti-LC3 antibody. In addition, we recommend using Western blotting control for anti-LC3 antibody (code no. PM036-PNY).
3. During the reaction of immunocytochemistry, the cells come off from the slide wells.
Never allow the slide to dry out.

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