Labiase

Cat#OZ-30EX

(from *Streptomyces fulvissimus* TU-6)

Labiase, produced by a submerged culture of *Streptomyces fulvissimus* TU-6, is a new enzyme preparation that lyses effectively cell walls of numerous lactic acid bacterium. This preparation is used as a tool for studies of cell walls structure of lactic acid bacterium and the preparation of plasmid DNA, intracellular enzyme and protoplast from lactic acid bacterium and so on.

SPECIFICATIONS:

Appearance.....Lyophilized powder
(containing lactose)
Activity..... ≥ 10 units/g
(β -N-Acetyl-D-glucosaminidase)
Other activities contained
Lysozyme..... $\geq 2 \times 10^4$ units/g
Optimum pH and temperature.....See Fig. 1,2
pH and thermal stabilitySee Fig. 1,2

ASSAY FOR β -N-ACETYL-D-GLUCOSAMINIDASE ACTIVITY: Unit Definition

One unit of enzyme catalyzes the release of 1 μ mole of *p*-nitrophenol from *p*-nitrophenyl N-acetyl- β -D-glucosaminide per minute at 37°C, pH4.0.

Method

Reaction mixture

Substrate.....5mM <i>p</i> -nitrophenyl N-acetyl- β -D-glucosaminide in 0.2M Citrate-phosphate buffer, pH 4.0	0.95 ml
Enzyme.....suitably diluted enzyme	0.05 ml
Total volume	1.00 ml

Procedure

After incubation for 10 minutes at 37°C, terminate by 2ml of 1M Na₂CO₃. Read A₄₀₅.

Calculation

$$\text{Enzyme unit} = \frac{A_{405} \times 3}{17.7 \times 10^3 \times 10} \times 10^3$$

The molar extinction coefficient of free *p*-nitrophenol under these conditions is 17,700.

ASSAY FOR LYTIC ACTIVITY TOWARD LACTIC ACID BACTERIUM: Method

Reaction mixture

Buffer.....50mM Citrate -sodium citrate buffer, pH 4.0	1ml
Substrate...Lactic acid bacterium suspension(A ₆₆₀ ≐5.0)	1ml
Enzyme.....500 μ g-protein/ml solution	1ml
Distilled water	2ml
Total volume	5ml

Procedure

(Before adding the enzyme solution preincubation is carried out for 5 minutes at 37°C.)After incubation for 2 hours at 37°C with gentle shaking, A₆₆₀ of the mixture is determined. As a reference, 1 ml of distilled water is used instead of enzyme solution.

Calculation

$$\text{Percentage decrease in } A_{660} = \frac{(d_0 - d_t) - (D_0 - D_t)}{d_0} \times 100$$

d(0 or t):A₆₆₀ of reaction mixture after 0 or 2 hours
D(0 or t):A₆₆₀ of reference mixture after 0 or 2 hours

Result

See Table 1. The extent of lysis of lactic acid bacterium cells by **Labiase** varies with lactic acid bacterium strain, growth stage of lactic acid bacterium, or cultural condition.

STORAGE:

Lyophilized preparation is stable for at least 1 year when stored at 4°C.


Fig.1 pH-stability and activity

Stability: 18 hours treatment at 25°C

Activity: 37°C


Fig.2 Thermal-stability and activity

Stability: 10 minutes treatment

Activity: pH4.0

Table 1. Lysis of lactic acid bacterium cells by Labiase

Strains	Lysis (%)
<i>Leuconostoc mesenteroides</i> IFO 3832	92.8
<i>Leuconostoc lactis</i> IFO 12455	99.6
<i>Streptococcus faecalis</i> IFO 3971	96.5
<i>Lactobacillus acidophilus</i> IFO 13951	97.9
<i>Lactobacillus plantarum</i> IAM 1041	96.8
<i>Lactobacillus casei subsp. casei</i> IFO 3533	59.4
<i>Lactobacillus casei subsp. ramosus</i> IAM 1118	91.8
<i>Lactobacillus brevis</i> IFO 3345	95.6
<i>Lactobacillus fermentum</i> IFO 3071	70.4
<i>Lactobacillus sp. (buchneri)</i> IFO 3961	89.4
<i>Lactobacillus fructivorans</i> IFO 13954	96.9
<i>Lactobacillus delbrueckii subsp. bulgaricus</i> JCM 1002	97.7
<i>Lactobacillus delbrueckii subsp. delbrueckii</i> JCM 1012	92.0
<i>Lactobacillus helveticus</i> JCM 1120	98.7
<i>Lactobacillus helveticus</i> JCM 1005	84.5
<i>Lactobacillus kefir</i> JCM 5818	93.0
<i>Lactobacillus sanfrancisco</i> JCM 5668	82.1
<i>Bifidobacterium bifidum</i> JCM 1255	62.9
<i>Streptococcus salivaris subsp. thermophilus</i> IFO 13957	94.0
<i>Streptococcus mutans</i> JCM 5705	73.9

Please note that the above data is an example. The values vary widely depending on the condition.

Manufactured by  OZEKI CORPORATION

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